

## High-Performance Cork for Vapor Shippers

### Models Affected

11037386*	CRYOSHIPPER MINI
10508967*	CRYOSHIPPER
11015195*	CRYOSHIPPER XC
10777411	IATA CRYOSHIPPER
21007700	CT-50 CRYOSHIPPER
21000341	CT-250 CRYOSHIPPER

\* including derivative part numbers for same tank



Figure 1 Regular Cork (Left) vs. High Performance Cork (Right)

Chart MVE introduces the new High-Performance Cork into select Vapor Shipper models. The new High-Performance Cork was designed to protect shipments that, due to unforeseen circumstances, are shipped on their side. Chart does NOT recommend laying a tank on its side, and we do not list the horizontal performance as a spec as this is another continuous improvement over the previous design. The High-Performance cork includes a felt-like hydrophobic material to prevent LN2 vapor from escaping the tank in non-upright positions.

When a vapor shipper is shipped on its side, it can reduce the tank's static hold time by 90% or more. For example, Cryoshipper XC may last only one day in the horizontal configuration with the previous cork but with the High-Performance Cork it extends the static hold time to approximately one week. Out-of-position orientation of the various models with the high-performance cork will see an improvement in holding time, with that improvement being lesser as the neck tube diameter decreases.

The High-Performance Cork now comes standard on vapor shippers (implemented Apr 9<sup>th</sup>, 2019). High-Performance Corks can be retrofitted with old vapor shipper models utilizing the part numbers below.




### High-Performance (HP) Lid Kit – Part Numbers

20552556	CRYOSHIPPER MINI HP LID KIT
20553315	CRYOSHIPPER, CRYOSHIPPER XC, IATA CRYOSHIPPER HP LID KIT
21007703	CT-50 CRYOSHIPPER HP LID KIT
21003472	CT-250 CRYOSHIPPER HP LID KIT

Always wipe condensate and water from the outside of the cork and the inside of the neck tube before inserting the cork into the dewar.

### High Performance Cork Temperature Mapping

\*NOTE: – Actual field results will vary.



**\*NOTE: – Actual field results will vary.**