MVE Operating Instruction Manual for Vapor Shipper Dewars

ADVANCED
QWick Charge Technology

Charges in under 2 hours

Table of Contents	Page
English	2-7
English (SI Version)	8-13
German	14-19
Spanish	20-25
Italian	26-31
French	32-37

 **CHART INC.**
2200 Airport Industrial Dr., Ste 500
Ball Ground, GA 30107
U.S.A
www.chartbiomed.com

Chart Australia Pty Ltd.
ABN 21 075 909 410
Sydney Business & Technology Centre
Unit 43 / 2 Railway Parade
Lidcombe, NSW 2141, Australia

Chart BioMedical Ltd.
Unit 6 – Ashville Way
Wokingham RG41 2PL
United Kingdom

EC	REP
----	-----

 Medical Product Services
Borngasse 20
35619 Braunfels, Germany

Ref 11562640 Rev P 05/16

A Chart Industries Company

MVE Vapor Shippers are ideally suited for the transportation of cryobiological materials. The liquid nitrogen is retained in an absorbent material and the cold nitrogen vapors maintain the Cryo Storage Area between -150°C and -190°C. The absorbent retains the liquid nitrogen and prevents accidental spills.

Liquid nitrogen is classified as "Dangerous Goods" by the Department of Transportation (DOT). At least three regulatory agencies publish dangerous goods regulations:

- DOT - Department of Transportation
- IATA - International Air Transportation Association
- ICAO - International Civil Aviation Organization

If used outside of the USA, please consult your applicable regulatory agencies. Liquid nitrogen is classified under the name "Nitrogen, refrigerated liquid" and has a UN 1977 code number. The name classification and UN Code Number must be indicated on the shipping carton in addition to a "Gas nonflammable" label. Packing and labeling requirements can be obtained from the above organizations. Also, contact your carrier for any variances in the rules as they may apply specifically to that company and/or the ultimate destination of the shipment.

Be aware that the regulations undergo periodic revisions. If you need current information prior to the shipping of your cryobiological specimen, please contact your air carrier for the most recent regulations.

GENERAL DESCRIPTION:

The vapor series dewar is a vacuum insulated container of aluminum with fiberglass neck construction providing the highest efficiency possible in nitrogen vapor storage. Use the container for inert fluids only. Liquid oxygen is not compatible with this unit and must not be stored inside the Vapor Shipper.

These high quality vacuum insulated units are constructed of durable material, compatible with the divergent temperature extremes and broad applications of cryobiology. The life expectancy of this product is five (5) years.

Upon receipt of the container, examine both the container and packaging for any evidence of damage during shipping. It is normal for some dust or white powder-like residue to accumulate at the bottom of the unit due to shipping vibrations. If desired, the owner may wipe or vacuum it out. Contact the shipping agent if any damage is apparent. Watch after the first fill for any signs of vacuum loss, such as excessive frost or sweating on the outside jacket. (Some frost near the top just after filling is normal).

A) SAFETY

To avoid injury by frostbite, use extreme care whenever handling liquid nitrogen, liquid nitrogen storage or transfer vessels or any objects which have come in contact with liquid nitrogen.

- Leave no area of skin exposed.
- Always wear proper safety attire over clothing: face shield, cryogenic gloves, cryogenic apron.
- Never overfill vapor shippers with liquid nitrogen. Overfilling the tank may cause immediate or premature vacuum failure to occur.
- Always keep vapor shippers in upright position.
- Do not tightly seal liquid nitrogen container or prevent nitrogen gas from escaping.
- Use extreme care to prevent spilling and splashing liquid nitrogen during transfer.
- Immediately remove any clothing or safety attire on which liquid nitrogen has spilled.
- Get immediate medical attention for any frostbite injuries due to liquid nitrogen.
- Never ship Vapor Shipper on its side or upside down. This can lead to vacuum failure and loss of product inside the Vapor Shipper.

WARNING: The venting of nitrogen vapors will deplete oxygen in the air, possibly leading to asphyxiation or even death. Do not store or use container in areas that have poor ventilation.

Name and address of manufacturer

Authorized representative in the European Community

B) FILLING INSTRUCTIONS

Note: Fill the container with a funnel or transfer line when possible. Avoid spilling liquid nitrogen over the vacuum cap near the neck as this can shrink the seal and allow air to leak into the vacuum space causing premature vacuum failure. To ensure maximum performance from your MVE vapor dewar simply follow the listed steps just prior to shipping to final destination:

1. Open container that dewar is in and remove cork/cover by lifting straight up (do not twist).
2. Fill unit to bottom of neck tube.
 - a) Follow established safety practices and procedures for transferring LN₂.
 - b) Transfer using LN₂ hose with phase separator or pouring container using a funnel.
3. Replace cork/cover and allow unit to stand for 2 hours (cooling down unit).
4. Pour off excess liquid just prior to shipment.
5. Weigh unit and record.

C) PREVENTIVE MAINTENANCE

The owner may follow two methods of maintenance, Continuous or Annual.

For Continuous Maintenance:

Follow steps 1 - 4 under Filling Instructions above but allow the unit to stand undisturbed for 2 hours. Verify no major frost or condensation on outside of container, which would indicate either a weak or no vacuum.

For Annual Maintenance:

1. Weigh empty unit and record.
2. Fill unit to bottom of neck tube.

Note: Follow established safety practices and procedures for transferring LN₂.

3. Transfer LN₂ using cryogenic hose with a phase separator, pouring container using a funnel or other approved filling device.
4. Replace cork/cover and allow unit to stand for 24 hours.
5. Weigh unit.
6. Allow filled unit to sit for another 24 (+/-0.25) hours.
7. Weigh second time.
8. Calculate evaporation rate. The difference between the first weight and the second weight is the daily evaporation rate in lbs. This figure roughly signifies the normal evaporation rate (N.E.R.) (first weight - second weight) x 0.5606 = liters/day.
9. Pour out LN₂. Take weight to determine 'charged weight'.

NOTE: During this time verify no major frost or condensation on outside of container, which would indicate either a weak or no vacuum.

10. Calculate Hold Time. Subtract the Empty Weight from the Charged Weight and divide by 1.782; then divide by NER. [(Charged Weight - Empty Weight) ÷ 1.782] ÷ NER). Remember that the Empty Weight must be taken before filling, while the inside of the Vapor Shipper is at room temperature.

MVE vapor shippers were primarily designed as vapor shipping containers; however, they can also be used for immersion of samples. A sharp blow to the outer vessel can damage the neck tube or start a vacuum leak. Use caution and common sense in handling the container. Use the following weight table as a general guide to determine if your vapor shipper is fully charged.

Model	Empty Weight		Suggested Wt. Charged		Static Hold Time Days
	lbs	(kg)	lbs	(kg)	
SC 2/1V	6	(2.7)	8	(3.6)	8
SC 4/2V	10	(4.5)	17.1	(7.8)	13
SC 4/3V	11.6	(5.3)	19.1	(8.7)	21
Cryoshipper XC/IATA	33	(15.0)	49	(22.2)	14
Cryoshipper	26.1	(11.8)	39	(17.7)	10
Cryomoover	30.5	(13.8)	40.9	(18.6)	12
Mini-Moover	9	(4.1)	13.7	(6.2)	14
XC 20/3V	24.1	(10.9)	34	(15.4)	16
SC 20/12V	25.7	(11.7)	41.3	(18.7)	85
Doble Models					
Doble 11	14.5	(6.6)	18.5	(8.4)	17
Doble 20	22.7	(10.3)	29.7	(13.5)	21
Doble 28	33	(15.0)	47	(21.3)	24
Doble 34	34	(15.4)	50	(22.7)	20
Doble 47	40	(18.1)	53.5	(24.3)	21
Doble 47-10	40	(18.1)	53.5	(24.3)	21

D) SHIPPING INSTRUCTIONS

Fill the dewar per Section 'B' and allow unit to stand for 2 hours (cooling down unit) with refills as needed.

Prior to packaging dewar for shipment the entire contents of liquid nitrogen must be removed. This is usually done by pouring out excess liquid nitrogen until no liquid is visible on the bottom of the inner dewar. Place in bio samples, package dewar and ship. The plastic shipping container provided by MVE/Chart is recommended to help keep dewar in upright position. **SHIPPING THE DEWAR IN ANY POSITION OTHER THAN UPRIGHT COULD REDUCE THE HOLDING TIME TO LESS THAN 10% OF NORMAL HOLD TIME. DO NOT SHIP UNIT ON ITS SIDE OR UPSIDE DOWN.** Remember that all MVE vapor shippers can also be used for liquid nitrogen storage as well, so it is imperative that all liquid nitrogen be removed so dewar remains classified as a vapor shipper. If liquid nitrogen is visible in the bottom of the inner it then becomes a liquid shipper and the exception status is void. The liquid inside is now classified as hazardous material.

E) UTILIZING THE QWICK CHARGE TECHNOLOGY

To utilize the QWick Charge Technology for faster charging times, do not allow inside of the dewar to warm up above -150°C between shipments. If the inside of the dewar does not warm up above -150°C, it will not need to be cooled down again before re-charging. Refill the dewar to the bottom of the neck with liquid nitrogen. Allow dewar to stand until liquid nitrogen stops boiling rapidly. After liquid nitrogen has stopped boiling rapidly, the dewar is fully charged. Pour off excess liquid nitrogen before shipment. The charged weight may be verified against the table in Section 'C' to ensure the dewar is fully charged. Please note that the Charged Weight is not an indication of unit performance.

Sanitizing and Decontaminating:

Any cleaning solution that does not react with aluminum or stainless steel can be used in the sanitation process of these dewars. In most cases, any household detergent or mild soap solution is suitable. Other cleaners and disinfectants that can be used safely include hydrogen peroxide, chlorine/water mixture, and denatured alcohol. The generally accepted practice of using 10% chlorine bleach with 90% water solution still holds as the best method for decontamination. **NOTE: DO NOT USE ANY PERTROLEUM BASED CLEANING SOLUTION.** It is important that all surfaces being sanitized are thoroughly rinsed and that all cleaner solution residues are removed after cleaning. For Vapor Shipper units, it is recommended that the unit is filled to its full capacity with the cleaning solution mixture, agitated, and then thoroughly rinsed. Allow the unit to dry completely before putting into service. It is suggested that the unit is inverted to drain and dry completely.

LETTER EXPLAINING EXCEPTION

This is concerning the applicability of the USA Federal Hazardous Material Regulations to the shipment of refrigerated samples in the "Dry Shipper" container. A "Dry Shipper" package consists of an outer container that is lined with an absorbent material. The container is charged with nitrogen refrigerated liquid which is absorbed into the container lining. The charged, completed package serves as refrigerated container for the shipment of samples.

In consideration of the above, consultation with the Research and Special Programs Administration of the DOT has determined that the use of nitrogen refrigerated liquid charged "dry shipper" containers for the shipment of samples fall within the regulation exception provided in 49CFR 173.320 paragraph (a) of the section states the requirements of this subchapter do not apply to atmospheric gases and helium when used in the operation of the process system such as refrigeration system. Paragraph (c) of 173.320 pertains to air transport of same refrigeration system. For exception status of air shipments please refer to IATA-Dangerous Goods Regulations for nitrogen refrigerated liquid. This falls in the class of 2.2 nonflammable gas, packing instructions 202 with special provisions A-152. If shipping outside of the USA, please consult your applicable regulatory agencies. For answers to questions regarding shipping regulations contact a Chart, AI-Cryobiological Tech Service Representative.

Limited Warranty Provisions:

Purchaser agrees that before this limited warranty shall become effective, Purchaser shall fully inspect each New Product, Repair/Blemished Product, or Service Part within ten (10) days of delivery and before such Product is put to use. Purchaser also agrees to operate the New Product, Repair/Blemished Product, or Service Part in accordance with CHART's operating instructions and that failure to do so shall void this limited warranty. Purchaser further agrees that any claim for breach of warranty must be made in writing within 60 days of discovery of a purported defect. CHART will not be responsible for any alleged breach of warranty, which, as a result of CHART's inspection, CHART determines to have arisen from a cause not covered by this limited warranty. In this case, CHART will charge the purchaser a nominal fee to repair the unit.

This limited warranty does not apply to: (A) Normal routine service items; (B) Repair or replacement necessitated by misuse, abuse, accident, or repairs made by persons other than CHART or persons not authorized by CHART, (C) Use of external equipment or parts with the New Product, Repair/Blemished Product, or Service Parts other than those approved by CHART, (D) Defects caused by effects of normal wear and tear; and (E) Acts of God, or other causes not within the control of CHART.

If Purchaser believes that a New Product, Repair/Blemished Product, or Service Part does not comply with the limited warranty stated above, Purchaser should contact CHART at the address stated above, describing the problem and providing proof of the date of purchase. If directed by CHART, Purchaser shall return the New Product, Repair/Blemished Product, or Service Part freight prepaid, properly packaged in a CHART approved shipping container and properly identified by a Return Material Authorization Number issued by CHART. New Products, Repair/Blemished Product, or Service Parts returned without a Return Material Authorization Number will be refused and returned at Purchaser's expense.

The remedies available for any breach of this limited warranty are limited to repair or replacement of the defective New Product, Repair/Blemished Product, Service Part, or refund of the purchase price, at the sole discretion of CHART. CHART warrants that replacement or repaired New Product, Repair/Blemished Product, or Service Part shall be free from defects in material and workmanship for the duration of the unexpired portion of the original warranty or ninety (90) days from the date of re-shipment to Purchaser, whichever is longer. PURCHASER'S RECOVERY FROM CHART FOR ANY CLAIM SHALL NOT EXCEED PURCHASER'S PURCHASE PRICE FOR THE NEW PRODUCT, REPAIR/BLEMISHED PRODUCT, OR SERVICE PART GIVING RISE TO SUCH CLAIM, IRRESPECTIVE OF THE NATURE OF THE CLAIM, WHETHER IN CONTRACT, TORT, WARRANTY, NEGLIGENCE, STRICT LIABILITY OR OTHERWISE. CHART SHALL NOT BE LIABLE FOR AND PURCHASER SHALL INDEMNIFY, DEFEND AND HOLD CHART HARMLESS FROM ANY CLAIMS BASED ON CHART'S COMPLIANCE WITH PURCHASER'S DESIGNS, SPECIFICATIONS OR INSTRUCTIONS, OR MODIFICATION OF ANY PRODUCTS BY PARTIES OTHER THAN CHART, OR USE IN COMBINATION WITH OTHER PRODUCTS.

PURCHASER SHALL NOT IN ANY EVENT BE ENTITLED TO, AND CHART SHALL NOT BE LIABLE FOR, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY NATURE INCLUDING, WITHOUT LIMITATION, BUSINESS INTERRUPTION COSTS, REMOVAL AND/OR REINSTALLATION COSTS, REPROCUREMENT COSTS, LOSS OF PROFIT OR REVENUE, LOSS OF DATA, PROMOTIONAL OR MANUFACTURING EXPENSES, OVERHEAD, INJURY TO REPUTATION OR LOSS OF CUSTOMERS, EVEN IF CHART HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

EXCEPT FOR THIS LIMITED WARRANTY, CHART HAS MADE NO WARRANTIES OR REPRESENTATIONS, EXPRESSED OR IMPLIED, AND HEREBY DISCLAIMS ANY OTHER WARRANTIES, INCLUDING, BUT NOT LIMITED TO, IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. NO REPRESENTATION OR STATEMENT OF CHART MAY CHANGE OR ALTER THIS LIMITED WARRANTY.

Any claims for breach of this limited warranty shall be governed by Georgia law and without regard to conflict of law rules thereof and must be brought in a state or federal court in Georgia.

Some states do not allow limitations on implied warranties or on incidental or consequential damages, so the above limitations may not apply. This limited warranty gives the Purchaser specific legal rights. Purchaser may also have other rights, which vary from state to state.

Limited Warranty Statement

Chart Inc., 2200 Airport Industrial Dr., Suite 500, Ball Ground, GA 30107

General Terms:

"**New Product**" – Chart Inc. ("CHART") warrants to the original purchaser ("Purchaser") that each new Vapor Shipper Dewar (collectively, the "New Products"), shall be free from defects in materials and workmanship for a period of two (2) years from the date of shipment, except as provided in the Limited Warranty Provisions below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship as for a period of three (3) years from the date of shipment, except as provided in the Limited Warranty Provisions below.

"**Repair/Blemished Product**" – CHART warrants to the Purchaser that all repaired equipment and factory blemished equipment (collectively, the "Repair/Blemished Products") shall be free of defects in materials and workmanship for a period of (90) days from the date of shipment, except as provided below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship for a period of (1) year from date of shipment, except as provided below.

"**Service Part**" – CHART warrants to the Purchaser that all Service Parts and Accessory Items (each a "Service Part") are warranted for (90) days from date of shipment to be free of defects in material and workmanship, except as provided below.

**Supplement to Chart Inc.
Limited Warranty Statement for Sales in Australia (AI)**

General Terms:

"New Product" – Chart Inc. ("CHART") warrants to the original purchaser only (meaning the party from whom CHART accepts the purchase order), whether such purchaser is a wholesaler, distributor or end user ("Purchaser"), that each new Liquid Nitrogen Dewar (collectively, the "New Products"), shall be free from defects in materials and workmanship for a period of two (2) years from the date of shipment, except as provided below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship as listed below, except as provided in the Limited Warranty Provisions below: Spectrum and CryoSystem Series, ET11, ET24, A-200, Workhorse 2003: 3 years / XC, SC, and Lab Series, ET 23-JLM: 5 years.

"Repair/Blemished Product" – CHART warrants to the Purchaser that all repaired equipment and factory blemished equipment (collectively, the "Repair/Blemished Products") shall be free of defects in materials and workmanship for a period of (90) days from the date of shipment, except as provided below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship for a period of (1) year from date of shipment, except as provided below.

"Service Part" - CHART warrants to the Purchaser that all Service Parts and Accessory Items (each a "Service Part") are warranted for (90) days from date of shipment to be free of defects in material and workmanship, except as provided below.

Sales of Product in Australia:

Where the Purchaser is the end user of the Product, the following additional information is provided in relation to the warranty against defects (the Warranty) provided by CHART (as set out in the Limited Warranty Statement).

1. In circumstances where a New Product, Repair/Blemished Product or Service Part does not comply with the Warranty, CHART will at its expense repair or replace that defective New Product, Repair/Blemished Product or Service Part, or refund the purchase price at the sole discretion of CHART. CHART warrants that the replacement or repaired New Product, Repair/Blemished Product or Service Part shall be free from defects in materials and workmanship for the duration of the unexpired portion of the Warranty, or 90 days from the date of reshipment to the purchaser, whichever is longer.
2. The Warranty is subject to the Limited Warranty Provisions set out below and all exclusions and limitations set forth in the CHART Limited Warranty Statement.
3. To claim the Warranty, the Purchaser must contact CHART's affiliate in Australia for further instruction, including the return of the defective New Product Repair/Blemished Product or Service Part (at the Purchaser's expense), to the following address:

Chart Australia Pty Ltd.
Unit 43/2 Railway Parade
Lidcombe, NSW 2141, Australia
Attention: Customer Service Department
Phone number: (02) 9749-4333
Email address: customerservice.australia@chartindustries.com

4. The Purchaser must otherwise comply with the requirements set out in the Limited Warranty Statement when making a claim under the Warranty.
5. Notwithstanding the Limited Warranty Statement, our Products come with guarantees that cannot be excluded under the Australian Consumer Law (ACL). You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the Products repaired or replaced if the Products fail to be of acceptable quality and the failure does not amount to a major failure.
6. The provider of the Warranty is CHART. Relevant contact details are set out below for CHART (to claim the Warranty, please use the contact details provided in paragraph 3 above):
 - (a) Customer Service Department
 - (b) Business Address:
Chart Inc.
2200 Airport Industrial Drive, Suite 500
Ball Ground, GA 30107 (USA)
 - (c) Phone number: 770-721-7759
 - (d) Email address: customerservice.usa@chart-ind.com

Repair of Products Sold in Australia

Where the Purchaser is the end user of the Product, the following additional information is provided in relation to the repair of Products by CHART:

The Products may contain "user generated data" as defined under the Australian Consumer Law (ACL). Repair of Products may result in the loss of user generated data.

MVE Vapor Shippers are ideally suited for the transportation of cryobiological materials. The liquid nitrogen is retained in an absorbent material and the cold nitrogen vapors maintain the Cryo Storage Area between -150°C and -190°C. The absorbent retains the liquid nitrogen and prevents accidental spills.

Liquid nitrogen is classified as "Dangerous Goods" by the Department of Transportation (DOT). At least three regulatory agencies publish dangerous goods regulations:

- DOT - Department of Transportation
- IATA - International Air Transportation Association
- ICAO - International Civil Aviation Organization

If used outside of the USA, please consult your applicable regulatory agencies. Liquid nitrogen is classified under the name "Nitrogen, refrigerated liquid" and has a UN 1977 code number. The name classification and UN Code Number must be indicated on the shipping carton in addition to a "Gas non-flammable" label. Packing and labeling requirements can be obtained from the above organizations. Also, contact your carrier for any variances in the rules as they may apply specifically to that company and/or the ultimate destination of the shipment.

Be aware that the regulations undergo periodic revisions.

If you need current information prior to the shipping of your cryobiological specimen, please contact your air carrier for the most recent regulations.

GENERAL DESCRIPTION:

The Vapor series dewar is a vacuum insulated container of aluminum with fiberglass neck construction providing the highest efficiency possible in nitrogen vapor storage. Use the container for inert fluids only. Liquid oxygen is not compatible with this unit and must not be stored inside the Vapor Shipper.

These high quality vacuum insulated units are constructed of durable material, compatible with the divergent temperature extremes and broad applications of cryobiology. The life expectancy of this product is five (5) years.

Upon receipt of the container, examine both the container and packaging for any evidence of damage during shipping. It is normal for some dust or white powder-like residue to accumulate at the bottom of the unit due to shipping vibrations. If desired, the owner may wipe or vacuum it out. Contact the shipping agent if any damage is apparent. Watch after the first fill for any signs of vacuum loss, such as excessive frost or sweating on the outside jacket. (Some frost near the top just after filling is normal).

A) SAFETY

To avoid injury by frostbite, use extreme care whenever handling liquid nitrogen, liquid nitrogen storage or transfer vessels or any objects which have come in contact with liquid nitrogen.

- Leave no area of skin exposed.
- Always wear proper safety attire over clothing: face shield, cryogenic gloves, cryogenic apron.
- Never overfill vapor shippers with liquid nitrogen. Overfilling the tank may cause immediate or premature vacuum failure to occur.
- Always keep vapor shippers in upright position.
- Do not tightly seal liquid nitrogen container or prevent nitrogen gas from escaping.
- Use extreme care to prevent spilling and splashing liquid nitrogen during transfer.
- Immediately remove any clothing or safety attire on which liquid nitrogen has spilled.
- Get immediate medical attention for any frostbite injuries due to liquid nitrogen.
- Never ship Vapor Shipper on its side or upside down. This can lead to vacuum failure and loss of product inside the Vapor Shipper.

WARNING: The venting of nitrogen vapors will deplete oxygen in the air, possibly leading to asphyxiation or even death. Do not store or use container in areas that have poor ventilation.

Name and address of manufacturer

Authorized representative in the European Community

B) FILLING INSTRUCTIONS

Note: Fill the container with a funnel or transfer line when possible. Avoid spilling liquid nitrogen over the vacuum cap near the neck as this can shrink the seal and allow air to leak into the vacuum space causing premature vacuum failure. To ensure maximum performance from your MVE vapor dewar simply follow the listed steps just prior to shipping to final destination:

1. Open container that dewar is in and remove cork/cover by lifting straight up (do not twist).
2. Fill unit to bottom of neck tube.
 - a) Follow established safety practices and procedures for transferring LN₂.
 - b) Transfer using LN₂ hose with phase separator or pouring container using a funnel.
3. Replace cork/cover and allow unit to stand for 2 hours (cooling down unit).
4. Pour off excess liquid just prior to shipment.
5. Weigh unit and record.

C) PREVENTIVE MAINTENANCE

The owner may follow two methods of maintenance, Continuous or Annual.

For Continuous Maintenance:

Follow steps 1 - 4 under Filling Instructions above but allow the unit to stand undisturbed for 2 hours. Verify no major frost or condensation on outside of container, which would indicate either a weak or no vacuum.

For Annual Maintenance:

1. Weigh empty unit and record.
2. Fill unit to bottom of neck tube.

Note: Follow established safety practices and procedures for transferring LN₂.

3. Transfer LN₂ using cryogenic hose with a phase separator, pouring container using a funnel or other approved filling device.
4. Replace cork/cover and allow unit to stand for 24 hours.
5. Weigh unit.
6. Allow filled unit to sit for another 24 (+/-0.25) hours.
7. Weigh second time.
8. Calculate evaporation rate. The difference between the first weight and the second weight is the daily evaporation rate in kgs. This figure roughly signifies the normal evaporation rate (N.E.R.) (first weight - second weight) x 1.235 = liters/day.

NOTE: During this time verify no major frost or condensation on outside of container, which would indicate either a weak or no vacuum.

9. Pour out LN₂. Take weight to determine 'charged weight'.
10. Calculate Hold Time. Subtract the Empty Weight from the Charged Weight and divide by 0.8083; then divide by NER. {[Charged Weight - Empty Weight] ÷ 0.8083} ÷ NER. Remember that the Empty Weight must be taken before filling, while the inside of the Vapor Shipper is at room temperature.

MVE vapor shippers were primarily designed as vapor shipping containers; however, they can also be used for immersion of samples. A sharp blow to the outer vessel can damage the neck tube or start a vacuum leak. Use caution and common sense in handling the container. Use the following weight table as a general guide to determine if your vapor shipper is fully charged.

Model	Empty Weight		Suggested Wt. Charged		Static Hold Time Days
	lbs	(kg)	lbs	(kg)	
SC 2/1V	6	(2.7)	8	(3.6)	8
SC 4/2V	10	(4.5)	17.1	(7.8)	13
SC 4/3V	11.6	(5.3)	19.1	(8.7)	21
Cryoshipper XC/IATA	33	(15.0)	49	(22.2)	14
Cryoshipper	26.1	(11.8)	39	(17.7)	10
Cryomoover	30.5	(13.8)	40.9	(18.6)	12
Mini-Moover	9	(4.1)	13.7	(6.2)	14
XC 20/3V	24.1	(10.9)	34	(15.4)	16
SC 20/12V	25.7	(11.7)	41.3	(18.7)	85
Doble Models					
Doble 11	14.5	(6.6)	18.5	(8.4)	17
Doble 20	22.7	(10.3)	29.7	(13.5)	21
Doble 28	33	(15.0)	47	(21.3)	24
Doble 34	34	(15.4)	50	(22.7)	20
Doble 47	40	(18.1)	53.5	(24.3)	21
Doble 47-10	40	(18.1)	53.5	(24.3)	21

MVE vapor shippers were primarily designed as vapor shipping containers; however, they can also be used for immersion of samples. A sharp blow to the outer vessel can damage the neck tube or start a vacuum leak. Use caution and common sense in handling the container. Use the following weight table as a general guide to determine if your vapor shipper is fully charged.

D) SHIPPING INSTRUCTIONS

Fill the dewar per Section 'B' and allow unit to stand for 2 hours (cooling down unit) with refills as needed.

Prior to packaging dewar for shipment the entire contents of liquid nitrogen must be removed. This is usually done by pouring out excess liquid nitrogen until no liquid is visible on the bottom of the inner dewar. Place in bio samples, package dewar and ship. The plastic shipping container provided by MVE/Chart is recommended to help keep dewar in upright position. **SHIPPING THE DEWAR IN ANY POSITION OTHER THAN UPRIGHT COULD REDUCE THE HOLDING TIME TO LESS THAN 10% OF NORMAL HOLD TIME. DO NOT SHIP UNIT ON ITS SIDE OR UPSIDE DOWN.** Remember that all MVE vapor shippers can also be used for liquid nitrogen storage as well, so it is imperative that all liquid nitrogen be removed so dewar remains classified as a vapor shipper. If liquid nitrogen is visible in the bottom of the inner it then becomes a liquid shipper and the exception status is void. The liquid inside is now classified as hazardous material.

E) UTILIZING THE QWICK CHARGE TECHNOLOGY

To utilize the QWick Charge Technology for faster charging times, do not allow inside of the dewar to warm up above -150°C between shipments. If the inside of the dewar does not warm up above -150°C, it will not need to be cooled down again before re-charging. Refill the dewar to the bottom of the neck with liquid nitrogen. Allow dewar to stand until liquid nitrogen stops boiling rapidly. After liquid nitrogen has stopped boiling rapidly, the dewar is fully charged. Pour off excess liquid nitrogen before shipment. The charged weight may be verified against the table in Section 'C' to ensure the dewar is fully charged. Please note that the Charged Weight is not an indication of unit performance.

Sanitizing and Decontaminating:

Any cleaning solution that does not react with aluminum or stainless steel can be used in the sanitation process of these dewars. In most cases, any household detergent or mild soap solution is suitable. Other cleaners and disinfectants that can be used safely include hydrogen peroxide, chlorine/water mixture, and denatured alcohol. The generally accepted practice of using 10% chlorine bleach with 90% water solution still holds as the best method for decontamination. **NOTE: DO NOT USE ANY PERTROLEUM BASED CLEANING SOLUTION.** It is important that all surfaces being sanitized are thoroughly rinsed and that all cleaner solution residues are removed after cleaning. For Vapor Shipper units, it is recommended that the unit is filled to its full capacity with the cleaning solution mixture, agitated, and then thoroughly rinsed. Allow the unit to dry completely before putting into service. It is suggested that the unit is inverted to drain and dry completely.

LETTER EXPLAINING EXCEPTION

This is concerning the applicability of the USA Federal Hazardous Material Regulations to the shipment of refrigerated samples in the "Dry Shipper" container. A "Dry Shipper" package consists of an outer container that is lined with an absorbent material. The container is charged with nitrogen refrigerated liquid which is absorbed into the container lining. The charged, completed package serves as refrigerated container for the shipment of samples.

In consideration of the above, consultation with the Research and Special Programs Administration of the DOT has determined that the use of nitrogen refrigerated liquid charged "dry shipper" containers for the shipment of samples fall within the regulation exception provided in 49CFR 173.320 paragraph (a) of the section states the requirements of this subchapter do not apply to atmospheric gases and helium when used in the operation of the process system such as refrigeration system. Paragraph (c) of 173.320 pertains to air transport of same refrigeration system. For exception status of air shipments please refer to IATA-Dangerous Goods Regulations for nitrogen refrigerated liquid. This falls in the class of 2.2 nonflammable gas, packing instructions 202 with special provisions A-152. If shipping outside of the USA, please consult your applicable regulatory agencies. For answers to questions regarding shipping regulations contact a Chart, AI-Cryobiological Tech Service Representative.

Limited Warranty Provisions:

Purchaser agrees that before this limited warranty shall become effective, Purchaser shall fully inspect each New Product, Repair/Blemished Product, or Service Part within ten (10) days of delivery and before such Product is put to use. Purchaser also agrees to operate the New Product, Repair/Blemished Product, or Service Part in accordance with CHART's operating instructions and that failure to do so shall void this limited warranty. Purchaser further agrees that any claim for breach of warranty must be made in writing within 60 days of discovery of a purported defect. CHART will not be responsible for any alleged breach of warranty, which, as a result of CHART's inspection, CHART determines to have arisen from a cause not covered by this limited warranty. In this case, CHART will charge the purchaser a nominal fee to repair the unit.

This limited warranty does not apply to: (A) Normal routine service items; (B) Repair or replacement necessitated by misuse, abuse, accident, or repairs made by persons other than CHART or persons not authorized by CHART, (C) Use of external equipment or parts with the New Product, Repair/Blemished Product, or Service Parts other than those approved by CHART, (D) Defects caused by effects of normal wear and tear; and (E) Acts of God, or other causes not within the control of CHART.

If Purchaser believes that a New Product, Repair/Blemished Product, or Service Part does not comply with the limited warranty stated above, Purchaser should contact CHART at the address stated above, describing the problem and providing proof of the date of purchase. If directed by CHART, Purchaser shall return the New Product, Repair/Blemished Product, or Service Part freight prepaid, properly packaged in a CHART approved shipping container and properly identified by a Return Material Authorization Number issued by CHART. New Products, Repair/Blemished Product, or Service Parts returned without a Return Material Authorization Number will be refused and returned at Purchaser's expense.

The remedies available for any breach of this limited warranty are limited to repair or replacement of the defective New Product, Repair/Blemished Product, Service Part, or refund of the purchase price, at the sole discretion of CHART. CHART warrants that replacement or repaired New Product, Repair/Blemished Product, or Service Part shall be free from defects in material and workmanship for the duration of the unexpired portion of the original warranty or ninety (90) days from the date of re-shipment to Purchaser, whichever is longer. PURCHASER'S RECOVERY FROM CHART FOR ANY CLAIM SHALL NOT EXCEED PURCHASER'S PURCHASE PRICE FOR THE NEW PRODUCT, REPAIR/BLEMISHED PRODUCT, OR SERVICE PART GIVING RISE TO SUCH CLAIM, IRRESPECTIVE OF THE NATURE OF THE CLAIM, WHETHER IN CONTRACT, TORT, WARRANTY, NEGLIGENCE, STRICT LIABILITY OR OTHERWISE. CHART SHALL NOT BE LIABLE FOR AND PURCHASER SHALL INDEMNIFY, DEFEND AND HOLD CHART HARMLESS FROM ANY CLAIMS BASED ON CHART'S COMPLIANCE WITH PURCHASER'S DESIGNS, SPECIFICATIONS OR INSTRUCTIONS, OR MODIFICATION OF ANY PRODUCTS BY PARTIES OTHER THAN CHART, OR USE IN COMBINATION WITH OTHER PRODUCTS.

PURCHASER SHALL NOT IN ANY EVENT BE ENTITLED TO, AND CHART SHALL NOT BE LIABLE FOR, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY NATURE INCLUDING, WITHOUT LIMITATION, BUSINESS INTERRUPTION COSTS, REMOVAL AND/OR REINSTALLATION COSTS, REPROCUREMENT COSTS, LOSS OF PROFIT OR REVENUE, LOSS OF DATA, PROMOTIONAL OR MANUFACTURING EXPENSES, OVERHEAD, INJURY TO REPUTATION OR LOSS OF CUSTOMERS, EVEN IF CHART HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

EXCEPT FOR THIS LIMITED WARRANTY, CHART HAS MADE NO WARRANTIES OR REPRESENTATIONS, EXPRESSED OR IMPLIED, AND HEREBY DISCLAIMS ANY OTHER WARRANTIES, INCLUDING, BUT NOT LIMITED TO, IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. NO REPRESENTATION OR STATEMENT OF CHART MAY CHANGE OR ALTER THIS LIMITED WARRANTY.

Any claims for breach of this limited warranty shall be governed by Georgia law and without regard to conflict of law rules thereof and must be brought in a state or federal court in Georgia.

Some states do not allow limitations on implied warranties or on incidental or consequential damages, so the above limitations may not apply. This limited warranty gives the Purchaser specific legal rights. Purchaser may also have other rights, which vary from state to state.

Limited Warranty Statement

Chart Inc., 2200 Airport Industrial Dr., Suite 500, Ball Ground, GA 30107

General Terms:

"New Product" – Chart Inc. ("CHART") warrants to the original purchaser ("Purchaser") that each new Vapor Shipper Dewar (collectively, the "New Products"), shall be free from defects in materials and workmanship for a period of two (2) years from the date of shipment, except as provided in the Limited Warranty Provisions below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship as for a period of three (3) years from the date of shipment, except as provided in the Limited Warranty Provisions below.

"Repair/Blemished Product" – CHART warrants to the Purchaser that all repaired equipment and factory blemished equipment (collectively, the "Repair/Blemished Products") shall be free of defects in materials and workmanship for a period of (90) days from the date of shipment, except as provided below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship for a period of (1) year from date of shipment, except as provided below.

"Service Part" – CHART warrants to the Purchaser that all Service Parts and Accessory Items (each a "Service Part") are warranted for (90) days from date of shipment to be free of defects in material and workmanship, except as provided below.

**Supplement to Chart Inc.
Limited Warranty Statement for Sales in Australia (AI)**

General Terms:

"New Product" – Chart Inc. ("CHART") warrants to the original purchaser only (meaning the party from whom CHART accepts the purchase order), whether such purchaser is a wholesaler, distributor or end user ("Purchaser"), that each new Liquid Nitrogen Dewar (collectively, the "New Products"), shall be free from defects in materials and workmanship for a period of two (2) years from the date of shipment, except as provided below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship as listed below, except as provided in the Limited Warranty Provisions below: Spectrum and CryoSystem Series, ET11, ET24, A-200, Workhorse 2003: 3 years / XC, SC, and Lab Series, ET 23-JLM: 5 years.

"Repair/Blemished Product" – CHART warrants to the Purchaser that all repaired equipment and factory blemished equipment (collectively, the "Repair/Blemished Products") shall be free of defects in materials and workmanship for a period of (90) days from the date of shipment, except as provided below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship for a period of (1) year from date of shipment, except as provided below.

"Service Part" - CHART warrants to the Purchaser that all Service Parts and Accessory Items (each a "Service Part") are warranted for (90) days from date of shipment to be free of defects in material and workmanship, except as provided below.

Sales of Product in Australia:

Where the Purchaser is the end user of the Product, the following additional information is provided in relation to the warranty against defects (the Warranty) provided by CHART (as set out in the Limited Warranty Statement).

1. In circumstances where a New Product, Repair/Blemished Product or Service Part does not comply with the Warranty, CHART will at its expense repair or replace that defective New Product, Repair/Blemished Product or Service Part, or refund the purchase price at the sole discretion of CHART. CHART warrants that the replacement or repaired New Product, Repair/Blemished Product or Service Part shall be free from defects in materials and workmanship for the duration of the unexpired portion of the Warranty, or 90 days from the date of reshipment to the purchaser, whichever is longer.
2. The Warranty is subject to the Limited Warranty Provisions set out below and all exclusions and limitations set forth in the CHART Limited Warranty Statement.
3. To claim the Warranty, the Purchaser must contact CHART's affiliate in Australia for further instruction, including the return of the defective New Product Repair/Blemished Product or Service Part (at the Purchaser's expense), to the following address:

Chart Australia Pty Ltd.
Unit 43/2 Railway Parade
Lidcombe, NSW 2141, Australia
Attention: Customer Service Department

Phone number: (02) 9749-4333
Email address: customerservice.australia@chartindustries.com

4. The Purchaser must otherwise comply with the requirements set out in the Limited Warranty Statement when making a claim under the Warranty.

5. Notwithstanding the Limited Warranty Statement, our Products come with guarantees that cannot be excluded under the Australian Consumer Law (ACL). You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the Products repaired or replaced if the Products fail to be of acceptable quality and the failure does not amount to a major failure.
6. The provider of the Warranty is CHART. Relevant contact details are set out below for CHART (to claim the Warranty, please use the contact details provided in paragraph 3 above):
 - (a) Customer Service Department
 - (b) Business Address:
Chart Inc.
2200 Airport Industrial Drive, Suite 500
Ball Ground, GA 30107 (USA)
 - (c) Phone number: 770-721-7759
 - (d) Email address: customerservice.usa@chart-ind.com

Repair of Products Sold in Australia

Where the Purchaser is the end user of the Product, the following additional information is provided in relation to the repair of Products by CHART:

The Products may contain "user generated data" as defined under the Australian Consumer Law (ACL). Repair of Products may result in the loss of user generated data.

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Medical Product Services, Borngasse 20, 35619 Braunfels, Germany

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IATA – International Air Transportation Association

ICAO – International Civil Aviation Organization

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Complément à la Déclaration de garantie limitée de Chart Inc. Limited pour les ventes en Australie (AI)