Guardian 5000  
Full Instruction manual

# Guardian 5000 Quick Start Guide

English

Español

Français

*Please download the user manual from [www.ohaus.com](http://www.ohaus.com).*

# Version History

Date	Version	Description
2024/7/25	A	<ul style="list-style-type: none"><li data-bbox="497 197 658 220">• New Release</li></ul>
2024/8/22	B	<ul style="list-style-type: none"><li data-bbox="497 244 927 266">• Added Run Dry Protection to the Menu</li><li data-bbox="497 272 938 295">• Updated the product specification sheet</li><li data-bbox="497 301 729 323">• Fixed content errors</li></ul>

# Contents

<b>1. Introduction</b>	<b>1</b>
1.1. Safety Information	1
1.2. Intended Use	2
<b>2. Installation</b>	<b>3</b>
2.1. Unpacking	3
2.2. Selecting the location	3
2.3. Connecting Power	3
2.4. Switch On / Off	3
<b>3. Product Overview</b>	<b>4</b>
3.1. Product Structure	4
3.2. Interface	5
3.3. Display	5
<b>4. Operations</b>	<b>7</b>
4.1. Basic Operations	7
<b>5. Menu</b>	<b>9</b>
5.1. Enter or Exit the Menu	9
5.2. Menu Map	9
5.3. Syetem - SYS	10
<b>6. Maintenance</b>	<b>12</b>
<b>7. Technical Data</b>	<b>14</b>
7.1. Conditions	14
7.2. Specifications	14
<b>8. Compliance</b>	<b>16</b>

# 1. Introduction

## 1.1. Safety Information

### Definition of Signal Warnings and Symbols

<b>WARNING</b>	For a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided.
<b>CAUTION</b>	For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or minor or medium injuries if not avoided.
<b>Attention</b>	For important information about the product. May lead to equipment damage if not avoided.
<b>Note</b>	For useful information about the product.

### Warning Symbols

Caution

Explosion  
HazardElectrical  
shockCaution, hot  
surfaceProtective  
Conductor  
TerminalAlternating  
Current

### Safety Precautions

**WARNING!** The protection provided by the unit may be impaired if used with accessories not provided or recommended by the manufacturer, or used in a manner not specified by the manufacturer.

- Always operate unit on a level surface for best performance and maximum safety.
- DO NOT lift unit by the top plate.
- Disconnect unit from the power supply prior to maintenance and servicing.
- Spills should be removed promptly, after the unit has cooled down.
- Alkalis spills, hydrofluoric acid or phosphoric acid spills may damage the unit and lead to thermal failure.
- DO NOT immerse the unit for cleaning.
- DO NOT operate the unit at high temperatures without a vessel/sample on the top plate.
- DO NOT operate the unit if it shows signs of electrical or mechanical damage.
- Protective earthing of the equipment is achieved via connection of the provided power cord to a compatible grounded power outlet.

**CAUTION!** To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplugging from the wall outlet.

WARNING! unit is not explosion proof. Use caution when unit is on or when heating volatile materials.

WARNING! DO NOT use the unit in explosive atmospheres or with materials that could cause a hazardous environment from processing. Keep in mind the material blink point relative to the target temperature that has been set.

CAUTION! The top plate can reach up-to 550°C, DO NOT touch the heated surface. Use caution at all times. Keep the unit away from explosive vapors and clear of papers, drapery and other flammable materials. Keep the power cord away from the heater plate.

CAUTION! Beware of the following risks when heating:

- flammable materials
- low boiling point combustible substances
- glass breakage as a result of mechanical shaking power
- incorrect container size
- too much medium
- unsafe condition of container

Earth Ground - Protective Conductor Terminal

Alternating Current

## 1.2. Intended Use

This instrument is intended for use in laboratories, pharmacies, schools, businesses and light industry. It must only be used for processing materials as described in these operating instructions. Any other type of use and operation beyond the limits of technical specifications, without written consent from OHAUS, is considered as not intended.

This instrument complies with current industry standards and the recognized safety regulations; however, it can constitute a hazard in use.

If the instrument is not used according to these operating instructions, the intended protection provided by the instrument may be impaired.

## 2. Installation

### 2.1. Unpacking

Check that there's no damage occurred during shipment. Notify the carrier for the damages found at the time of unpacking.

Packing list:

- Hotplate-Stirrer / Hotplate / Stirrer Unit
- Power Cord
- Stir Bar (Not applicable for e-G52HP07C)

### 2.2. Selecting the location

Select a location that meets the following requirement to position the instrument:

- the location must be sturdy, flat and level.
- away from explosive vapors
- the surface on which the instrument is placed can withstand the typical heat produced by the instrument
- avoid locations that is difficult to disconnect the power cord during use.

### 2.3. Connecting Power

The instrument is supplied with a 3-conductor grounded power cord. It should be plugged into a matching standard grounded outlet.

If the cord supplied does not meet your needs, please use an approved power cord that has ratings equal or exceeding those of the originally provided cord and that complies with the local/national regulations of the country in which the equipment is to be used.

Attention:

Replacement of the plug must be made by a qualified electrician.

### 2.4. Switch On / Off


After the power is connected, the instrument will be ready for switch on. Flip the **Standby Switch** to switch On or Off the instrument.

### 3. Product Overview

This chapter gives a product overview of its construction, interface, control panel and display.

#### 3.1. Product Structure

#	Controls	Functions
A	<b>Left Knob<sup>1</sup></b>	Control temperature and Menu (and Timer for e-G52HP07C)
B	<b>Right Knob<sup>2</sup></b>	Control speed and Timer (and Menu for e-G52ST07C)
C	<b>Display Screen</b>	Show operation status, parameters and menu settings
D	<b>Caution, Hot Surface Mark</b>	Hot surface warning
E	<b>Top Heating Plate</b>	Heat medium
F	<b>Accessory Rod Mounting Hole x 2</b>	Mount rod holders
G	<b>Standby Switch</b>	Switch ON / Off the instrument

**Note:**

- 1 Not applicable to e-G52ST07C
- 2 Not applicable to e-G52HP07C


### 3.2. Interface

#	Description
H	External Temperature Probe Port
I	RS232 Port
J	Power Entry Module

### 3.3. Display

Heating Display		General Display		Stirring Display	
A	Real-Time Temperature Display	R	Timer Display	L	Stirring Speed Display
B	Heating Indicator			M	Stirring Indicator
C	Set Temperature Display			N	Set Speed Icon
D	°C Icon			O	Speed Unit - RPM Icon
E	Set Temperature Icon			P	Soft Acceleration Rate Icon

Heating Display		General Display	Stirring Display	
F	Single Point Calibration Icon		Q	Fast Acceleration Rate Icon
G	External Probe Connection Icon			
H	Soft Heating Rate Icon			
I	Fast Heating Rate Icon			
J	Set Temperature Limit Icon (SmartHeat™)			
K	Hot Top Indicator Icon			

Attention:

The **Hot Top Indicator Icon** will illuminate when the top plate temperature is  $\geq 40^{\circ}\text{C}$ .

This icon will continuous to light up when the instrument is switched Off.

Attention: RUN DRY

The Run Dry alarm will be displayed and heating will be shut off, if a probe-controlled heating is running when probe tip not immersed in the heating medium, or sample level goes down due to evaporation.

Turn off the standby switch to clear the alarm.

 **Note:**

RUN DRY protection is by default Disable.

To enable the function, enter Menu > SYS > RUN DRY > EN

## 4. Operations

### 4.1. Basic Operations

#### Control Definition

Control	Definition
Short press	Press for > 0.05 seconds, release prior to 0.3 seconds.
Long press	Press for > 0.3 seconds, release prior to 1.5 seconds
Press & Hold	Press and hold for > 1.5 seconds
Rotate	Clockwise rotation to increment, counterclockwise rotation to decrement

#### Button Operations

Button	Task	Operation
Standby Switch	Turn on the instrument	Press the Switch
	Turn off the instrument	Press the Switch
	Adjust temperature	Rotate the Left Knob
	Start heating	Long press the Left Knob until you hear a beep sound or see the green light blinking.
	Adjust temperature when heating is on	Rotate the Left Knob, the Set Temperature is adjusted in real-time. (if Real-Time Adjustment is disabled, short press the knob to confirm the adjustment)
	Stop heating	Long press the Left Knob and release the knob when you hear a beep sound, or release in 1.5 second. When heating is stopped, the green light will fades.
	Adjust stirring speed	Rotate the Right Knob
	Start Stirring	Long press the Right Knob until you hear a beep sound or see the green light blinking.
	Adjust speed when stirring is on	Rotate the Right Knob, the Set Speed is adjusted in real-time. (if Real-Time Adjustment is disabled, short press the knob to confirm the adjustment)

<b>Button</b>	<b>Task</b>	<b>Operation</b>
	Stop stirring	Long press the Right Knob and release the knob when you hear a beep sound, or release in 1.5 second. When heating is stopped, the green light will fades.

## 5. Menu

### 5.1. Enter or Exit the Menu

#### Enter the Menu

**Note:**

Heating and Stirring need to be turned off to enter the Menu.

1. Press and hold the knob until Menu **Menu** is displayed.  
**Menu** would be displayed for 2 seconds, and then the first sub-menu will be displayed for configuration.
2. Rotate the knob to select sub-menu. Press the knob to enter the sub-menu.
3. Configure the settings
4. Press the knob to confirm the adjustment, and back to the Menu.

#### Exit the Menu

1. If in the sub-menu, rotate the knob to **BACK**, then press the knob to back to the Menu.
2. From the Menu, rotate the knob to Exit **EXIT**, then press the knob.

### 5.2. Menu Map

Menu	Sub-Menu
TEMP RATE	<ul style="list-style-type: none"> <li>• STD</li> <li>• FAST </li> <li>• SOFT </li> </ul>
TLIM	Max °C - 0 °C adjustable
SPC	<ul style="list-style-type: none"> <li>• NEW</li> <li>• CLR</li> </ul>
RPM RATE	<ul style="list-style-type: none"> <li>• STD</li> <li>• FAST </li> <li>• SOFT </li> </ul>

Menu	Sub-Menu
<b>SYS</b>	<ul style="list-style-type: none"> <li>• Beeper <b>BEEP</b></li> <li>• Timer Start Setting <b>TMDE</b></li> <li>• Power Recovery <b>PWRR</b></li> <li>• Reset <b>RSET</b></li> <li>• Real-time Adjustment <b>RTA</b></li> <li>• Run Dry Protection <b>RUN DRY</b></li> <li>• System Version <b>V 1.0 1<sup>1</sup></b></li> </ul> <div style="background-color: #e0f2f1; padding: 10px; margin-top: 10px;"> <p> <b>Note:</b> 1 Example display of system version</p> </div>

### 5.3. System - SYS

SYS Sub-menu	Definition	Options
<b>BEEP</b>	Set to enable or disable beeper.	<ul style="list-style-type: none"> <li>• Enable <b>EN</b> (default)</li> <li>• Disable <b>DIS</b></li> </ul>
<b>TMDE</b>	Timer Start Setting controls the initiation of the timer in count or countdown situations	<ul style="list-style-type: none"> <li>• <b>STD</b> <b>STD</b> (default): timer start when heating or stirring is on</li> <li>• Temperature dependent <b>TEMP</b>: the timer will start once the heater temperature or probe temperature reaches the set temperature in count or countdown situations</li> </ul>
<b>PWRR</b>	If power has been disconnected, Power Recovery allows the unit to automatically start heater and stirrer functions again when power is returned.	<ul style="list-style-type: none"> <li>• Enable <b>ON</b></li> <li>• Disable <b>OFF</b> (default)</li> </ul>
<b>RSET</b>	Restore the instrument to factory settings.	<ul style="list-style-type: none"> <li>• <b>YES</b> (long press the knob to confirm selection)</li> <li>• <b>NO</b></li> </ul>
<b>RTA</b>	Enable RTA (Real-time Adjustment), temperature and stirring speed will be adjusted in real-time. Disable RTA, the classic Guardian adjustment mode will be activated automatically.	<ul style="list-style-type: none"> <li>• Enable <b>ENEN</b> (default)</li> <li>• Disable <b>DIS</b></li> </ul>

<b>SYS Sub-menu</b>	<b>Definition</b>	<b>Options</b>
<b>RUN DRY</b>	Enable or disable RUN DRY protection.	<ul style="list-style-type: none"><li>• Enable <b>EN</b></li><li>• Disable <b>DIS</b> (default)</li></ul>
<b>V 1.0 1</b>	Display of the current system version	--

**Note:**

1 Example display of system version

## 6. Maintenance

### Cleaning

**WARNING:** Electric Shock Hazard. Disconnect the equipment from the power supply before cleaning. Make sure that no liquid enters the interior of the instrument.

**Attention:** Do not use solvents, harsh chemicals, ammonia or abrasive cleaning agents.

The housing may be cleaned with a cloth dampened with a mild detergent if necessary.

### Error Codes

Error Code	Description / Cause of Error	How to Fix
<b>N/A</b>	Unit fails to power on Missing or blown fuse	Not fixable by user, please contact OHAUS.
<b>E1</b>	Plate RTD Disconnected	Not fixable by user, please contact OHAUS.
<b>E2</b>	Plate RTD Error	Not fixable by user, please contact OHAUS.
<b>E3</b>	Stir Error	Flip the standby switch off and back on. If the problem persists, please contact OHAUS.
<b>E4</b>	Probe RTD Disconnected	Flip the standby switch off and back on. If the problem persists, please contact OHAUS.
<b>E5</b>	Probe RTD Error	<ul style="list-style-type: none"> <li>Flip the standby switch off, and remove the probe from the unit. Then, flip the standby switch on to return the normal operating mode.</li> <li>If this does not solve the issue, replace the external probe or the main PCBA, and then try again.</li> <li>If the problem persists, please contact OHAUS.</li> </ul>
<b>E7</b>	User probe error. An External Probe is connected to the unit while heating is on.	Flip the standby switch off and back on. If the problem persists, please contact OHAUS.
<b>E8</b>	Plate over temperature	Not fixable by user, please contact OHAUS.
<b>E9</b>	Heating not working	Not fixable by user, please contact OHAUS.
<b>E10</b>	Triac Fault	Not fixable by user, please contact OHAUS.
<b>AC Err</b>	Mains frequency out of the range 40~55Hz (rating 50Hz) or 55~70Hz (rating 60Hz)	Regulate mains frequency within the range.

## **Technical Support Information**

For technical issues, please speak to an Authorized Ohaus Service Agent. Please visit our website [www.ohaus.com](http://www.ohaus.com) to find the Ohaus office nearest you.

## 7. Technical Data

### 7.1. Conditions

The technical data is valid under the following conditions:

Indoor Use Only

Altitude: 0 to 2000m

Operating temperature: 5 to 40°C

Storage temperature: -20 to 65°C

Operating humidity: 20 to 80% relative humidity, non-condensing

Storage Humidity: 20 to 80% relative humidity, non-condensing

Electrical Supply: 100 - 120V ~, 1-15A, or 220 - 240V ~, 1-10A. (Depending on the models)

Voltage fluctuations: Mains supply voltage fluctuations up to  $\pm 10\%$  of the nominal voltage.

Overvoltage category (Installation category): II

Pollution Degree: 2

### 7.2. Specifications

Guardian 5000-2nd Generation					
Model No.	e-G52HSRDA	e-G52HS10C	e-G52HS07C	e-G52HP07C	e-G52ST07C
Function	Heating and Stirring		Heating		Stirring
Top Plate Size	$\varnothing$ 5.3 in / 13.5 cm	10 x 10 in / 26 x 26 cm	7 x 7 in / 18 x 18 cm		
Top Plate Material	Ceramic Coated Aluminum	Ceramic			
Control	Digital				
Display	115 x 30 mm White Character LCD				
Timer	1 minute to 99 hours 59 minutes				
Temperature Range	Ambient + 5°C – 360°C	Ambient + 5°C – 500°C	Ambient + 5°C – 550°C		--
Set Temperature Increments	0.5°C				--
Temperature Readability	0.5°C				--
Temperature Stability <sup>1</sup>	+/-1% Plate, below 100°C +/-1°C +/-0.5% Probe, below 100°C +/-0.5°C				--
Temperature Calibration (SPC) Set Points	3 points for Plate 3 points for Probe				--
Probe Connection	PT100, PT1000 (Class A)				--

Guardian 5000-2nd Generation					
Model No.	e-G52HSRDA	e-G52HS10C	e-G52HS07C	e-G52HP07C	e-G52ST07C
Probe Temperature Measurement Accuracy	$\pm 0.2^{\circ}\text{C} + \text{PT1000 (Class A) Tolerance}^2$				--
Heating Power	<ul style="list-style-type: none"> <li>600W at 120V</li> <li>656W at 230V (600W at 220V)</li> </ul>	<ul style="list-style-type: none"> <li>1650W at 120V</li> <li>2066W at 230V (2136W at 240V)</li> </ul>	<ul style="list-style-type: none"> <li>1582W at 120V (1050W at 100V)</li> <li>1050W at 230V</li> </ul>	--	
SmartHeat™	User is able to adjust Set Temperature upper limit of the Top Plate, from 0°C to Rated Max Temperature				--
Speed Range	50 - 1800 rpm			--	50 - 1800 rpm
Speed Stability <sup>1</sup>	+/-2%			--	+/-2%
Stirring Capacity	20 L H2O	22 L H2O	20 L H2O	--	20 L H2O
SmartRate™	User is able to select the speed or temperature ramping rate				
Data Logging & Remote Control	built-in RS232 optional RS232-USB interface kit optional RS232-Ethernet interface kit				
IP Rating	IP32				
Housing	Painted Aluminum				
Dimensions (L x W x H)	259 x 175 x 108 mm 10.21 x 6.87 x 4.26 inch	402 x 287 x 118 mm 15.84 x 11.28 x 4.65 inch	329 x 221 x 116 mm 12.94 x 8.69 x 4.56 inch		
Working Environment	41°F – 104°F, 80%RH, non-condensing / 5°C – 40°C, 80%RH, non-condensing				
Net Weight	5.7 lb / 2.6 kg	14.8 lb / 6.7 kg	10.1 lb / 4.6 kg	8.6 lb / 3.9 kg	9.9 lb / 4.5 kg
Electrical (50/60Hz)	100-120V, 8A 220-240V, 4A	100-120V, 15A 220-240V, 10A	100-120V, 15A 220-240V, 7A		100-120V, 1A 220-240V, 1A


**Note:**

<sup>1</sup> Conditions permitting. Variations in temperature and speed measurement processes, vessel, ambient and sample will impact actual performance. To improve temperature accuracy of the system use the Single Point Calibration feature.

<sup>2</sup> PT1000 Class A Tolerance:  $\pm(0.15^{\circ}\text{C} + 0.002 \times T^{\circ}\text{C})$ , T is Set Temperature

## 8. Compliance

Compliance to the following standards is indicated by the corresponding mark on the product.

Mark	Standard
	This product complies with the applicable harmonized standards of EU Directives 2011/65/EU (RoHS), 2014/30/EU (EMC) and 2014/35/EU (LVD) The EU Declaration of Conformity is available online at <a href="http://www.ohaus.com/ce">www.ohaus.com/ce</a> .
	This product complies with the applicable statutory standards of the Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012, UK Electromagnetic Compatibility Regulations 2016 and Electrical Equipment (Safety) Regulations 2016. The UK Declaration of Conformity is available online at <a href="http://www.ohaus.com/uk-declarations">www.ohaus.com/uk-declarations</a> .
	This product complies with the EU Directive 2012/19/EU (WEEE). Please dispose this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to <a href="http://www.ohaus.com/weee">www.ohaus.com/weee</a> .
	EN 61326-1
	CAN/CSA-C22.2 No. 61010-1, CAN/CSA-C22.2 No. 61010-2-010, CAN/CSA-C22.2 No. 61010-2-051
	UL 61010-1, UL 61010-2-010, UL 61010-2-051

### ISED Canada Compliance Statement:

CAN ICES-003(A) / NMB-003(A)

### ISO 9001 Registration

The management system governing the production of this product is ISO 9001 certified.

## 8.1. FCC Supplier Declaration of Conformity

### Unintentional Radiator per 47CFR Part B

Trade Name: OHAUS CORPORATION

Model: Guardian 5000, e-52 series

### Party issuing Supplier's Declaration of Conformity:

Ohaus Instruments (Changzhou) Co., Ltd.

C Block, 6 Zhengqiang Road, Xinbei District, Changzhou

Jiangsu 213022,

China

Phone: +86 519 85287270

**Responsible Party – U.S. Contact Information:**

Ohaus Corporation  
8 Campus Drive, Suite 105  
Parsippany, NJ 07054  
United States  
Phone: +1 973 377 9000  
Web: [www.ohaus.com](http://www.ohaus.com)

**FCC Compliance Statement:**

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications.

Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

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P/N 30910711 B © 2024 Ohaus Corporation, all rights reserved. / todos los derechos reservados./ tous droits réservés.