**STARTER 400D
Portable DO Meter
Instruction Manual**

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TABLE OF CONTENTS

1	INTRODUCTION	1
1.1	Definition of Signal Warnings and Symbols	1
1.2	Safety Precautions	2
1.3	Display and controls	3
2	INSTALLATION	5
2.1	Package contents	5
2.2	Installing the batteries	5
2.3	Installing the electrode clip	6
2.4	Installing the IP54 cover caps	6
2.5	Integrated stand for table top use	6
3	SETUP	8
3.1	Setting Auto/Manual Endpoint	7
3.2	Setting Salinity Value	7
3.3	Setting Temperature Unit	8
3.4	Setting the Date and Time	8
3.5	Setting Power State	8
3.6	Setting Factory Reset	8
3.7	Setting New Probe Cap Parameter	9
4	OPERATION	10
4.1	Calibration	10
4.1.1	Zero-oxygen solution	10
4.1.2	Performing 1-point calibration	10
4.1.3	Performing 2-point calibration	11
4.2	Sample measurement	11
4.3	Temperature measurement	11
4.4	Barometer pressure measurement	11
4.5	Using the memory	12
4.5.1	Storing a reading	12
4.5.2	Recalling from memory	12
4.5.3	Clearing the memory	12
5	MAINTENANCE	13
5.1	Error message	13
5.2	Meter maintenance	13
5.3	Probe maintenance	14
6	TECHNICAL DATA	15
6.1	Specifications	15

1 INTRODUCTION

😊 Thank you for choosing OHAUS. Please read the manual completely before using the **STARTER 400D** portable Optical DO (**D**issolved **O**xygen) meter to avoid incorrect operation.

The STARTER 400D has an excellent performance/price ratio and is designed with many useful features. Optical DO probe Other accessories are also available from OHAUS.

Starter 400D offers many practical features such as:

- Ergonomic, multifunctional design for ambidextrous measurements or bench top use
- Meter can measure barometric pressure automatically to compensate the DO value; Salinity value can be set to compensate the DO value.
- IP54 protection against dust and water
- Ohaus' renowned user friendly software

1.1 Definition of Signal Warnings and Symbols

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results.

Signal Words

WARNING	For a hazardous situation with medium risk, possibly resulting in injuries or death if not avoided.
CAUTION	For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or injuries if not avoided.
Attention Note	For important information about the product For useful information about the product

Warning Symbols

General Hazard

Explosion Hazard

Corrosive Hazard

1.2 Safety Precautions

CAUTION Read all safety warnings before installing, making connections, or servicing this equipment. Failure to comply with these warnings could result in personal injury and/or property damage. Retain all instructions for future reference.

- Use the equipment only in dry locations.
- Dry off any liquid spills immediately. The instrument is not watertight.
- Use only approved accessories and peripherals.
- Operate the equipment only under ambient conditions specified in these instructions.
- Disconnect the equipment from the power supply when cleaning.
- Do not operate the equipment in hazardous or unstable environments.
- Service should only be performed by authorized personnel.


WARNING Do not use this instrument in an environment subject to explosion hazards. The housing of the instrument is not gas tight (explosion hazard due to spark formation, corrosion caused by the ingress of gases).

CAUTION When using chemicals and solvents, comply with the instructions of the producer and the general lab safety rules.

1.3 Display and controls

Display

31	Time: HH:MM:SS	
2	Date: YYYY-MM-DD	
3	 calibration data: $ K-1 \leq 0.3$ and $ B < 0.1$ Probe condition is very good	 calibration data: $0.3 \leq K-1 \leq 0.4$ and $ B < 0.1$ Probe condition is acceptable.
4	Measurement icon - measurement (calibration) is running	
5	Calibration icon - calibration is running	
6	Power state icon, include power always on icon , means meter will not automatically turn off; and auto-off icon which means after 10 minutes no operation meter will turn off(see 4.4); power state icon position is at the top-left of the battery status icon.	
7	Battery status icon – the charged status of the battery	
8	DO reading: mg/L and %	
9	Temperature during measurement or in calibration process	
10	Barometric pressure during measurement or in calibration process	
11	Data log: memory number	
12	ppt for salinity entry	

Controls

Button	Press & release	Press & hold for 2 seconds 
	<ul style="list-style-type: none"> - Start or endpoint measurement (manual endpoint mode) - Confirm setting, confirm entered value for setting 	<ul style="list-style-type: none"> - switch between auto endpoint and manual endpoint mode
	<ul style="list-style-type: none"> - Start calibration 	<ul style="list-style-type: none"> - Review the latest calibration data
	<ul style="list-style-type: none"> - Meter turn on - Back to measurement screen 	<ul style="list-style-type: none"> - Meter turn off
	<ul style="list-style-type: none"> - Store current reading to memory - Increase value during setting - Scroll up through the memory 	<ul style="list-style-type: none"> - Recall stored data
	<ul style="list-style-type: none"> - Enter setup screen - Scroll down through the memory 	
		<ul style="list-style-type: none"> - Reset factory setting

2 INSTALLATION

Carefully unpack the box.

2.1 Package contents

The ST400D kit has the following items:

ST400D	Units
STARTER 400D meter	1
STDO21 (with 1.0m cable)	1
AAA battery	4
Electrode Clip	1
IP54 cover caps	1 set
Wrist Strap	1


We also offer ST400D-G kit which includes ST400D kit and a portable bag.

Accessories:

Model	Description	P/N
STDO21	Optical DO probe	30219117


DO Zero Oxygen Chemical	30059257
Portable Bag for portable meters	30031635

2.2 Installing the batteries

- Remove the battery cover screw using a coin or screwdriver and remove the battery cover
- Insert the batteries in the battery compartment as shown
- Replace the battery cover and tighten the battery cover screw

2.3 Installing the electrode clip


The electrode clip is an electrode holder that can be placed next to the display on either side of the housing.

- Attach the clip by inserting the tabs into the recess.
- Slide the clip forward to lock in position. Slide the shaft of the probe into the clip from above.


2.4 Installing the IP54 cover caps

Use the small plastic tool as shown:

2.5 Integrated stand for table top use

Starter 400D has an integrated table stand as shown, which can be pushed out to use the meter on a lab bench or desk. To close the stand, push the middle edge of the stand toward the meter.

3 SETUP

Press  to enter the setup interface. We have 7 setting items, they are:


- 1) Auto/Manual Endpoint
- 2) Salinity Value
- 3) Temperature Unit
- 4) Date and Time
- 5) Power State
- 6) Reset Factory Setting
- 7) New Probe Cap Parameter

3.1 Setting Auto/Manual Endpoint

When in calibration, ST400D is only work at Auto End Point Mode.

The default endpoint mode of ST400D for measurement is Auto End Point Mode. When users want to judge the endpoint manually, you can shift to Manual Endpoint Mode.

End point Mode:

- When in Manual Mode, to manually reach and lock a pH measurement value, you need to press button- when reading is stable, then the measurement value freezes; And meter displays .
- When in Auto End Point Mode, the meter determines when the reading is stable then locks the measurement or calibration value automatically, the reading freezes; the meter displays .

Users can also change the endpoint mode by shortcut -- Press and hold button- to shift the **End Point Mode**.

3.2 Setting Salinity Value

Since salinity has an influence on the dissolved oxygen value, it is important to input the correct salinity value. The meter can compensate salinity for measurements in solutions with salinity from 0 to 50 ppt. The default value for salinity is 0.0 ppt.

- Press button- or button- to adjust the salinity value.
- Press button- to confirm the value or press button- to return to the Setting screen.

Normally seawater salinity is 35ppt, which means 35g of salt per 1000g seawater.

3.3 Setting Temperature Unit

Ohaus optical DO probe STDO21 include a built-in temperature probe, the meter will be in Automatic Temperature Compensation (ATC) mode automatically. You can select the temperature unit according to your preferred: °C or °F.

Note: °C = 5/9 (°F -32)

The default setting is 25 °C (77°F).

3.4 Setting the Date and Time

The Date format is yy-mm-dd. Press button- or button- to adjust the value of each numerical digit. Press button- to confirm current value then goes to next numerical digit.

3.5 Setting Power State

We have 2 power states: "Auto OFF" or "Always ON".

When "Auto OFF" is selected, which means the meter will turn off the powder without any press operation after 10 minutes and auto-off icon  will display at the top-left of the battery status icon.

"Always ON" means power will not turn off automatically, and the power always-on icon  displays.


Ohaus suggest you to select "Auto OFF" so to save power of batteries.

3.6 Setting Factory Reset


Press button- and button- together to enter the Reset Factory Setting screen.

3.7 Setting New Probe Cap Parameter

When replace a new probe cap, user must input the parameter of this new cap into the meter. Or the measurement value will be incorrect.

Press button-, button- and button- together to enter the new probe cap parameter setting screen. Input the parameter of the new cap. Then do calibration before any measurement.

4 OPERATION

WARNING Do not operate the equipment in hazardous environments. The equipment is not explosion protected.

WARNING When using chemicals and solvents, comply with the instructions of the chemical producer and the general lab safety rules.

4.1 Calibration

STARTER 400D allows you to perform 1-point or 2-point calibrations; Ohaus suggest 1-point calibration(100% DO calibration) for samples which DO value is not less than 2 mg/L.

Only when measure samples which DO value is lower than 2 mg/L, Ohaus suggest 2-point calibration.

4.1.1 Zero-oxygen solution

Preparing a zero-oxygen solution: put approximately 20ml of pure or deionized water into the bottle of DO Zero Oxygen Chemical (30059257, there are about 5g powder in it), stir for at least 5 minutes.

This zero-oxygen solution is used for 2-point calibration.

4.1.2 Performing 1-point calibration

When performing calibration, the meter only uses Auto End Point Mode.

1 point calibration can be performed in the air or oxygen-saturated water. There will be a little difference between calibration in the air and in oxygen-saturated water.

Before calibration you need to rinse the DO probe with pure water and carefully wipe the water off with a tissue.

- Place the probe in the air (or in oxygen-saturated water).
- Press button  , the calibration icon  appear on the display.
- The meter reaches endpoints automatically, 100% will display. Measurement icon  disappears.

Once you have completed your 1-point calibration, you have 3 options new:

- Press button-  to finish 1-point calibration then return to measurement screen.
- Press button-  to reject calibration and return to measurement screen.
- Press button-  to complete a 2-point calibration (after preparing the zero oxygen solution and putting STDO21 into it).

4.1.3 Performing 2-point calibration

- Perform the 1-point calibration.
- Place the probe in the zero oxygen solution and press button-  .

The meter reaches the endpoint automatically. 0% will display then return to measurement screen.

Pressing and holding button-  will display the calibration data.

4.2 Sample measurement

Dissolved oxygen readings can display unit of mg/L or %.

Place the DO probe in the sample. Press  to start the measurement.

When the meter reaches endpoint, the measurement reading will freeze on the screen.

4.3 Temperature measurement

For better accuracy, we recommend the use of the temperature probe. If a temperature probe is used, **ATC** and the sample temperature are displayed.

STDO21 has a built-in temperature probe.

4.4 Barometer pressure measurement

Barometric pressure influences the DO value.

STARTER 400D has a built-in barometric pressure sensor which can measure the barometric pressure automatically. This meter uses barometric pressure unit **KPa** (kilopascal).

4.5 Using the memory

4.5.1 Storing a reading

The STARTER 400D can store up to 99 measurement readings. Press button-  when the measurement reaches end point. ^{Data}Log 02 indicates that is the 2nd measurement reading had been stored. The starter 400D can store up to 99 measurement readings.

4.5.2 Recalling from memory

Press and hold button-  to recall the stored values from memory.

Press button-  or button-  to scroll through the stored results. Press button-  to go back to the measurement screen.

4.5.3 Clearing the memory

Continue pressing button-  to scroll through the stored results until "Clear All" appears.

Press button-  to confirm the deletion or press button-  to return to the measurement mode without deleting all the data in the memory.

5 MAINTENANCE

5.1 Error message

Low battery	If low battery, the battery icon will blink to show need new battery. Then turn off after 10 seconds.
Communication fail	When do calibration or new cap parameter input, meter will detect if the probe is connected properly. If not, meter displays " Fail, please confirm the STDO21 probe is connected"
Measured temperature out of range (< -5°C or >55 °C)	Meter will displays "Temperature Out of Range!" Please keep the sample temperature within the range for measurement.
Calibration out of range	Meter displays "Calibration Fail". $ K-1 > 0.4$ or $ B > 0.1$ You may need to do a proper calibration again, if calibration data is still not good, you may need to replace a new probe cap or probe.
Recall memory fail	Meter displays "No Data" When recall memory while no data memorized yet.
The current data set has already been stored once	Meter displays "Already saved!" A measurement can only be stored once. Perform a new measurement to store a new data set.

Any more problems, you can contact Ohaus.

5.2 Meter maintenance

The STARTER 400D does not require any maintenance other than replacement of depleted batteries. To clean, use a damp cloth.

Attention: Do not disassemble the instrument. There are no user serviceable parts inside.

The housing is made of acrylonitrile butadiene styrene (ABS). This material is susceptible to damage by some organic solvents, such as toluene, xylene and methyl ethyl ketone (MEK). Any spillage should be immediately wiped off.

5.3 Probe maintenance

The DO probe measures the oxygen based on luminescence principle. The optical DO probe measures dissolved oxygen by emitting a blue light of the proper wavelength that causes the sensing element to luminesce. Oxygen dissolved in the sample continually passes through the diffusion layer to the sensing layer, affecting the luminescence lifetime. STDO21 can measure the lifetime of the luminescence as it is affected by the presence of oxygen and compared with the reference. This technology has been approved by EPA.

When measurement, optical probe did not consume the DO in the sample, that means the sample do not need stir

Optical probe STDO21 do not need to replace the membrane or any maintenance of anode or cathode. The optical probe only requires that its sensor cap be replaced every 12-18 months. Please input the new cap parameter when new probe cap in installed, then do calibration before any measurement.

When storing the probe for a longer period, make sure it is placed in the wetting protection cover.

Calibration frequency: Ohaus suggest verifying the calibration regularly according to your measurement requirement. STDO21 only has very little calibration drift when compared to electrochemical probes. This advantage results in less frequent calibrations. In fact, optical probes are so stable they are capable of holding their calibration for many months.

Please also read the STDO21 probe manuals before using it.

6 TECHNICAL DATA

6.1 Specifications

Ambient conditions

- Altitude: Up to 2000 m
- Specified Temperature range: 5°C to 40°C
- Humidity: maximum relative humidity 80 % for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C
- Installation category: N/A
- Pollution degree: N/A
- Operability is assured at ambient temperatures between 5°C to 40°C

Model: ST400D

Measurement Range	%	0.0 to 200.0%
	mg/L	0.00 to 20.00mg/L (ppm)
	Temp.	0 to 50 °C
Measurement Resolution	%	0.1%
	mg/L	0.01 mg/L
	Temp.	0.1 °C
Accuracy		± 0.2 mg/L (<8 mg/L) ± 0.3 mg/L (8 to 20 mg/L); ± 0.3 °C
Barometric	Range	50.0 to 115.0 kPa
	Resolution	0.1 kPa
	Accuracy	1.5 kPa
Temp compensation		ATC
Salinity compensation		0.0 to 40.0 ppt
Calibration		1 point or 2 points
Memory		99 sets, last calibration data
Power		4 x AAA(LR03) batteries, 12 hrs
Size/Weight		90 W x 150 D x 35 H mm
Weight		0.16kg
Display		liquid crystal
IP rating		IP54
T input		Mini-Din
Materials		ABS

Compliance

This product conforms to the EMC Directive 2004/108/EC and the Low Voltage Directive 2006/95/EC. The Declaration of Conformity is available online at www.ohaus.com/ce

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements. Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related. Disposal instructions in Europe are available online at www.ohaus.com/weee. Thank you for your contribution to environmental protection.

FCC Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

ISO 9001 Registration

In 1994, OHAUS Corporation, USA, was awarded a certificate of registration to ISO 9001 by Bureau Veritas Quality International (BVQI), confirming that the OHAUS quality management system is compliant with the ISO 9001 standards requirements. On June 21, 2012, OHAUS Corporation, USA, was re-registered to the ISO 9001:2008 standard.

LIMITED WARRANTY

Ohaus products are warranted against defects in materials and workmanship from the date of delivery through the duration of the warranty period. During the warranty period Ohaus will repair, or, at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to Ohaus.

This warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials, has foreign material penetrating to the inside of the product, or as a result of service or modification by other than Ohaus. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by Ohaus Corporation. Ohaus Corporation shall not be liable for any consequential damages.

As warranty legislation differs from state to state and country to country, please contact Ohaus or your local Ohaus dealer for further details

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PN 30208731 B © Ohaus Corporation 2015, all rights reserved/ todos los