




Accessories

Sliding and Fixed Baskets



BSI-SLD BSKT-FULL-5.2



Compatible Refrigerator Size	5.2 ft ³
Max Quantity	3
Basket type	Sliding
Width x Depth x Height	16 7/8" x 13 5/8" x 3"
Max Load	30 lbs.


*Factory install only

BSI-BSKT-BASE-FIXED-5.2



Compatible Refrigerator Size	5.2 ft ³
Max Quantity	1
Basket type	Fixed
Width x Depth x Height	18 3/8" x 8 1/2" x 6 1/4"
Max Load	44 lbs.



BSI-SLD BSKT-FULL-10.5



Compatible Refrigerator Size	10.5 ft ³
Max Quantity	6
Basket type	Sliding
Width x Depth x Height	16 7/8" x 15 3/4" x 3"
Max Load	35 lbs.

*Factory install only

BSI-BSKT-BASE-FIXED-10.5



Compatible Refrigerator Size	10.5 ft ³
Max Quantity	1
Basket type	Fixed
Width x Depth x Height	18 3/8" x 10 5/8" x 6 1/4"
Max Load	55 lbs.



Accessories

Sliding and Fixed Baskets

BSI-SLD BSKT-FULL-5.2




BSI-BSKT-BASE-FIXED-5.2



BSI-SLD BSKT-FULL-10.5



BSI-BSKT-BASE-FIXED-10.5



BSI-HC-UCFS-0504
 BSI-HC-UCFS-0504G
 †BSI-HC-RFC7
 †BSI-HC-RFC7A
 †BSI-HC-RFC7S
 †BSI-HC-RFC7SA
 †BSI-HC-RFC1020
 †BSI-HC-RFC1020G
 †BSI-HC-RFC1030
 †BSI-HC-RFC1030G
 †BSI-HC-RFC1040
 †BSI-HC-RFC1040G
 CRT-BSI-HC-UCFS-0504G
 CRT-BSI-HC-UCFS-0504
 PH-BSI-NSF-UCFS-0504
 PH-BSI-NSF-UCFS-0504G
 PH-BSI-HC-UCFS-0504G
 PH-BSI-HC-UCFS-0504
 †PH-BSI-HC-RFC7
 †PH-BSI-HC-RFC7A
 †PH-BSI-HC-RFC7S
 †PH-BSI-HC-RFC7SA
 †PH-BSI-HC-RFC7A-CAD
 †PH-BSI-HC-RFC7SA-CAD
 †PH-BSI-HC-RFC1020
 †PH-BSI-HC-RFC1020G
 †PH-BSI-HC-RFC1030
 †PH-BSI-HC-RFC1030G
 †PH-BSI-HC-RFC1040
 †PH-BSI-HC-RFC1040G

BSI-HC-UCBI-0404*
 BSI-HC-UCBI-0404G*
 BSI-HC-UCBI-0404SS*
 BSI-HC-UCBI-0404-ADA*
 BSI-HC-UCBI-0404G-ADA*
 BSI-HC-UCFS-0504
 BSI-HC-UCFS-0504G
 †BSI-HC-RFC7
 †BSI-HC-RFC7A
 †BSI-HC-RFC7S
 †BSI-HC-RFC7SA
 †BSI-HC-RFC9*
 †BSI-HC-RFC9G*
 †BSI-HC-RFC9SS*
 †BSI-HC-RFC1020
 †BSI-HC-RFC1020G
 †BSI-HC-RFC1030
 †BSI-HC-RFC1030G
 †BSI-HC-RFC1040
 †BSI-HC-RFC1040G
 CRT-BSI-HC-UCBI-0404*
 CRT-BSI-HC-UCBI-0404G*
 CRT-BSI-HC-UCFS-0504G
 CRT-BSI-HC-UCFS-0504
 PH-BSI-NSF-UCBI-0404*
 PH-BSI-NSF-UCBI-0404G*
 PH-BSI-NSF-UCBI-0404SS*
 PH-BSI-NSF-UCBI-0404-ADA*
 PH-BSI-NSF-UCBI-0404G-ADA*
 PH-BSI-NSF-UCFS-0504
 PH-BSI-NSF-UCFS-0504G
 PH-BSI-HC-UCBI-0404*
 PH-BSI-HC-UCBI-0404G*
 PH-BSI-HC-UCBI-0404SS*
 PH-BSI-HC-UCBI-0404-ADA*
 PH-BSI-HC-UCBI-0404G-ADA*
 PH-BSI-HC-UCFS-0504G
 PH-BSI-HC-UCFS-0504
 †PH-BSI-HC-RFC7
 †PH-BSI-HC-RFC7A
 †PH-BSI-HC-RFC7S
 †PH-BSI-HC-RFC7SA
 †PH-BSI-HC-RFC7A-CAD
 †PH-BSI-HC-RFC7SA-CAD
 †PH-BSI-HC-RFC9*
 †PH-BSI-HC-RFC9G*
 †PH-BSI-HC-RFC9SS*
 †PH-BSI-HC-RFC1020
 †PH-BSI-HC-RFC1020G
 †PH-BSI-HC-RFC1030
 †PH-BSI-HC-RFC1030G
 †PH-BSI-HC-RFC1040
 †PH-BSI-HC-RFC1040G

BSI-HC-10PG
 BSI-HC-10PS
 †BSI-HC-RFC12G
 †BSI-HC-RFC12
 †BSI-HC-RFC12A
 †BSI-HC-RFC12GA
 CRT-BSI-HC-10PG
 CRT-BSI-HC-10PS
 PH-BSI-NSF-10PG
 PH-BSI-NSF-10PS
 PH-BSI-HC-10PG
 PH-BSI-HC-10PS
 †PH-BSI-HC-RFC12G
 †PH-BSI-HC-RFC12
 †PH-BSI-HC-RFC12A
 †PH-BSI-HC-RFC12A-CAD
 †PH-BSI-HC-RFC12GA
 †PH-BSI-HC-RFC12GA-CAD

BSI-HC-10PG
 BSI-HC-10PS
 †BSI-HC-RFC12G
 †BSI-HC-RFC12
 †BSI-HC-RFC12A
 †BSI-HC-RFC12GA
 CRT-BSI-HC-10PG
 CRT-BSI-HC-10PS
 PH-BSI-NSF-10PG
 PH-BSI-NSF-10PS
 PH-BSI-HC-10PG
 PH-BSI-HC-10PS
 †PH-BSI-HC-RFC12G
 †PH-BSI-HC-RFC12
 †PH-BSI-HC-RFC12A
 †PH-BSI-HC-RFC12A-CAD
 †PH-BSI-HC-RFC12GA
 †PH-BSI-HC-RFC12GA-CAD

*Model available in left hinge
 †Compatible with refrigerator compartment only