



## Silver Series Manual Defrost -45°C Chest Freezers

The BSI Silver Series manual defrost -45°C chest freezers are designed for laboratory & pharmaceutical applications that require more rigorous temperature control and construction features. These durable, easy to use units feature an intelligent microprocessor controller with digital display and adjustable temperature control.

### **Standard Product Features:**

- Digital temperature display
- Intelligent microprocessor controller with digital display and adjustable temperature control
- Visual and audible hi or low temperature, sensor failure, door open alarms
- Keyed door locks
- Magnetically sealed door with rapid cool-down mode after door opening for quick temperature recovery
- Cabinet constructed of white powder-coated 304 steel with CFC-free, high density insulation (2.4" thick)
- Easy swivel casters
- Doors: 1 Lid
- Defrost: Manual
- Built and shipped from ISO 9001:2008 manufacturing and distribution facilities
- Each cabinet undergoes three quality control checks prior to shipping
- One wire basket (B2X-7-C45 and B2X-11-C45 only)
- Certification: UL / cUL
- Warranty: Two-year parts & compressor



**B2X-3-C45**



**B2X-7-C45**



**B2X-11-C45**

Model	CAPACITY (cu. ft.)	Exterior Dimensions (HxWxD)	Interior Dimensions (HxWxD)	Temperature Range (Setpoint)	Shipping Weight
B2X-3-C45	3.5	33.5" x 27.25" x 22"	24.25" x 20.5" x 13.5"	-10 °C to -45°C C	102 lbs
B2X-7-C45	7	35" x 37.13" x 30.5"	24" x 29.6" x 19.1"	Preset -40 °C	113.5 lbs
B2X-11-C45	11	33.5" x 50.25" x 30.5"	24" x 42.38" x 19"		176 lbs