

# KEY ACCESSORIES

## REFRIGERATED INCUBATORS

*Models: 7001-10/25/33/50/75*



### One-Handed Door Latch

---

Secure yet easy chamber access. Pull handle to open, door latches automatically. Available with or without integral lock, for 25 & 33 ft<sup>3</sup> units.



### Extended Temperature Range

---

Extends unit operation down to 0°C, includes defrost cycle for 0-4.5°C



### Additional Shelving

---

Perforated, solid, wire, and heavy-duty variants. All shelves made of polished type 304 stainless steel and include shelf brackets. Factory-installed kits available to swap all shelves for another type.



### Outputs / Recorders

---

Circular temp/RH chart recorders in 6", 10" and 12" sizes, plus analog outputs and digital USB datalogging.