




MVE BL-7

The new MVE BL-7 (Biologicistic) is a secure cryogenic shipping option for your biological samples. The MVE BL-7 incorporates creative engineering and simplistic design to allow shipping in any orientation without sacrificing temperature or hold time. Improved welded construction provides the BL-7 with consistent, high quality design with every unit. **The MVE BL-7 allows for savings on dry ice, packaging, shipping cost, and disposal.**

Features include:

- Charges in 2 hours
- Operates as a dry shipper
- -150° C temperature
- 7-day holding time
- Safe shipping in any orientation
- Welded construction for added integrity
- Safe for international shipping; IATA Dangerous Goods Regulations Exempt (Special Provision A152)



	BL-7
Maximum Storage Capacity	
Number of Canisters	1
Number of 1/2 cc Straws in Bulk	115
Number of 1/4 cc Straws in Bulk	270
Number of 1.2 & 2.0 ml vials	15
Performance	
LN2 Capacity w/o Inventory (min) <i>L est.</i>	2.1
Static Holding Time (days)	7
Unit Dimensions	
Neck Opening <i>in. (mm)</i>	1.97 (50)
Overall Height <i>in. (mm)</i>	15.0 (381)
Overall Width <i>in. (mm)</i>	12.64 (321)
Overall Depth <i>in. (mm)</i>	12.64 (321)
Inner Canister Height <i>in. (mm)</i>	5.00 (127)
Inner Canister Diameter <i>in. (mm)</i>	1.50 (38)
Weight Empty <i>lb. (kg)</i>	9.26 (4.2)
Weight Charged <i>lb. (kg)</i>	13.01 (5.9)

TWO Year Vacuum Warranty

Conforms to MDD 93/42/EEC, the Medical Device Directive for the EU.

Visit www.chartbiomed.com for more information.



MVE BL-7

Safe, Secure, Cold.
No matter where it travels.