



**Class 10
Polyurethane Clean
Room Chair -
Medium Bench
Height with Medium
Back & Casters in
Black Polyurethane**

SKU: LRC43526

Details

Class 10 Clean Room Black Polyurethane Medium Bench Height Chair with Medium Back in Chrome is adjustable from 20.25-28 inches from the top of the seat to the floor and features ergonomic adjustments including seat height and locking back tilt, chrome adjustable foot ring, and cast-aluminum base with chrome-hooded, dual-wheel, self-braking, non-marring casters suitable for hard floors in chemical, soil, and puncture resistant polyurethane.

Specifications

Adjustment	Pneumatic Adjustment: 20.25-28 in
Seat	Width: 16in, Depth: 17in, Thickness: 1.5in
Back	Width: 15.5in, Height: 15in, Thickness: 1in
Footring	20 in diameter round tube chrome foot ring with internal lock
Color	Black

Seat Material	Rugged polyurethane that is self-skinning, chemical- and puncture-resistant, easily cleaned, and soil resistant
Feet	Chrome-hooded dual-wheel, self-braking, non-marring casters (suitable for both hard and carpeted floors)
Base	Polished Cast Aluminum Base
Base Width	24 in
Base Depth	24 in
Weight (lbs)	35
Weight Capacity	300 lbs
Warranty	Warranty: 15 year parts, lifetime on pneumatic lift