


Weight load capacity for BC and Deep Ledge series.  
 150 lbs. (68kg) per shelf.  
 2-shelf model: 300 lbs. (136kg)  
 3-shelf model: 400 lbs. (181kg)



Gray



Blue



Black



Slate Blue  
(with Microban product protection)

**BC Series Utility Carts — 12.28 (Available through July, 2013 —see pages 84-85)**



- Polymer shelves are easy to clean and corrosion proof.
- 7/16" (11mm) ship's edge lip around shelf perimeter contains spills and helps prevent items from sliding off during transport.
- Available colors are gray, black, blue and slate blue. Slate blue contains Microban antimicrobial product protection.
- Available in 2-shelf and 3-shelf styles.
- Center shelf option adjusts at 1" (25mm) increments.

Width/Length/Height		Description	Approx. Pkd. Wt.		Cat. No.*	Cat. No. Slate Blue with Microban
(in.)	(mm)		(lbs.)	(kg)		
18x28x33 <sup>1</sup> / <sub>4</sub>	457x710x845	2-shelf unit	27 <sup>1</sup> / <sub>2</sub>	12.4	LAB42-7261CB	LABBM42-7261CB
18x28x33 <sup>1</sup> / <sub>4</sub>	457x710x845	3-shelf unit	33	14.9	LAB43-7261CB	LABBM43-7261CB
21 <sup>1</sup> / <sub>2</sub> x33 <sup>3</sup> / <sub>4</sub> x33 <sup>1</sup> / <sub>4</sub>	546x857x845	2-shelf unit	34 <sup>1</sup> / <sub>2</sub>	15.6	LAB42-0302CB	LABBM42-0302CB
21 <sup>1</sup> / <sub>2</sub> x33 <sup>3</sup> / <sub>4</sub> x33 <sup>1</sup> / <sub>4</sub>	546x857x845	3-shelf unit	43	19.5	LAB43-0302CB	LABBM43-0302CB
27x39 <sup>1</sup> / <sub>2</sub> x33 <sup>1</sup> / <sub>4</sub>	686x1003x845	2-shelf unit	40 <sup>1</sup> / <sub>2</sub>	18.3	LAB42-6362CB	
27x39 <sup>1</sup> / <sub>2</sub> x33 <sup>1</sup> / <sub>4</sub>	686x1003x845	3-shelf unit	52	23.5	LAB43-6362CB	

\*Add color designation to part number when ordering: Gray (G), Black (BL) and Blue (BU).



3-Shelf Deep Ledge Utility Cart

**Deep Ledge Utility Carts — 12.28**

Specially designed with a 2<sup>3</sup>/<sub>4</sub>" (70mm) deep ledge to contain product and spills.

- Polymer shelves are easy to clean and corrosion proof.
- Available colors are gray, black, blue and slate blue. Slate blue contains Microban antimicrobial product protection.
- Available in 2-shelf and 3-shelf models.
- Center shelf option adjusts at 1" (25mm) increments.

Width/Length/Height		Description	Approx. Pkd. Wt.		Cat. No.*	Cat. No. Slate Blue with Microban
(in.)	(mm)		(lbs.)	(kg)		
21 <sup>1</sup> / <sub>2</sub> x33 <sup>3</sup> / <sub>4</sub> x33 <sup>1</sup> / <sub>4</sub>	546x857x845	2-shelf unit	29	31.1	LABD2-0302CB	LABBMD2-0302CB
21 <sup>1</sup> / <sub>2</sub> x33 <sup>3</sup> / <sub>4</sub> x33 <sup>1</sup> / <sub>4</sub>	546x857x845	3-shelf unit	37 <sup>1</sup> / <sub>2</sub>	17.0	LABD3-0302CB	LABBMD3-0302CB
27x39 <sup>1</sup> / <sub>2</sub> x33 <sup>1</sup> / <sub>4</sub>	685x1003x845	2-shelf unit	35	15.9	LABD2-6362CB	
27x39 <sup>1</sup> / <sub>2</sub> x33 <sup>1</sup> / <sub>4</sub>	685x1003x845	3-shelf unit	46 <sup>1</sup> / <sub>2</sub>	21.1	LABD3-6362CB	

\*Add color designation to part number when ordering: Gray (G), Black (BL) and Blue (BU).



**Accessories and Replacement Casters — 12.28**

Description	Fits Cart	Width/Height/Depth		Approx. Pkd. Wt.		BC Series Cat. No.	Deep Ledge Cat. No.
		(in.)	(mm)	(lbs.)	(kg)		
Wastebasket + Holder	BC1627	17 <sup>3</sup> / <sub>8</sub> x27 <sup>1</sup> / <sub>2</sub> x13 <sup>1</sup> / <sub>4</sub>	440x700x337	5.5	2.5	LAB1BWCB	LABD1BWCB
Wastebasket + Holder	BC2030	21 <sup>1</sup> / <sub>4</sub> x27 <sup>1</sup> / <sub>2</sub> x13 <sup>1</sup> / <sub>4</sub>	540x700x337	8.0	3.8	LAB2BWCB	LABD2BWCB
Utility Bin + Holder	BC1627	17 <sup>1</sup> / <sub>2</sub> x7 <sup>3</sup> / <sub>4</sub> x14 <sup>1</sup> / <sub>2</sub>	445x197x368	3.8	1.7	LAB1BUCB	LABD1BUCB
Utility Bin + Holder	BC2030	21 <sup>1</sup> / <sub>4</sub> x7 <sup>3</sup> / <sub>4</sub> x14 <sup>1</sup> / <sub>2</sub>	540x197x368	6.0	2.7	LAB2BUCB	LABD2BUCB
Wastebasket Only		14 <sup>3</sup> / <sub>8</sub> x15 <sup>3</sup> / <sub>8</sub> x10 <sup>3</sup> / <sub>8</sub>	365x391x264	2.8	1.3	LAB222FM	LAB222FM
Utility Bin Only		16 <sup>3</sup> / <sub>8</sub> x6x11	416x152x279	2.0	0.9	LAB1BU	LAB1UB
Replacement Caster Kit — All						LAB4-M4CBPR	LAB4-M4CBPR

\*MICROBAN® and the MICROBAN® symbol are registered trademarks of the Microban Products Company, Huntersville, NC.