



MINI VORTEX MIXERS

Crystal Technology and Industries' Mini Vortex Mixer is the ideal Vortex Mixer for every lab with its compact size, ease of use, and its economic price point. The 3,000 rpm fixed speed ensures proper mixing and accurate results.



Mini Vortex Mixer
VM-03UA



Mini Vortex Mixer
VM-03GUA



Mini Vortex Mixer
VM-03RUA



Mini Vortex Mixer
VM-03PUA

PRODUCT FEATURES:

- Powerful vortexing for tubes up to 30mm
- Extremely compact and easy operation
- Vortex head is made of durable material to ensure long life
- Refined steel base is designed to ensure product stability
- Choice of different colored working head

Specifications

Model	VM-03UA	VM-03RUA	VM-03GUA	VM-03PUA
Product Name	Mini Vortex Mixer			
Cap Color	Blue	Red	Green	Purple
Speed	3000 RPM			
Speed Accuracy	±10%			
Oscillation	Circular Oscillation			
Orbital Diameter	Ø4.5mm (Ø0.18 in)			
Capacity	Max 80mL			
Max Tube Diameter	≤Ø30mm (≤Ø1.2 in)			
Operating Environment	Ambient Temperature: 5°C~40°C Relative Humidity: ≤ 80%			
Power Supply	AC110V, 60Hz			
Max. Power	5W			
Dimensions DiaxH	Ø96 x 85 mm (Ø3.8 x 3.4 in)			
Net Weight	0.6 kg (1.2 lb)			
Protection Level	IP21			
Certifications	CE, RoHS, WEEE			

