



## **Labsit Chair 9123, Desk Height, Freewheeling Casters, Black Vinyl, Armless**

**Article #: LRC-LBS23\***

*\*In-Stock Quick Ship Item*

---

### **Features**

Labsit is a progressive chair design for all types of laboratories. The intelligent inherent flex of the seat shell offers comfort and ergonomic support with minimum adjustment. Lightweight and easy to maneuver, Labsit is a winning solution for lab stations, ESD areas and classrooms with an unbeatable value and performance.

---

### **Specifications**

<b>Seat Adjustment:</b>	Standard Desk Height Range 17-3/4" - 25-1/2"
<b>Seat Depth:</b>	18-1/8"
<b>Shell Colors:</b>	In-Stock (Blue); Factory Orders (Orange, Green, White, Black)
<b>Seat Material:</b>	Black Vinyl
<b>Castors:</b>	Freewheeling
<b>Base:</b>	Black Plastic Base
<b>Cleanroom Rating:</b>	ISO Level 5
<b>Bio-Safety BSL:</b>	Level 2
<b>BIFMA Rating:</b>	Complies with Level e3 Platinum
<b>SEFA Rating:</b>	Complies with Scientific Equipment & Furniture Assn. - SEFA12
<b>Weight Capacity:</b>	300 lbs.
<b>Product Weight:</b>	25 lbs.
<b>Warranty:</b>	10-Year Warranty
<b>Original Design:</b>	German Engineering and Production
<b>ESD Capable:</b>	Option Available by Request