



**Labsit Chair 9123,
Desk Height,
Freewheeling Casters,
Black Vinyl, With Arms**

Article #: LRC-LBS23A

**In-Stock Quick Ship Item*

Features

Labsit is a progressive chair design for all types of laboratories. The intelligent inherent flex of the seat shell offers comfort and ergonomic support with minimum adjustment. Lightweight and easy to maneuver, Labsit is a winning solution for lab stations, ESD areas and classrooms with an unbeatable value and performance.

Specifications

Seat Adjustment:	Standard Desk Height Range 17-3/4" - 25-1/2"
Seat Depth Range:	18-1/8"
Shell Colors:	In-Stock (Blue); Factory Orders (Orange, Green, White, Black)
Seat Material:	Vinyl Black
Castors:	Freewheeling
Base:	Black Base
Cleanroom Rating:	ISO Level 5
Bio-Safety BSL:	Level 2
BIFMA Rating:	Complies with Level e3 Platinum
SEFA Rating:	Complies with Scientific Equipment & Furniture Assn SEFA12
Weight Capacity:	300 lbs.
Product Weight:	25 lbs.
Warranty:	10-Year Warranty
Original Design:	German Engineering and Production
ESD Capable:	Option Available by Request