

TW 35LD

LN₂ DEWAR



For safe storage and dispensing of small amounts of liquid nitrogen.



| SPECIFICATIONS | 35LD |
|------------------------------|------|
| LN ₂ Capacity (L) | 35 |
| Static Hold Time (Days) | 152 |
| Evap. Rate (Liters/Day) | 0.23 |
| DIMENSIONS | |
| Necktube (in) | 2.5 |
| Overall Dia. (in) | 18.8 |
| Height (in) | 26.3 |
| WEIGHT | |
| Empty (Lbs) | 35.1 |
| Full (Lbs) | 97.4 |

FEATURES

- 35 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA