



LD Series

Storage Dewar

TW 4LD

LN₂ DEWAR



For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	4LD
LN ₂ Capacity (L)	4
Static Hold Time (Days)	10
Evap. Rate (Liters/Day)	0.40
DIMENSIONS	
Necktube (in)	1.2
Outside Dia. (in)	7.6
Height (in)	17
WEIGHT	
Empty (Lbs)	6.6
Full (Lbs)	13.7

FEATURES

- 4 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA

TW 5LD

LN₂ DEWAR



For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	5LD
LN ₂ Capacity (L)	5
Static Hold Time (Days)	6
Evap. Rate (Liters/Day)	0.77
DIMENSIONS	
Necktube (in)	5.6
Overall Dia. (in)	7.6
Height (in)	17.5
WEIGHT	
Empty (Lbs)	6.9
Full (Lbs)	15.8

FEATURES

- 5 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA

TW 10LD

LN₂ DEWAR



For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	10LD
LN ₂ Capacity (L)	10
Static Hold Time (Days)	45
Evap. Rate (Liters/Day)	0.22
DIMENSIONS	
Necktube (in)	2.0
Overall Dia. (in)	11.4
Height (in)	23.5
WEIGHT	
Empty (Lbs)	14.5
Full (Lbs)	32.3

FEATURES

- 10 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA

TW 25LD

LN₂ DEWAR



For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	25LD
LN ₂ Capacity (L)	25
Static Hold Time (Days)	109
Evap. Rate (Liters/Day)	0.23
DIMENSIONS	
Necktube (in)	2.5
Overall Dia. (in)	15.6
Height (in)	25.8
WEIGHT	
Empty (Lbs)	23.2
Full (Lbs)	67.7

FEATURES

- 25 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA

TW CLASSIC-25

LN₂ DEWAR



For safe storage and dispensing of small accounts of liquid nitrogen.



SPECIFICATIONS	CLASSIC-25
LN ₂ Capacity (L)	25
Static Hold Time (Days)	119
Evap. Rate (Liters/Day)	0.21
DIMENSIONS	
Necktube (in)	2.0
Overall Dia. (in)	15.5
Height (in)	22.9
WEIGHT	
Empty (Lbs)	19
Full (Lbs)	63.5

FEATURES

- 25 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA

TW 35LD

LN₂ DEWAR



For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	35LD
LN ₂ Capacity (L)	35
Static Hold Time (Days)	152
Evap. Rate (Liters/Day)	0.23
DIMENSIONS	
Necktube (in)	2.5
Overall Dia. (in)	18.8
Height (in)	26.3
WEIGHT	
Empty (Lbs)	35.1
Full (Lbs)	97.4

FEATURES

- 35 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA

TW 50LD

LN₂ DEWAR



For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	50LD
LN ₂ Capacity (L)	50
Static Hold Time (Days)	122
Evap. Rate (Liters/Day)	0.41
DIMENSIONS	
Necktube (in)	2.5
Overall Dia. (in)	18.8
Height (in)	32.4
WEIGHT	
Empty (Lbs)	38.7
Full (Lbs)	127.7

FEATURES

- 50 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA