



MetroMax i® with Open Grid Shelves



MetroMax i® with Solid Shelves

MetroMax i® — Stem Caster Carts

- Corrosion proof shelves and posts; removable polymer mats.
- Standard Units consist of four shelves, MX63UP posts 63" (1600mm), wedge connectors and choice of resilient rubber or polyurethane casters.
- Approximate cart height — 67⁵/₁₆" (1710mm). Weight capacity of cart: 900 lbs. (408kg) evenly distributed.

MetroMax i® Open Grid

Nominal Width/Length (in.)	Nominal Width/Length (mm)	Approx. Pkd. Wt.		Cat. No. 2-Swivel 2-Brake Resilient Rubber Casters		Cat. No. 2-Swivel 2-Brake Polyurethane Casters	
		(lbs.)	(kg)				
18x36	457x914	74	33	LAB3XGB633X	LAB3XGE633X		
18x48	457x1220	86	38	LAB3XGB653X	LAB3XGE653X		
18x60	457x1524	98	44	LAB3XGB663X	LAB3XGE663X		
24x36	610x914	82	37	LAB3XGB635X	LAB3XGE635X		
24x48	610x1220	98	44	LAB3XGB655X	LAB3XGE655X		
24x60	610x1524	114	51	LAB3XGB665X	LAB3XGE665X		

MetroMax i® Solid Shelves

Approx. Pkd. Wt.		Cat. No. 2-Swivel 2-Brake Polyurethane Casters
(lbs.)	(kg)	
94	42	LAB3XFE633X
106	48	LAB3XFE653X
118	53	LAB3XFE663X
102	46	LAB3XFE635X
118	53	LAB3XFE655X
134	60	LAB3XFE665X



MetroMax Q® with Open Grid Shelves

MetroMax Q® — Stem Caster Cart, Open Grid Shelves

- Quick adjust shelves with removable polymer shelf mats.
- Standard Units consist of four open grid shelves, MQ63UPE posts 63" (1600mm), wedge connectors, and choice of resilient rubber or polyurethane casters.
- Approximate cart height — 67⁵/₁₆" (1710mm). Weight capacity of cart: 900 lbs. (408kg) evenly distributed.

Nominal Width/Length (in.)	Nominal Width/Length (mm)	Approx. Pkd. Wt.		Cat. No. 2-Swivel 2-Brake Resilient Rubber Casters		Cat. No. 2-Swivel 2-Brake Polyurethane Casters	
		(lbs.)	(kg)				
18x36	457x914	68	31	LAB3GB633Q	LAB3GE633Q		
18x48	457x1220	84	38	LAB3GB653Q	LAB3GE653Q		
18x60	457x1524	100	45	LAB3GB663Q	LAB3GE663Q		
21x36	530x914	74	34	LAB3GB634Q	LAB3GE634Q		
21x48	530x1220	90	41	LAB3GB654Q	LAB3GE654Q		
21x60	530x1524	106	48	LAB3GB664Q	LAB3GE664Q		
24x36	610x914	80	36	LAB3GB635Q	LAB3GE635Q		
24x48	610x1220	96	44	LAB3GB655Q	LAB3GE655Q		
24x60	610x1524	112	51	LAB3GB665Q	LAB3GE665Q		



MetroMax i® with a Solid Bottom Shelf

Stem Caster Carts with Solid Bottom MetroMax i® Shelves

- 5-tier models with Microban® antimicrobial product protection.
- Standard Units consist of 5 shelves (4 open grid, 1 MetroMax i solid), 74" (1880mm) posts, wedge connectors, and polyurethane casters.
- Approximate cart height — 79⁵/₁₆" (2015mm). Weight capacity of cart: 900 lbs. (408kg) evenly distributed.

Nominal Width/Length (in.)	Nominal Width/Length (mm)	Approx. Pkd. Wt.		i 2-Swivel 2-Brake Polyurethane Casters Cat. No.		Q 2-Swivel 2-Brake Polyurethane Casters Cat. No.	
		(lbs.)	(kg)				
18x36	457x914	82	37	LAB3XGE733X5	85	39	LAB3GE733Q5
18x48	457x1220	109	49	LAB3XGE753X5	111	50	LAB3GE753Q5
18x60	457x1524	121	55	LAB3XGE763X5	123	56	LAB3GE763Q5
24x36	610x914	102	46	LAB3XGE755X5	104	47	LAB3GE735Q5
24x48	610x1220	123	56	LAB3XGE755X5	125	57	LAB3GE755Q5
24x60	610x1524	141	63	LAB3XGE765X5	144	65	LAB3GE765Q5

✓ Indicates antimicrobial product.