

Upright Freezer Racks New

Modifiable Racks for 2" and 3" Standard Boxes

- ◆ Removable clip can be adjusted so that the rack can hold both standard 2" and 3" boxes.
- ◆ Constructed of corrosion resistant stainless steel.



Part Number	Dimension L x H x W	Capacity (Boxes)			Shipping Weight(lb/Kg)
		2"	3"	2"and 3"	
UFMC-33-2-3	16 ⁷ / ₁₆ x 6 ⁷ / ₈ x 5 ¹ / ₂ inch 418.2 x 168.8 x 140 mm	9	6	3 + 3	5.5/2.5
UFMC-34-2-3	16 ⁷ / ₁₆ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 418.2 x 225 x 140 mm	12	6	6 + 3	6.7/3.0
UFMC-35-2-3	16 ⁷ / ₁₆ x 11 x 5 ¹ / ₂ inch 418.2 x 280 x 140 mm	15	6	9 + 3	7.8/3.6
UFMC-43-2-3	21 ⁹ / ₁₆ x 6 ⁷ / ₈ x 5 ¹ / ₂ inch 557 x 168.8 x 140 mm	12	8	4 + 4	7.1/3.2
UFMC-44-2-3	21 ⁹ / ₁₆ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 557 x 225 x 140 mm	16	8	8 + 4	8.6/3.9
UFMC-45-2-3	21 ⁹ / ₁₆ x 11 x 5 ¹ / ₂ inch 557 x 280 x 140 mm	20	8	12 + 4	10.1/4.6
UFMC-53-2-3	26 ⁹ / ₁₆ x 6 ⁷ / ₈ x 5 ¹ / ₂ inch 681 x 168.8 x 140 mm	15	10	5 + 5	8.7/3.9
UFMC-54-2-3	26 ⁹ / ₁₆ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 681 x 225 x 140 mm	20	10	10 + 5	10.6/4.8
UFMC-55-2-3	26 ⁹ / ₁₆ x 11 x 5 ¹ / ₂ inch 681 x 280 x 140 mm	25	10	15 + 5	12.2/5.6

For Mixed Storage of 2", 3" New and 3.75" Standard Boxes

- ◆ Constructed of corrosion resistant stainless steel.



Part Number	Dimension L x H x W	Capacity (Boxes)			Shipping Weight(lb/Kg)
		2"	3"	2"and 3"	
UF-306	16 ³ / ₁₆ x 6 ³ / ₄ x 5 ¹ / ₂ Inch 415 x 159 x 140 mm	3x3 =9	3x2 =6	3x1 + 3x1	3.3/1.5
UF-308	16 ³ / ₁₆ x 8 ³ / ₄ x 5 ¹ / ₂ Inch 415 x 210 x 140 mm	3x4 =12	3x2 =6	3x1 + 3x2	3.7/1.7
UF-312L	16 ³ / ₁₆ x 12 ³ / ₈ x 5 ¹ / ₈ Inch 415 x 315 x 143 mm	3x6 =16	3x4 =12	3x3 + 3x2	4.9/2.2
UF-406	21 ³ / ₈ x 6 ³ / ₄ x 5 ¹ / ₂ Inch 550 x 159 x 140 mm	4x3 =12	4x2 =8	4x1 + 4x1	4.2/1.9
UF-408	21 ³ / ₈ x 8 ³ / ₄ x 5 ¹ / ₂ Inch 550 x 210 x 140 mm	4x4 =16	4x2 =8	4x1 + 4x2	4.8/3.9
UF-412L	21 ³ / ₈ x 12 ³ / ₈ x 5 ¹ / ₈ Inch 550 x 315 x 143 mm	4x6 =24	4x3 =12	4x3 + 4x2	6.3/2.9
UF-506	26 ³ / ₄ x 6 ³ / ₄ x 5 ¹ / ₂ Inch 679 x 159 x 140 mm	5x3 =15	5x2 =10	5x1 + 5x1	8.7/3.9
UF-508	26 ³ / ₄ x 8 ³ / ₄ x 5 ¹ / ₂ Inch 679 x 210 x 140 mm	5x4 =15	5x2 =10	5x1 + 5x2	5.1/2.4
UF-512L	26 ³ / ₄ x 12 ³ / ₈ x 5 ¹ / ₈ Inch 679 x 315 x 143 mm	5x6 =30	5x3 =15	5x3 + 5x2	7.7/3.5