



# NEW Model NU-99100J



## Ultralow Under-Counter Freezer

**-86°C ULTRA-LOW TEMPERATURE FREEZER**

The NuAire Mini-Blizzard -86°C Ultralow Freezer features reliable performance for the storage of stem cells, blood & tissue samples, and vaccine vials while occupying the smallest space among our ultralow freezers. The NU-99100J model now features uniform temperature control up to 1°C, a simplified LED display screen, a quadruple gasketed outer door, 2 interior steel mini-doors, and produces visual/audible alarms. The quiet compressor motor utilizes non-toxic Hydrocarbon refrigerants and operates under 47 dB(A) while continuously maximizing energy efficiency.



### Storage Capacity: 100 L (3.5 ft<sup>3</sup>)



**1 Box = 195 Vials**

Box Width: 9" [229 mm]  
Box Depth: 9" [229 mm]  
Box Height: 1.6" [41 mm]



**1 Vial = 5 Doses**

*\*Note: 2 Doses are needed to inoculate 1 person*

*\*Note: Vaccine storage capacity is calculated by manufacturer's provided measurements, actual capacity may vary.*

Model	Configuration	Storage Capacity	Interior Dimensions (W x D x H)	Exterior Dimensions (W x D x H)	Temperature Range	Net Weight	Electrical Configurations
NU-99100J	Under-Counter	3.5 cu ft [100 L]	13 x 18.9 x 24.8 in [330 x 481 x 630 mm]	30.2 x 26 x 31.9 in [770 x 660 x 810 mm]	-40°C to -86°C	238 lbs [108 kg]	120 VAC/60 Hz 240 VAC/50 Hz

More Information or Request a Quote at [www2.nuaire.com/blizzard-ultralow-freezer-nu-99100j](http://www2.nuaire.com/blizzard-ultralow-freezer-nu-99100j)

**3 YEAR WARRANTY**  
PARTS + LABOR  
United States & Canada

**4 YEAR WARRANTY**  
PARTS  
Worldwide

