

Medical-Grade Refrigerators and Freezer Manufacturers - Limited Storage Capacity

Manufacturer	ROAM Mobile	Helmer Traditional Stationary Undercounters	Thermo/Revco Traditional Stationary Undercounters	NorLake Traditional Stationary Undercounters	Follett Traditional Stationary Undercounters
Mobility	Yes, carried or rolled	None, Stationary	None, Stationary	None, Stationary	None, Stationary
Size (Height)	13.5" to 22"	34"	33 3/8"	31 5/8"	34"
Weight	26 lbs - 95 lbs	189 lbs	212 lbs	125 lbs	205 lbs
Carry with handles	Yes	No	No	No	No
Storage Capacity (liters)	3 liters to 100 liters	142 liters	133 liters	119 liters	127 liters
Battery Power Capable	Yes, Optional	No	No	No	No
Payload Risk (Blood, Plasma, Vaccines, Pharm, Reagents)	Very Low (w/Optional Battery)	Medium, if power outage and using generator	Medium, if power outage and using generator	Medium, if power outage and using generator	Medium, if power outage and using generator
Data Logging	Yes, Optional	Yes, Optional	Yes, Optional	Yes, Optional	Yes, Optional
Warranty	Yes	Yes	Yes	Yes	Yes
Audible/Visual Alarms	Yes	Yes	Yes	Yes	Yes
Shipping Damage	<1%	3-5%	3-5%	3-5%	3-5%
Outdoor Use	Yes	No	No	No	No
Service	Ship to San Diego	Local Heating and Air Professionals	Local Heating and Air Professionals	Local Heating and Air Professionals	Local Heating and Air Professionals
Heat Capable	Yes	No	No	No	No
Made in USA	Yes	Yes	Yes	Yes	Yes
115v Electrical Power Plug	Yes	Yes	Yes	Yes	Yes