



**Polyurethane
Tractor Sit-Stand
Stool - Medium
Bench Height with
Chrome Foot Ring &
Casters in Blue
Polyurethane**

SKU: BSI45152

Details

Blue Polyurethane Tractor Sit-Stand Medium Bench Height Stool is adjustable from 19-26.25 inches from the top of the seat to the floor and features ergonomic adjustments such as seat height and locking seat tilt, chrome adjustable foot ring, and black, nylon-reinforced fiberglass base with dual-wheel, self-braking, non-marring casters suitable for hard floors in chemical, soil, and puncture resistant polyurethane.

Specifications

Adjustment	Pneumatic Adjustment: 19-26.5 in
Seat	Width: 17in, Depth: 14in, Thickness: 2.5in
Footring	18 in diameter round tube chrome foot ring with knob lock
Color	Blue
Seat Material	Rugged polyurethane that is self-skinning, chemical- and puncture-resistant, easily cleaned, and soil resistant

Feet	Black dual-wheel, self-braking, non-marring casters (suitable for both hard and carpeted floors)
Base	Nylon Reinforced Fiberglass Base
Base Width	25
Base Depth	25
Weight (lbs)	27
Weight Capacity	300 lbs
Warranty	Warranty: 15 year parts, lifetime on pneumatic lift

