

















# Accessories

|   |   |   |   |   |  |   |   |   |   |   |  |
|---|---|---|---|---|--|---|---|---|---|---|--|
|  |  |  |  |  |  |  |  |  |  |  |  |
| <b>Cat. No.</b>   | <b>18900001</b>   | <b>18900002</b>   | <b>18900003</b>   | <b>18900004</b>   | <b>18900005</b>  | <b>189000048</b>  | <b>189000049</b>  | <b>189000080</b>  | <b>189000081</b>  | <b>189000082</b>  | <b>189000083</b>   |
| <b>Descriptions</b>   | Blue carrying plate, use with color quarter pies                                    | Red quarter pie, 11 holes, 4 ml reaction vessel, ±15.2mm, 20mm depth                | Purple quarter pie, 4 holes, 20 ml reaction vessel, ±28mm, 24mm depth               | Blue quarter pie, 4 holes, 30 ml reaction vessel, ±28mm, 30mm depth                 | Black quarter pie, 4 holes, 40 ml reaction vessel, ±28mm, 43mm depth               | Green quarter pie, 6 holes, 8ml reaction vessel, ±17.8mm, 26mm depth              | Golden quarter pie, 4 holes, 16ml reaction vessel, ±21.6mm, 31.7mm depth          | Reaction block for 50ml round bottom flask (one flask capacity)                   | Reaction block for 100ml round bottom flask (one flask capacity)                  | Reaction block for 250ml round bottom flask (one flask capacity)                  | Reaction block for 500ml round bottom flask (one flask capacity) |



MS-H-Pro<sup>+</sup> with PT1000 and 1/4 color pies

| <b>Models</b>         |  |  |  |  |  |  |  |  |  |  |  |
|-----------------------|--|--|--|--|--|--|--|--|--|--|--|
| MS7-H550-Pro          |  |  |  |  |  |  |  |  |  |  |  |
| MS7-H550-S            |  |  |  |  |  |  |  |  |  |  |  |
| MS-H-Pro <sup>+</sup> |  |  |  |  |  |  |  |  |  |  |  |
| MS-H-S                |  |  |  |  |  |  |  |  |  |  |  |
| MS-H280-Pro           |  |  |  |  |  |  |  |  |  |  |  |

|  |  |  |                  |  |
|--|--|---|--|--|
| <b>Cat. No.</b>  | <b>18900019</b>  | <b>18900016</b>   | <b>189000136</b>   | <b>189000148</b>   |
| <b>Descriptions</b>  | F101 protective cover, for hotplate models, temperature resistance of 135°C          | PTI 000-A, temperature sensor for digital hotplate models, length of 230mm          | PTI 000-B, temperature sensor with glass coated, length of 230mm, used for digital hotplate models | Support clamp of PT1000, use with MS-H280-PRO                                      |
| <b>Models</b>  |  |   |  |  |
| MS7-H550-Pro   |  |   |  |  |
| MS7-H550-S   |  |   |  |  |
| MS-H-Pro <sup>+</sup>  |  |   |  |  |
| MS-H-S   |  |   |  |  |
| MS-H280-Pro  |  |   |  |  |



MS-H-Pro<sup>+</sup> with PT1000 and reaction block