


Material : Stainless Steel


Upright Freezer Racks

For Standard 2" Boxes

Racks are for standard 2" cardboard (white or yellow), polycarbonate or metal (stainless steel or aluminum) boxes. High quality Upright Freezer Racks, including handles, are constructed of corrosion resistant stainless steel and offer convenient storage of samples in an organized format. Security lock devices are also available for each rack. Using the internal dimensions obtained for your freezer shelf space (see page 1, fig 1) determine the height (H), length (L) and width (W) of each rack you select for your upright freezer.



Rack Only



Rack w/Cardboard Box



Rack w/Metal Box



Rack w/Polycarbonate Plastic Box



w/DP Lock Device Option



Yellow 2" Cardboard Box



White 2" Cardboard Box



Polycarbonate Plastic Box















Metal Box
S-Stainless Steel
A-Aluminum

[Check page 81 for Boxes]

Material : Stainless Steel

Upright Freezer Racks

For Standard 2"Boxes

Part Number	Description	Dimension L x H x W	Capacity (Box)	Lock Device Option***
UF-332	Rack Only	16 ⁵ / ₁₆ x 6 ⁷ / ₈ x 5 ¹ / ₂ inch 415 x 168 x 140 mm	3 x 3 =9 	DP-1607
UF-332-X'-XX''	Rack w/2" box and divider			
UF-342	Rack Only	16 ⁵ / ₁₆ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 415 x 225 x 140 mm	3 x 4 =12 	DP-1609
UF-342-X'-XX''	Rack w/2" box and divider			
UF-352	Rack Only	16 ⁵ / ₁₆ x 11 x 5 ¹ / ₂ inch 415 x 280 x 140 mm	3 x 5 =15 	DP-1611
UF-352-X'-XX''	Rack w/2" box and divider			
UF-362	Rack Only	16 ⁵ / ₁₆ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 415 x 321 x 140 mm	3 x 6 =18 	DP-1613
UF-362-X'-XX''	Rack w/2" box and divider			
UF-432	Rack Only	21 ¹ / ₈ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 550 x 168 x 140 mm	4 x 3 =12 	DP-2107
UF-432-X'-XX''	Rack w/2" box and divider			
UF-442	Rack Only	21 ¹ / ₈ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 550 x 225 x 140 mm	4 x 4 =16 	DP-2109
UF-442-X'-XX''	Rack w/2" box and divider			
UF-452	Rack Only	21 ¹ / ₈ x 11 x 5 ¹ / ₂ inch 550 x 280 x 140 mm	4 x 5 =20 	DP-2111
UF-452-X'-XX''	Rack w/2" box and divider			
UF-462	Rack Only	21 ¹ / ₈ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 550 x 321 x 140 mm	4 x 6 =24 	DP-2113
UF-462-X'-XX''	Rack w/2" box and divider			
UF-472	Rack Only	21 ¹ / ₈ x 15 ³ / ₁₀ x 5 ¹ / ₂ inch 550 x 391 x 140 mm	4 x 7 =28 	DP-2115
UF-472-X'-XX''	Rack w/2" box and divider			
UF-532	Rack Only	26 ³ / ₄ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 679 x 168 x 140 mm	5 x 3 =15 	DP-2707
UF-532-X'-XX''	Rack w/2" box and divider			
UF-542	Rack Only	26 ³ / ₄ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 679 x 225 x 140 mm	5 x 4 =20 	DP-2709
UF-542-X'-XX''	Rack w/2" box and divider			
UF-552	Rack Only	26 ³ / ₄ x 11 x 5 ¹ / ₂ inch 679 x 280 x 140 mm	5 x 5 =25 	DP-2711
UF-552-X'-XX''	Rack w/2" box and divider			

* To order, replace X with box type: Y-Yellow cardboard box; W-White cardboard box; P-Plastic Polycarbonate box; PP-Plastic Polypropylene Box; S-Stainless Steel Box or A-Aluminum Box.

** To order, replace XX with 16,25,36,49,64,81,100 or 144 cell divider. Plastic Polycarbonate and Polypropylene Boxes only have 81 and 100 cell dividers.

*** Lock device is ordered separately.

UFD-462 and UFD-462-X'-XX'' does not accept Plastic Polycarbonate nor Plastic Polypropylene Boxes.