


Material : Stainless Steel

Chest Freezer Racks For 25-Place Small Slide Boxes



L1 x W1 x H1: 5.5" x 3.3" x 1.4"
140 x 83 x 34mm

Racks listed are for chest freezers to store 25-Place Small Slide Boxes. Racks are made of corrosion resistant stainless steel. All racks include stainless steel locking rods and handles. Security locks are also available for each rack. Using the internal dimensions obtained for your chest freezer (see page 1 , fig 2) determine the maximum height (H), length (L) and width (W) and the number of racks you need for your chest freezer.



Rack Only

Part Number	Dimension H x D x W	Capacity	Lock Device Option*
CFS-25-4	155x152x97mm 6 ¹ / ₁₀ x 6x 3 ⁴ / ₅ inch	4	DR25-06
CFS-25-6	232x152x97mm 9 ¹ / ₁₀ x 6x 3 ⁴ / ₅ inch	6	DR25-09
CFS-25-8	309x152x97mm 12 ¹ / ₅ x 6x 3 ⁴ / ₅ inch	8	DR25-12
CFS-25-10	386x152x97mm 15 ¹ / ₅ x 6x 3 ⁴ / ₅ inch	10	DR25-15
CFS-25-12	463x152x97mm 18 ¹ / ₅ x 6x 3 ⁴ / ₅ inch	12	DR25-18
CFS-25-14	540x152x97mm 21 ³ / ₁₀ x 6x 3 ⁴ / ₅ inch	14	DR25-21
CFS-25-16	617x152x97mm 24 ³ / ₁₀ x 6x 3 ⁴ / ₅ inch	16	DR25-24
CFS-25-18	694x152x97mm 27 ³ / ₁₀ x 6x 3 ⁴ / ₅ inch	18	DR25-27
CFS-25-20	771x152x97mm 30 ³ / ₅ x 6x 3 ⁴ / ₅ inch	20	DR25-30