


Upright Freezer Rack



For Matrix

1.0ml Screw-Top Latch Rack (H=2¹/₈ inch / 54.8mm)

H: 2⁵/₁₆" - 2⁷/₁₆" (54.8-61.8mm)

For Micronic

Roborack-96 with High Cover (H=2¹/₄ inch / 56.1mm)

For Fluidx

96-Format Rack for Tube with Protective Jacket with Tall Profile Lid (H=2³/₁₆ inch / 55mm)


96-Format Rack Tack for 1.0ml and 1.4ml Tube with Tall Profile Lid (H=2⁷/₁₆ inch / 61.8mm)

24-Format Rack (H=2¹/₄ inch / 56mm)

For Abgene

1.2ml 2D Barcoded Tube in Stackable Rack (H=2¹/₄ inch / 56.13mm)

1.2ml 2D Barcoded Storage Tube Rack (H=2¹/₄ inch / 56.13mm)



- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.

Part Number	Dimension L x H x W	Capacity	Shipping Weight(lb/Kg)
UFMX-332	10 ¹¹ / ₁₆ x 7 ² / ₅ x 5 ⁷ / ₁₆ inch 271 x 191 x 138 mm	3 x 3 = 9	9.5/4.3
UFMX-342	10 ¹¹ / ₁₆ x 10 x 5 ⁷ / ₁₆ inch 271 x 254 x 138 mm	3 x 4 = 12	12.5/5.7
UFMX-352	10 ¹¹ / ₁₆ x 12 ¹ / ₂ x 5 ⁷ / ₁₆ inch 271 x 317 x 138 mm	3 x 5 = 15	15.3/6.9
UFMX-432	14 ³ / ₁₆ x 7 ¹ / ₂ x 5 ⁷ / ₁₆ inch 361 x 191 x 138 mm	4 x 3 = 12	12.7/5.8
UFMX-442	14 ³ / ₁₆ x 10 x 5 ⁷ / ₁₆ inch 361 x 254 x 138 mm	4 x 4 = 16	16.5/7.5
UFMX-452	14 ³ / ₁₆ x 12 ¹ / ₂ x 5 ⁷ / ₁₆ inch 361 x 317 x 138 mm	4 x 5 = 20	20.1/9.1
UFMX-532	17 ³ / ₄ x 7 ¹ / ₂ x 5 ⁷ / ₁₆ inch 451 x 191 x 138 mm	5 x 3 = 15	15.6/7.1
UFMX-542	17 ³ / ₄ x 10 x 5 ⁷ / ₁₆ inch 451 x 254 x 138 mm	5 x 4 = 20	20.2/9.2
UFMX-552	17 ³ / ₄ x 12 ¹ / ₂ x 5 ⁷ / ₁₆ inch 451 x 317 x 138 mm	5 x 5 = 25	24.8/11.3
New UFMX-562	17 ³ / ₄ x 15 x 5 ⁷ / ₁₆ inch 451 x 380 x 138 mm	5 x 6 = 30	29.5/13.4
UFMX-632	21 ¹ / ₄ x 7 ¹ / ₂ x 5 ⁷ / ₁₆ inch 540 x 191 x 138 mm	6 x 3 = 18	18.6/8.4
UFMX-642	21 ¹ / ₄ x 10 x 5 ⁷ / ₁₆ inch 540 x 254 x 138 mm	6 x 4 = 24	24.2/11.0
UFMX-652	21 ¹ / ₄ x 12 ¹ / ₂ x 5 ⁷ / ₁₆ inch 540 x 317 x 138 mm	6 x 5 = 30	29.5/13.5
New UFMX-662	21 ¹ / ₄ x 15 x 5 ⁷ / ₁₆ inch 540 x 380 x 138 mm	6 x 6 = 36	35.0/16.0