


Upright Freezer Rack



For Matrix

1.0ml Screw-Top Latch Rack (H=2¹/₈ inch / 54.8mm)

H: 2⁵/₁₆" - 2⁷/₁₆" (54.8-61.8mm)

For Micronic

Roborack-96 with High Cover (H=2¹/₄ inch / 56.1mm)

For Fluidx

96-Format Rack for Tube with Protective Jacket with Tall Profile Lid (H=2³/₁₆ inch / 55mm)


96-Format Rack Tack for 1.0ml and 1.4ml Tube with Tall Profile Lid (H=2⁷/₁₆ inch / 61.8mm)

24-Format Rack (H=2¹/₄ inch / 56mm)

For Abgene

1.2ml 2D Barcoded Tube in Stackable Rack (H=2¹/₄ inch / 56.13mm)

1.2ml 2D Barcoded Storage Tube Rack (H=2¹/₄ inch / 56.13mm)




- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.

Part Number	Dimension L x H x W	Capacity	Shipping Weight (lb/Kg)
UFMX-332	10 ¹¹ / ₁₆ x 7 ² / ₅ x 5 ⁷ / ₁₆ inch 271 x 191 x 138 mm	3 x 3 = 9	9.5/4.3
UFMX-342	10 ¹¹ / ₁₆ x 10 x 5 ⁷ / ₁₆ inch 271 x 254 x 138 mm	3 x 4 = 12	12.5/5.7
UFMX-352	10 ¹¹ / ₁₆ x 12 ¹ / ₂ x 5 ⁷ / ₁₆ inch 271 x 317 x 138 mm	3 x 5 = 15	15.3/6.9
UFMX-432	14 ³ / ₁₆ x 7 ¹ / ₂ x 5 ⁷ / ₁₆ inch 361 x 191 x 138 mm	4 x 3 = 12	12.7/5.8
UFMX-442	14 ³ / ₁₆ x 10 x 5 ⁷ / ₁₆ inch 361 x 254 x 138 mm	4 x 4 = 16	16.5/7.5
UFMX-452	14 ³ / ₁₆ x 12 ¹ / ₂ x 5 ⁷ / ₁₆ inch 361 x 317 x 138 mm	4 x 5 = 20	20.1/9.1
UFMX-532	17 ³ / ₄ x 7 ¹ / ₂ x 5 ⁷ / ₁₆ inch 451 x 191 x 138 mm	5 x 3 = 15	15.6/7.1
UFMX-542	17 ³ / ₄ x 10 x 5 ⁷ / ₁₆ inch 451 x 254 x 138 mm	5 x 4 = 20	20.2/9.2
UFMX-552	17 ³ / ₄ x 12 ¹ / ₂ x 5 ⁷ / ₁₆ inch 451 x 317 x 138 mm	5 x 5 = 25	24.8/11.3
New UFMX-562	17 ³ / ₄ x 15 x 5 ⁷ / ₁₆ inch 451 x 380 x 138 mm	5 x 6 = 30	29.5/13.4
UFMX-632	21 ¹ / ₄ x 7 ¹ / ₂ x 5 ⁷ / ₁₆ inch 540 x 191 x 138 mm	6 x 3 = 18	18.6/8.4
UFMX-642	21 ¹ / ₄ x 10 x 5 ⁷ / ₁₆ inch 540 x 254 x 138 mm	6 x 4 = 24	24.2/11.0
UFMX-652	21 ¹ / ₄ x 12 ¹ / ₂ x 5 ⁷ / ₁₆ inch 540 x 317 x 138 mm	6 x 5 = 30	29.5/13.5
New UFMX-662	21 ¹ / ₄ x 15 x 5 ⁷ / ₁₆ inch 540 x 380 x 138 mm	6 x 6 = 36	35.0/16.0

Upright Freezer Racks

For 96 Deep-Well Microtiter Plates



- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.
- ◆ Security lock devices are available for each rack.




Part Number	Dimension L x H x W	Capacity (Plate)	Lock Device Option*	Shipping Weight (lb/kg)
New UFDP-34L	11 x 7 ⁷ / ₈ x 5 ¹ / ₂ inch 280.5 x 187.1 x 140mm	3 x 4 = 12		DP-1107 8.5/3.9
New UFDP-35L	11 x 9 ³ / ₁₆ x 5 ¹ / ₂ inch 280.5 x 233.5 x 140mm	3 x 5 = 15		DP-1109 10.3/4.7
New UFDP-36L	11 x 11 x 5 ¹ / ₂ inch 280.5 x 279.5 x 140mm	3 x 6 = 18		DP-1111 12.5/5.5
New UFDP-44L	14 ¹ / ₈ x 7 ³ / ₈ x 5 ¹ / ₂ inch 373.3 x 187.1 x 140mm	4 x 4 = 16		DP-1507 11.1/5.1
New UFDP-45L	14 ¹¹ / ₁₆ x 9 ³ / ₁₆ x 5 ¹ / ₂ inch 373.3 x 233.5 x 140mm	4 x 5 = 20		DP-1509 13.6/6.2
New UFDP-46L	14 ¹ / ₈ x 11 x 5 ¹ / ₂ inch 373.3 x 279.5 x 140mm	4 x 6 = 24		DP-1511 16.0/7.3
New UFDP-54L	18 ³ / ₈ x 7 ³ / ₈ x 5 ⁷ / ₁₆ inch 466.1 x 187.1 x 140mm	5 x 4 = 20		DP-1807 13.7/6.2
UFDP-55L	18 ³ / ₈ x 9 ³ / ₁₆ x 5 ¹ / ₂ inch 466 x 233 x 140mm	5 x 5 = 25		DP-1809 16.8/7.6
UFDP-56L	18 ³ / ₈ x 11 x 5 ¹ / ₂ inch 466 x 280 x 140mm	5 x 6 = 30		DP-1811 19.8/9.0
New UFDP-64L	22 x 7 ⁷ / ₈ x 5 ¹ / ₂ inch 559 x 187.1 x 140mm	6 x 4 = 24		DP-2207 16.4/7.4
UFDP-65L	22 x 9 ³ / ₁₆ x 5 ¹ / ₂ inch 559 x 233 x 140mm	6 x 5 = 30		DP-2209 20.0/9.1
UFDP-66L	22 x 11 x 5 ¹ / ₂ inch 559 x 280 x 140mm	6 x 6 = 36		DP-2211 23.6/10.7
New UFDP-74L	25 ¹ / ₁₆ x 7 ³ / ₈ x 5 ¹ / ₂ inch 651.7 x 187.1 x 140mm	7 x 4 = 28		DP-2507 19.0/8.6
New UFDP-75L	25 ¹¹ / ₁₆ x 9 ³ / ₁₆ x 5 ¹ / ₂ inch 651.7 x 233.5 x 140mm	7 x 5 = 35		DP-2509 23.2/10.5
New UFDP-76L	25 ¹ / ₁₆ x 11 x 5 ¹ / ₂ inch 651.7 x 279.5 x 140mm	7 x 6 = 42		DP-2511 27.4/12.4


*Lock device is ordered separately.

Upright Freezer Drawer Racks

For 96 Deep-Well Microtiter Plates



- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer racks allow easy access to fewer microtiter plates at a time thus protecting other samples in the rack from warming up.
- ◆ Security lock devices are available for each rack.




Part Number	Description	Dimension L x H x W	Capacity	Lock Device Option*	Shipping Weight(lb/Kg)
UFDDP-33	Rack w/3 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	10 ³ / ₁₆ x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 277 x 154 x 138mm	3x3=9	DD-05	8.5/3.9
UFDDP-34	Rack w/4 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	10 ⁷ / ₁₆ x 8 ⁷ / ₁₆ x 5 ⁷ / ₁₆ inch 277 x 205 x 138mm	3x4=12	DD-07	11.1/5.0
UFDDP-35	Rack w/5 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	10 ⁹ / ₁₆ x 10 ¹ / ₈ x 5 ⁷ / ₁₆ inch 277 x 256 x 138mm	3x5=15	DR-10	13.7/6.2
New UFDDP-36	Rack w/6 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	10 ³ / ₁₆ x 12 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 277 x 307 x 138mm	3x6=18	DD-12	16.3/7.4
New UFDDP-37	Rack w/7 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	10 ⁵ / ₁₆ x 14 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 277 x 358 x 138mm	3x7=21	DD-14	19.2/8.7
UFDDP-43	Rack w/3 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	14 ¹ / ₂ x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 367 x 154 x 138mm	4x3=12	DD-05	10.9/5
UFDDP-44	Rack w/4 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	14 ¹ / ₂ x 8 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 367 x 205 x 138mm	4x4=16	DD-07	13.6/6.2
UFDDP-45	Rack w/5 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	14 ¹ / ₂ x 10 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 367 x 256 x 138mm	4x5=20	DR-10	17.6/8.0
New UFDDP-46	Rack w/6 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	14 ¹ / ₂ x 12 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 367 x 307 x 138mm	4x6=24	DD-12	20.9/9.5
New UFDDP-47	Rack w/7 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	14 ¹ / ₂ x 14 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 367 x 358 x 138mm	4x7=28	DD-14	27.0/12.6
UFDDP-53	Rack w/3 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	18 x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 457 x 154 x 138mm	5x3=15	DD-05	12.3/5.3
UFDDP-54	Rack w/4 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	18 x 8 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 457 x 205 x 138mm	5x4=20	DD-07	17.4/7.9
UFDDP-55	Rack w/5 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	18 x 10 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 457 x 256 x 138mm	5x5=25	DR-10	21.5/9.8
New UFDDP-56	Rack w/6 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	18 x 12 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 457 x 307 x 138mm	5x6=30	DD-12	25.5/11.6
New UFDDP-57	Rack w/7 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	18 x 14 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 457 x 358 x 138mm	5x7=35	DD-14	30.5/13.8
UFDDP-63	Rack w/3 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	21 ¹ / ₂ x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 547 x 154 x 138mm	6x3=18	DD-05	15.7/7.2
UFDDP-64	Rack w/4 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	21 ¹ / ₂ x 8 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 547 x 205 x 138mm	6x4=24	DD-07	20.5/9.3
UFDDP-65	Rack w/5 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	21 ¹ / ₂ x 10 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 547 x 256 x 138mm	6x5=30	DR-10	25.3/11.5
New UFDDP-66	Rack w/6 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	21 ¹ / ₂ x 12 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 547 x 307 x 138mm	6x6=36	DD-12	30.2/17.3
New UFDDP-67	Rack w/7 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	21 ¹ / ₂ x 14 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 547 x 358 x 138mm	6x7=42	DD-14	49.6/22.5
UFDDP-73	Rack w/3 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	25 x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 635 x 154 x 138mm	7x3=21	DD-05	12.1/8.2
UFDDP-74	Rack w/4 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	25 x 8 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 635 x 205 x 138mm	7x4=28	DD-07	23.7/10.8
UFDDP-75	Rack w/5 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	25 x 10 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 635 x 256 x 138mm	7x5=35	DR-10	29.2/13.3
New UFDDP-76	Rack w/6 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	25 x 12 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 635 x 307 x 138mm	7x6=42	DD-12	34.8/15.4
New UFDDP-77	Rack w/7 drawers (3 compartments per drawer) for 96Deep Well Microtiter Plates	25 x 14 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 635 x 358 x 138mm	7x7=49	DD-14	39.2/17.8

*Lock device is ordered separately.


Material : Stainless Steel

Upright Freezer Drawer Racks For MATRIX Boxes



L1 x W1 x H1: 5" x 3.4" x 2.3"
127.8 x 85.48 x 58mm

Listed racks are for 96-Well MATRIX racks and are constructed of corrosion resistant stainless steel, including handles and identification card holders. Each rack compartment holds one 96-Well MATRIX rack. Drawer style rack design allows easy access to fewer racks at a time thus protecting other samples in the rack from warming up. Security lock devices are also available for each rack. Using the internal dimensions obtained for your freezer shelf space (see page 1, fig 1) determine the height (H), length (L) and width (W) of each rack you select for your upright freezer. Racks fit all standard types of upright freezers.




Racks Only

Part Number	Dimension L x H x W		Capacity (MATRIX Racks)	Lock Device Option*
UFDMX-332	272.5x202x137	10 ⁷ / ₁₀ x 8x 5 ² / ₃ inch	9	DDMX-06
UFDMX-342	272.5x268.5x137	10 ⁷ / ₁₀ x10 ¹ / ₃ x 5 ² / ₃ inch	12	DDMX-09
UFDMX-352	272.5x335x137	10 ⁷ / ₁₀ x13 ¹ / ₃ x 5 ² / ₃ inch	15	DDMX-11
UFDMX-432	361x202x137	14 ¹ / ₃ x 8x 5 ² / ₃ inch	12	DDMX-06
UFDMX-442	361x268.5x137	14 ¹ / ₃ x10 ¹ / ₃ x 5 ² / ₃ inch	16	DDMX-09
UFDMX-452	361x335x137	14 ¹ / ₃ x13 ¹ / ₃ x 5 ² / ₃ inch	20	DDMX-11
UFDMX-532	450x202x137	17 ⁷ / ₁₀ x8x 5 ² / ₃ inch	15	DDMX-06
UFDMX-542	450x268.5x137	17 ⁷ / ₁₀ x10 ¹ / ₃ x 5 ² / ₃ inch	20	DDMX-09
UFDMX-552	450x335x137	17 ⁷ / ₁₀ x13 ¹ / ₃ x 5 ² / ₃ inch	25	DDMX-11
UFDMX-632	538.5x202x137	21 ¹ / ₃ x8x 5 ² / ₃ inch	18	DDMX-06
UFDMX-642	538.5x268.5x137	21 ¹ / ₃ x10 ¹ / ₃ x 5 ² / ₃ inch	24	DDMX-09
UFDMX-652	538.5x335x137	21 ¹ / ₃ x13 ¹ / ₃ x 5 ² / ₃ inch	30	DDMX-11


Material : Stainless Steel

Upright Freezer Drawer Racks

For Micronic LOBO Racks



Listed racks are for 96-Well Micronic LOBO racks and are constructed of corrosion resistant stainless steel, including handles and identification card holders. Each rack compartment holds one 96-Well Micronic LOBO rack. Drawer style rack design allows easy access to fewer racks at a time thus protecting other samples in the rack from warming up. Security lock devices are also available for each rack. Using the internal dimensions obtained for your freezer shelf space (see page 1, fig 1) determine the height (H), length (L) and width (W) of each rack you select for your upright freezer. Racks fit all standard types of upright freezers.



Part Number	Description	Dimension L x H x W	Capacity (LOBO Racks)	Lock Device Option*
UFDDP-33	Rack w/ 3 drawers (3 compartments per drawer) for Micronic LOBO Racks	10 ¹⁵ / ₁₆ x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 277 x 154 x 138mm	9	DD-05
UFDDP-34	Rack w/ 4 drawers (3 compartments per drawer) for Micronic LOBO Racks	10 ¹⁵ / ₁₆ x 8 ¹ / ₈ x 5 ⁷ / ₁₆ inch 277 x 205 x 138mm	12	DD-07
UFDDP-35	Rack w/ 5 drawers (3 compartments per drawer) for Micronic LOBO Racks	10 ¹⁵ / ₁₆ x 10 ¹ / ₄ x 5 ⁷ / ₁₆ inch 277 x 256 x 138mm	15	DD-10
UFDDP-43	Rack w/ 3 drawers (4 compartments per drawer) for Micronic LOBO Racks	14 ¹ / ₂ x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 367 x 154 x 138mm	12	DD-05
UFDDP-44	Rack w/ 4 drawers (4 compartments per drawer) for Micronic LOBO Racks	14 ¹ / ₂ x 8 ¹ / ₈ x 5 ⁷ / ₁₆ inch 367 x 205 x 138mm	16	DD-07
UFDDP-45	Rack w/ 5 drawers (4 compartments per drawer) for Micronic LOBO Racks	14 ¹ / ₂ x 10 ¹ / ₈ x 5 ⁷ / ₁₆ inch 367 x 256 x 138mm	20	DD-10
UFDDP-53	Rack w/ 3 drawers (5 compartments per drawer) for Micronic LOBO Racks	18 x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 457 x 154 x 138mm	15	DD-05
UFDDP-54	Rack w/ 4 drawers (5 compartments per drawer) for Micronic LOBO Racks	18 x 8 ¹ / ₈ x 5 ⁷ / ₁₆ inch 457 x 205 x 138mm	20	DD-07
UFDDP-55	Rack w/ 5 drawers (5 compartments per drawer) for Micronic LOBO Racks	18 x 10 ¹ / ₈ x 5 ⁷ / ₁₆ inch 457 x 256 x 138mm	25	DD-10
UFDDP-63	Rack w/ 3 drawers (6 compartments per drawer) for Micronic LOBO Racks	21 ¹ / ₂ x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 547 x 154 x 138mm	18	DD-05
UFDDP-64	Rack w/ 4 drawers (6 compartments per drawer) for Micronic LOBO Racks	21 ¹ / ₂ x 8 ¹ / ₈ x 5 ⁷ / ₁₆ inch 547 x 205 x 138mm	24	DD-07
UFDDP-65	Rack w/ 5 drawers (6 compartments per drawer) for Micronic LOBO Racks	21 ¹ / ₂ x 10 ¹ / ₈ x 5 ⁷ / ₁₆ inch 547 x 256 x 138mm	30	DD-10
UFDDP-73	Rack w/ 3 drawers (7 compartments per drawer) for Micronic LOBO Racks	25 x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 635 x 154 x 138mm	21	DD-05
UFDDP-74	Rack w/ 4 drawers (7 compartments per drawer) for Micronic LOBO Racks	25 x 8 ¹ / ₈ x 5 ⁷ / ₁₆ inch 635 x 205 x 138mm	28	DD-07
UFDDP-75	Rack w/ 5 drawers (7 compartments per drawer) for Micronic LOBO Racks	25 x 10 ¹ / ₈ x 5 ⁷ / ₁₆ inch 635 x 256 x 138mm	35	DD-10

*Lock device is ordered separately.