

# SUPER ADJUSTABLE 2™ SUPER ERECTA® SHELVING



## Shelving Starter and Add-On Units

Save time and money by eliminating adjacent posts and replacing with "S" hooks.

- Starter Units consist of four posts and indicated number of Super Adjustable 2™ Super Erecta shelves.
- Add-On Units consist of two posts, indicated number of shelves, and two "S" Hooks per shelf, for attaching shelf to adjacent unit.



**High Res**

**Replacement "S" Hook**  
2 are required for each storage level.

For additional wire shelves, order from page 11.


Add-on units with "S" hooks can only be attached to shelves supported by two or more posts.



**High Res**

**Security "S" Hook**  
2 are required for each storage level.

Chrome. Bag of 2  
Black. Bag of 2



Indicates antimicrobial product.

## Shelving Starter and Add-On Units — Easy Order Guide

| 5 Shelves<br>74" (1880mm) Posts | Chrome-Plated        |           |                  |           |                  |           | Metroseal 3 with Microban® |            |                  |            |                  |            |             |
|---------------------------------|----------------------|-----------|------------------|-----------|------------------|-----------|----------------------------|------------|------------------|------------|------------------|------------|-------------|
|                                 | 18" (457mm) Wide     |           | 21" (530mm) Wide |           | 24" (610mm) Wide |           | 18" (457mm) Wide           |            | 21" (530mm) Wide |            | 24" (610mm) Wide |            |             |
|                                 | Starter              | Add-On    | Starter          | Add-On    | Starter          | Add-On    | Starter                    | Add-On     | Starter          | Add-On     | Starter          | Add-On     |             |
| <b>24"</b><br>(610mm)           | <b>Cat. No.</b>      | LABC713A5 | LABC713AA5       | LABC714A5 | LABC714AA5       | LABC715A5 | LABC715AA5                 | LAB3K713A5 | LAB3K713AA5      | LAB3K714A5 | LAB3K714AA5      | LAB3K715A5 | LAB3K715AA5 |
| <b>Long</b>                     | <b>Wt. lbs. (kg)</b> | 51 (23)   | 43 (20)          | 57 (26)   | 48 (22)          | 62 (28)   | 53 (24)                    | 51 (23)    | 43 (20)          | 57 (26)    | 48 (22)          | 62 (28)    | 53 (24)     |
| <b>30"</b><br>(760mm)           | <b>Cat. No.</b>      | LABC723A5 | LABC723AA5       | LABC724A5 | LABC724AA5       | LABC725A5 | LABC725AA5                 | LAB3K723A5 | LAB3K723AA5      | LAB3K724A5 | LAB3K724AA5      | LAB3K725A5 | LAB3K725AA5 |
| <b>Long</b>                     | <b>Wt. lbs. (kg)</b> | 58 (26)   | 50 (23)          | 64 (29)   | 56 (25)          | 71 (32)   | 62 (28)                    | 58 (26)    | 50 (23)          | 64 (29)    | 56 (25)          | 71 (32)    | 62 (28)     |
| <b>36"</b><br>(914mm)           | <b>Cat. No.</b>      | LABC733A5 | LABC733AA5       | LABC734A5 | LABC734AA5       | LABC735A5 | LABC735AA5                 | LAB3K733A5 | LAB3K733AA5      | LAB3K734A5 | LAB3K734AA5      | LAB3K735A5 | LAB3K735AA5 |
| <b>Long</b>                     | <b>Wt. lbs. (kg)</b> | 65 (29)   | 57 (26)          | 73 (33)   | 65 (29)          | 81 (37)   | 73 (33)                    | 65 (29)    | 57 (26)          | 73 (33)    | 65 (29)          | 81 (37)    | 73 (33)     |
| <b>42"</b><br>(1066mm)          | <b>Cat. No.</b>      | LABC743A5 | LABC743AA5       | LABC744A5 | LABC744AA5       | LABC745A5 | LABC745AA5                 | LAB3K743A5 | LAB3K743AA5      | LAB3K744A5 | LAB3K744AA5      | LAB3K745A5 | LAB3K745AA5 |
| <b>Long</b>                     | <b>Wt. lbs. (kg)</b> | 70 (32)   | 62 (28)          | 81 (37)   | 73 (33)          | 89 (40)   | 81 (37)                    | 70 (32)    | 62 (28)          | 81 (37)    | 73 (33)          | 89 (40)    | 81 (37)     |
| <b>48"</b><br>(1219mm)          | <b>Cat. No.</b>      | LABC753A5 | LABC753AA5       | LABC754A5 | LABC754AA5       | LABC755A5 | LABC755AA5                 | LAB3K753A5 | LAB3K753AA5      | LAB3K754A5 | LAB3K754AA5      | LAB3K755A5 | LAB3K755AA5 |
| <b>Long</b>                     | <b>Wt. lbs. (kg)</b> | 78 (35)   | 70 (32)          | 88 (40)   | 80 (36)          | 98 (44)   | 90 (41)                    | 78 (35)    | 70 (32)          | 88 (40)    | 80 (36)          | 98 (44)    | 90 (41)     |
| <b>60"</b><br>(1524mm)          | <b>Cat. No.</b>      | LABC763A5 | LABC763AA5       | LABC764A5 | LABC764AA5       | LABC765A5 | LABC765AA5                 | LAB3K763A5 | LAB3K763AA5      | LAB3K764A5 | LAB3K764AA5      | LAB3K765A5 | LAB3K765AA5 |
| <b>Long</b>                     | <b>Wt. lbs. (kg)</b> | 102 (46)  | 93 (42)          | 109 (49)  | 101 (46)         | 122 (55)  | 113 (51)                   | 102 (46)   | 93 (42)          | 109 (49)   | 101 (46)         | 122 (55)   | 113 (51)    |
| <b>72"</b><br>(1829mm)          | <b>Cat. No.</b>      | LABC773A5 | LABC773AA5       | LABC774A5 | LABC774AA5       | LABC775A5 | LABC775AA5                 | LAB3K773A5 | LAB3K773AA5      | LAB3K774A5 | LAB3K774AA5      | LAB3K775A5 | LAB3K775AA5 |
| <b>Long</b>                     | <b>Wt. lbs. (kg)</b> | 117 (53)  | 108 (49)         | 122 (55)  | 113 (51)         | 149 (68)  | 141 (64)                   | 117 (53)   | 108 (49)         | 122 (55)   | 113 (51)         | 149 (68)   | 141 (64)    |