


Upright Freezer Drawer Racks

For Standard 2" Boxes

- ◆ Racks are for standard 2" cardboard (white or yellow), plastic (polypropylene or polycarbonate) or metal (stainless or aluminum) boxes.
- ◆ Drawer Racks offer a more organized system for storage and retrieval of samples.
- ◆ Constructed of corrosion resistant stainless steel.
- ◆ Security lock devices are available for each rack.



Rack Only



Rack w/Cardboard Box



Rack w/Polycarbonate Plastic Box



Rack w/Metal Box



w/DD Lock Device Option


* UFD-262 / UFD-362 / UFD-462 / UFD-562 (h:51.6mm)



Yellow 2" Cardboard Box



White 2" Cardboard Box



Polycarbonate Plastic Box





















Metal Box
S-Stainless Steel
A-Aluminum

[Check page 106 for Boxes]

Upright Freezer Drawer Racks

For Standard 2" Boxes

Part Number	Description	Dimension L x H x W	Capacity (Box)	Lock Device Option***	Shipping Weight lb(Kg)
New UFD-222	Rack Only	11 ¹ / ₈ x 4 ³ / ₄ x 5 ¹ / ₂ inch 283 x 121 x 140 mm	2 x 2 =4 	DD-04	4.0(1.8)
UFD-222-X'-XX''	Rack w/2" box and divider				4.0(1.8)+4*ω'
New UFD-232	Rack Only	11 ¹ / ₈ x 7 ¹ / ₈ x 5 ¹ / ₂ inch 283 x 181 x 140 mm	2 x 3 =6 	DD-06	5.7(2.6)
UFD-232-X'-XX''	Rack w/2" box and divider				5.7(2.6)+6*ω'
New UFD-242	Rack Only	11 ¹ / ₈ x 9 ⁷ / ₁₆ x 5 ¹ / ₂ inch 283 x 240 x 140 mm	2 x 4 =8 	DD-09	7.3(3.3)
UFD-242-X'-XX''	Rack w/2" box and divider				7.3(3.3)+8*ω'
UFD-252	Rack Only	11 ¹ / ₈ x 11 ¹³ / ₁₆ x 5 ¹ / ₂ inch 283 x 300 x 140 mm	2 x 5 =10 	DD-11	9.2(4.2)
UFD-252-X'-XX''	Rack w/2" box and divider				9.2(4.2)+10*ω'
UFD-262	Rack Only	11 ¹ / ₈ x 14 ¹ / ₈ x 5 ¹ / ₂ inch 283 x 320 x 140 mm	2 x 6 =12 	DD-13	8.4(3.8)
UFD-262-X'-XX''	Rack w/2" box and divider				8.4(3.8)+12*ω'
UFD-332	Rack Only	16 ¹ / ₂ x 7 ¹ / ₈ x 5 ¹ / ₂ inch 419 x 181 x 140 mm	3 x 3 =9 	DD-06	7.5(3.4)
UFD-332-X'-XX''	Rack w/2" box and divider				7.5(3.4)+9*ω'
UFD-342	Rack Only	16 ¹ / ₂ x 9 ⁷ / ₁₆ x 5 ¹ / ₂ inch 419 x 240 x 140 mm	3 x 4 =12 	DD-09	9.5(4.3)
UFD-342-X'-XX''	Rack w/2" box and divider				9.5(4.3)+12*ω'
UFD-352	Rack Only	16 ¹ / ₂ x 11 ¹³ / ₁₆ x 5 ¹ / ₂ inch 419 x 300 x 140 mm	3 x 5 =15 	DD-11	11.4(5.2)
UFD-352-X'-XX''	Rack w/2" box and divider				11.4(5.2)+15*ω'
New UFD-362	Rack Only	16 ¹ / ₂ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 419 x 320 x 140 mm	3 x 6 =18 	DD-13	11.3(5.1)
UFD-362-X'-XX''	Rack w/2" box and divider				11.3(5.1)+18*ω'
UFD-432	Rack Only	22 x 7 ¹ / ₈ x 5 ¹ / ₂ inch 559 x 181 x 140 mm	4 x 3 =12 	DD-06	9.7(4.4)
UFD-432-X'-XX''	Rack w/2" box and divider				9.7(4.4)+12*ω'
UFD-442	Rack Only	22 x 9 ⁷ / ₁₆ x 5 ¹ / ₂ inch 559 x 240 x 140 mm	4 x 4 =16 	DD-09	12.3(5.6)
UFD-442-X'-XX''	Rack w/2" box and divider				12.3(5.6)+16*ω'
UFD-452	Rack Only	22 x 11 ¹³ / ₁₆ x 5 ¹ / ₂ inch 559 x 300 x 140 mm	4 x 5 =20 	DD-11	14.8(6.7)
UFD-452-X'-XX''	Rack w/2" box and divider				14.8(6.7)+20*ω'
UFD-462	Rack Only	22 x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 559 x 320 x 140 mm	4 x 6 =24 	DD-13	14.9(6.8)
UFD-462-X'-XX''	Rack w/2" box and divider				14.9(6.8)+24*ω'
UFD-472	Rack Only	22 x 16 ¹ / ₂ x 5 ¹ / ₂ inch 559 x 419 x 140 mm	4 x 7 =28 	DD-16	19.9(9.0)
UFD-472-X'-XX''	Rack w/2" box and divider				19.9(9.0)+28*ω'
UFD-532	Rack Only	26 ³ / ₄ x 7 ¹ / ₈ x 5 ¹ / ₂ inch 679 x 181 x 140 mm	5 x 3 =15 	DD-06	11.3(5.1)
UFD-532-X'-XX''	Rack w/2" box and divider				11.3(5.1)+15*ω'
UFD-542	Rack Only	26 ³ / ₄ x 9 ⁷ / ₁₆ x 5 ¹ / ₂ inch 679 x 240 x 140 mm	5 x 4 =20 	DD-09	14.5(6.6)
UFD-542-X'-XX''	Rack w/2" box and divider				14.5(6.6)+20*ω'
UFD-552	Rack Only	26 ³ / ₄ x 11 ¹³ / ₁₆ x 5 ¹ / ₂ inch 679 x 300 x 140 mm	5 x 5 =25 	DD-11	17.6(8.0)
UFD-552-X'-XX''	Rack w/2" box and divider				17.6(8.0)+25*ω'
New UFD-562	Rack Only	26 ³ / ₄ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 679 x 320 x 140 mm	5 x 6 =30 	DD-13	16.8(7.6)
UFD-562-X'-XX''	Rack w/2" box and divider				16.8(7.6)+30*ω'

* To order, replace X with box type: Y-Yellow cardboard box; W-White cardboard box; P-Plastic Polycarbonate box; PP-Plastic Polypropylene Box; S-Stainless Steel Box or A-Aluminum Box.

** To order, replace XX with 16,25,36,49,64,81,100 or 144 cell divider. Plastic Polycarbonate and Polypropylene Boxes only have 81 and 100 cell dividers.

*** Lock device is ordered separately.

****ω'=0.22lb(0.1kg) for cardboard box, 0.37lb(0.17kg) for plastic box, 0.84lb(0.38kg) for stainless steel box, 0.53lb(0.24kg) for aluminum box.