


Upright Freezer Drawer Rack

New



H: 1 1/4" - 1 1/2" (31.7-38.7mm)

For Micronic

Loborack-96 with Standard Cover (H=1 1/4 inch / 31.7mm)

Loborack-96 with High Cover (H=1 1/2 inch / 38.7mm)

For Fluidx

Low Profile 96-Format Rack for all 0.5ml and 0.75ml Tube with Standard Profile Lid (H=1 5/16 inch / 32.9mm)


96-Format Rack for Tube with Protective Jacket with Low Cover (H=1 1/2 inch / 38mm)




For Abgene

0.5ml 2D Barcoded Lock Storage Tube Rack (H=1 1/2 inch / 38.5mm)

0.6ml 2D Barcoded Storage Tube Rack (H=1 1/2 inch / 38.5mm)

- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer racks allow easy access to and storage of samples without having to remove the whole rack from the freezer.



Part Number	Dimension L x H x W	Capacity	Shipping Weight(lb/kg)
UFDML-43	14 ³ / ₁₆ x 5 ⁹ / ₁₆ x 5 ⁷ / ₁₆ inch 361 x 142 x 138 mm	4 x 3 = 12 	11.9/5.4
UFDML-44	14 ³ / ₁₆ x 7 ⁷ / ₁₆ x 5 ⁷ / ₁₆ inch 361 x 188 x 138 mm	4 x 4 = 16 	15.6/7.1
UFDML-45	14 ³ / ₁₆ x 9 ¹ / ₄ x 5 ⁷ / ₁₆ inch 361 x 235 x 138 mm	4 x 5 = 20 	19.2/8.7
UFDML-46	14 ³ / ₁₆ x 11 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 361 x 281 x 138 mm	4 x 6 = 24 	22.9/10.4
UFDML-53	17 ¹¹ / ₁₆ x 5 ⁹ / ₁₆ x 5 ⁷ / ₁₆ inch 449 x 142 x 138 mm	5 x 3 = 15 	14.6/6.6
UFDML-54	17 ¹¹ / ₁₆ x 7 ⁷ / ₁₆ x 5 ⁷ / ₁₆ inch 449 x 188 x 138 mm	5 x 4 = 20 	19.1/8.7
UFDML-55	17 ¹¹ / ₁₆ x 9 ¹ / ₄ x 5 ⁷ / ₁₆ inch 449 x 235 x 138 mm	5 x 5 = 25 	23.6/10.7
UFDML-56	17 ¹¹ / ₁₆ x 11 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 449 x 281 x 138 mm	5 x 6 = 30 	28.1/12.8
UFDML-63	21 ³ / ₁₆ x 5 ⁹ / ₁₆ x 5 ⁷ / ₁₆ inch 538 x 142 x 138 mm	6 x 3 = 18 	17.3/7.9
UFDML-64	21 ³ / ₁₆ x 7 ⁷ / ₁₆ x 5 ⁷ / ₁₆ inch 538 x 188 x 138 mm	6 x 4 = 24 	22.7/10.3
UFDML-65	21 ³ / ₁₆ x 9 ¹ / ₄ x 5 ⁷ / ₁₆ inch 538 x 235 x 138 mm	6 x 5 = 30 	28.0/12.7
UFDML-66	21 ³ / ₁₆ x 11 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 538 x 281 x 138 mm	6 x 6 = 36 	33.3/15.2