



Upright Freezer Drawer Racks

For 100-Cell Hinged Top Plastic Storage Boxes




[Check page 107 for Boxes]

- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer racks allow easy access to and storage of samples without having to remove the whole rack from the freezer.
- ◆ Security lock devices are available for each rack.



Rack Only



Rack w/100-Cell Hinged Plastic Storage Box



Rack w/DDHT Lock Device Option

Part Number	Description	Dimension L x H x W	Capacity (Box)	Lock Device Option**	Shipping Weight lb(Kg)
UFDHT-23	Rack Only	11 ⁵ / ₈ x 7 ³ / ₁₆ x 6 ¹ / ₂ inch	6	DDHT-06	6.2(2.8)
UFDHT-23-PX*	Rack w/100-Cell Hinged Plastic Storage Box	294 x 186 x 164 mm			6.2(2.8)+6*ω ¹
UFDHT-24	Rack Only	11 ³ / ₈ x 9 ¹¹ / ₁₆ x 6 ¹ / ₂ inch	8	DDHT-09	8.0(3.6)
UFDHT-24-PX*	Rack w/100-Cell Hinged Plastic Storage Box	294 x 247 x 164 mm			8.0(3.6)+8*ω ¹
UFDHT-25	Rack Only	11 ⁵ / ₈ x 12 ¹ / ₈ x 6 ¹ / ₂ inch	10	DDHT-11	9.9(4.5)
UFDHT-25-PX*	Rack w/100-Cell Hinged Plastic Storage Box	294 x 308 x 164 mm			9.9(4.5)+10*ω ¹
New UFDHT-26	Rack Only	11 ⁵ / ₈ x 14 ¹ / ₂ x 6 ¹ / ₂ inch	12	DDHT-14	11.7(5.3)
UFDHT-26-PX*	Rack w/100-Cell Hinged Plastic Storage Box	294 x 369 x 164 mm			11.7(5.3)+12*ω ¹
UFDHT-33	Rack Only	17 ¹ / ₄ x 7 ³ / ₁₆ x 6 ¹ / ₂ inch	9	DDHT-06	8.6(3.9)
UFDHT-33-PX*	Rack w/100-Cell Hinged Plastic Storage Box	437 x 186 x 164 mm			8.6(3.9)+9*ω ¹
UFDHT-34	Rack Only	17 ¹ / ₄ x 9 ¹¹ / ₁₆ x 6 ¹ / ₂ inch	12	DDHT-09	11(5.0)
UFDHT-34-PX*	Rack w/100-Cell Hinged Plastic Storage Box	437 x 247 x 164 mm			11(5.0)+12*ω ¹
UFDHT-35	Rack Only	17 ¹ / ₄ x 12 ¹ / ₈ x 6 ¹ / ₂ inch	15	DDHT-11	13.5(6.1)
UFDHT-35-PX*	Rack w/100-Cell Hinged Plastic Storage Box	437 x 308 x 164 mm			13.5(6.1)+15*ω ¹
New UFDHT-36	Rack Only	17 ¹ / ₄ x 14 ¹ / ₂ x 6 ¹ / ₂ inch	18	DDHT-14	16.1(7.3)
UFDHT-36-PX*	Rack w/100-Cell Hinged Plastic Storage Box	437 x 369 x 164 mm			16.1(7.3)+18*ω ¹
UFDHT-43	Rack Only	22 ⁷ / ₈ x 7 ³ / ₁₆ x 6 ¹ / ₂ inch	12	DDHT-06	11(5.0)
UFDHT-43-PX*	Rack w/100-Cell Hinged Plastic Storage Box	580 x 186 x 164 mm			11(5.0)+12*ω ¹
UFDHT-44	Rack Only	22 ⁷ / ₈ x 9 ¹¹ / ₁₆ x 6 ¹ / ₂ inch	16	DDHT-09	14.1(6.4)
UFDHT-44-PX*	Rack w/100-Cell Hinged Plastic Storage Box	580 x 247 x 164 mm			14.1(6.4)+16*ω ¹
New UFDHT-45	Rack Only	22 ⁷ / ₈ x 12 ¹ / ₈ x 6 ¹ / ₂ inch	20	DDHT-11	17(7.7)
UFDHT-45-PX*	Rack w/100-Cell Hinged Plastic Storage Box	580 x 308 x 164 mm			17(7.7)+20*ω ¹
New UFDHT-46	Rack Only	22 ⁷ / ₈ x 14 ¹ / ₂ x 6 ¹ / ₂ inch	24	DDHT-14	20.5(9.3)
UFDHT-46-PX*	Rack w/100-Cell Hinged Plastic Storage Box	580 x 369 x 164 mm			20.5(9.3)+24*ω ¹

* To order replace, X with color of box : O=orange, P=pink, Y=yellow, B=blue, G=green or N=natural.

** ω¹=0.221lb(0.1kg) for cardboard box, 0.371lb(0.17kg) for plastic box, 0.841lb(0.38kg) for stainless steel box, 0.531lb(0.24kg) for aluminum box.