



# Upright Freezer Racks

## For 100-Cell Hinged Top Plastic Storage Boxes




[ Check page 107 for Boxes ]


- ◆ Made of corrosion resistant stainless steel.
- ◆ Security lock devices are available for each rack.



Rack Only



Rack w/100-Cell Hinged Plastic Storage Box



Rack w/DP Lock Device Option

New

New

New

New

New

New

Part Number	Description	Dimension L x H x W	Capacity (Box)	Lock Device Option**	Shipping Weight lb(Kg)
UFHT-23	Rack Only	12 <sup>1</sup> / <sub>2</sub> x 6 <sup>15</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch	2 x 3 =6	DP-1307	4.0(1.8)
UFHT-23-PX*	Rack w/100-Cell Hinged Plastic Storage Box	318 x 176 x 149 mm			4.0(1.8)+6*ω <sup>1</sup>
UFHT-24	Rack Only	12 <sup>1</sup> / <sub>2</sub> x 9 <sup>1</sup> / <sub>4</sub> x 5 <sup>7</sup> / <sub>8</sub> inch	2 x 4 =8	DP-1309	5.1(2.3)
UFHT-24-PX*	Rack w/100-Cell Hinged Plastic Storage Box	318 x 235 x 149 mm			5.1(2.3)+8*ω <sup>1</sup>
UFHT-25	Rack Only	12 <sup>1</sup> / <sub>2</sub> x 11 <sup>9</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch	2 x 5 =10	DP-1312	5.9(2.7)
UFHT-25-PX*	Rack w/100-Cell Hinged Plastic Storage Box	318 x 293 x 149 mm			5.9(2.7)+10*ω <sup>1</sup>
UFHT-26	Rack Only	12 <sup>1</sup> / <sub>2</sub> x 13 <sup>13</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch	2 x 6 =12	DP-1313	6.8(3.1)
UFHT-26-PX*	Rack w/100-Cell Hinged Plastic Storage Box	318 x 351 x 149 mm			6.8(3.1)+12*ω <sup>1</sup>
UFHT-33	Rack Only	18 <sup>3</sup> / <sub>4</sub> x 6 <sup>15</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch	3 x 3 =9	DP-1907	5.7(2.6)
UFHT-33-PX*	Rack w/100-Cell Hinged Plastic Storage Box	476 x 176 x 149 mm			5.7(2.6)+9*ω <sup>1</sup>
UFHT-34	Rack Only	18 <sup>3</sup> / <sub>4</sub> x 9 <sup>1</sup> / <sub>4</sub> x 5 <sup>7</sup> / <sub>8</sub> inch	3 x 4 =12	DP-1909	7.0(3.2)
UFHT-34-PX*	Rack w/100-Cell Hinged Plastic Storage Box	476 x 235 x 149 mm			7.0(3.2)+12*ω <sup>1</sup>
UFHT-35	Rack Only	18 <sup>3</sup> / <sub>4</sub> x 11 <sup>9</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch	3 x 5 =15	DP-1912	8.4(3.8)
UFHT-35-PX*	Rack w/100-Cell Hinged Plastic Storage Box	476 x 293 x 149 mm			8.4(3.8)+15*ω <sup>1</sup>
UFHT-36	Rack Only	18 <sup>3</sup> / <sub>4</sub> x 13 <sup>13</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch	3 x 6 =18	DP-1913	9.7(4.4)
UFHT-36-PX*	Rack w/100-Cell Hinged Plastic Storage Box	476 x 351 x 149 mm			9.7(4.4)+18*ω <sup>1</sup>
UFHT-43	Rack Only	24 <sup>15</sup> / <sub>16</sub> x 6 <sup>15</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch	4 x 3 =12	DP-2407	7.3(3.3)
UFHT-43-PX*	Rack w/100-Cell Hinged Plastic Storage Box	634 x 177 x 149 mm			7.3(3.3)+12*ω <sup>1</sup>
UFHT-44	Rack Only	24 <sup>15</sup> / <sub>16</sub> x 9 <sup>1</sup> / <sub>4</sub> x 5 <sup>7</sup> / <sub>8</sub> inch	4 x 4 =16	DP-2409	9.2(4.2)
UFHT-44-PX*	Rack w/100-Cell Hinged Plastic Storage Box	634 x 235 x 149 mm			9.2(4.2)+15*ω <sup>1</sup>
UFHT-45	Rack Only	24 <sup>15</sup> / <sub>16</sub> x 11 <sup>9</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch	4 x 5 =20	DP-2412	10.6(4.8)
UFHT-45-PX*	Rack w/100-Cell Hinged Plastic Storage Box	634 x 293 x 149 mm			10.6(4.8)+20*ω <sup>1</sup>
UFHT-46	Rack Only	24 <sup>15</sup> / <sub>16</sub> x 13 <sup>13</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch	4 x 6 =24	DP-2413	12.1(5.5)
UFHT-46-PX*	Rack w/100-Cell Hinged Plastic Storage Box	634 x 351 x 149 mm			12.1(5.5)+24*ω <sup>1</sup>

\* To order replace , X with color of box : O-orange, P-pink, Y-yellow, B-blue, G-green or N-natural.

\*\* ω<sup>1</sup>= 0.41lb(0.18kg) for plastic box.