


Upright Freezer Racks

For 96 Deep-Well Microtiter Plates



- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.
- ◆ Security lock devices are available for each rack.



Part Number	Dimension L x H x W	Capacity (Plate)	Lock Device Option*	Shipping Weight (lb/kg)
New UFDP-34L	11 x 7 ⁷ / ₈ x 5 ¹ / ₂ inch 280.5 x 187.1 x 140mm	3 x 4 = 12	 DP-1107	8.5/3.9
New UFDP-35L	11 x 9 ³ / ₁₆ x 5 ¹ / ₂ inch 280.5 x 233.5 x 140mm	3 x 5 = 15	 DP-1109	10.3/4.7
New UFDP-36L	11 x 11 x 5 ¹ / ₂ inch 280.5 x 279.5 x 140mm	3 x 6 = 18	 DP-1111	12.5/5.5
New UFDP-44L	14 ¹ / ₈ x 7 ³ / ₈ x 5 ¹ / ₂ inch 373.3 x 187.1 x 140mm	4 x 4 = 16	 DP-1507	11.1/5.1
New UFDP-45L	14 ¹ / ₈ x 9 ³ / ₁₆ x 5 ¹ / ₂ inch 373.3 x 233.5 x 140mm	4 x 5 = 20	 DP-1509	13.6/6.2
New UFDP-46L	14 ¹ / ₈ x 11 x 5 ¹ / ₂ inch 373.3 x 279.5 x 140mm	4 x 6 = 24	 DP-1511	16.0/7.3
New UFDP-54L	18 ³ / ₈ x 7 ³ / ₈ x 5 ⁷ / ₁₆ inch 466.1 x 187.1 x 140mm	5 x 4 = 20	 DP-1807	13.7/6.2
UFDP-55L	18 ³ / ₈ x 9 ³ / ₁₆ x 5 ¹ / ₂ inch 466 x 233 x 140mm	5 x 5 = 25	 DP-1809	16.8/7.6
UFDP-56L	18 ³ / ₈ x 11 x 5 ¹ / ₂ inch 466 x 280 x 140mm	5 x 6 = 30	 DP-1811	19.8/9.0
New UFDP-64L	22 x 7 ⁷ / ₈ x 5 ¹ / ₂ inch 559 x 187.1 x 140mm	6 x 4 = 24	 DP-2207	16.4/7.4
UFDP-65L	22 x 9 ³ / ₁₆ x 5 ¹ / ₂ inch 559 x 233 x 140mm	6 x 5 = 30	 DP-2209	20.0/9.1
UFDP-66L	22 x 11 x 5 ¹ / ₂ inch 559 x 280 x 140mm	6 x 6 = 36	 DP-2211	23.6/10.7
New UFDP-74L	25 ¹ / ₁₆ x 7 ³ / ₈ x 5 ¹ / ₂ inch 651.7 x 187.1 x 140mm	7 x 4 = 28	 DP-2507	19.0/8.6
New UFDP-75L	25 ¹ / ₁₆ x 9 ³ / ₁₆ x 5 ¹ / ₂ inch 651.7 x 233.5 x 140mm	7 x 5 = 35	 DP-2509	23.2/10.5
New UFDP-76L	25 ¹ / ₁₆ x 11 x 5 ¹ / ₂ inch 651.7 x 279.5 x 140mm	7 x 6 = 42	 DP-2511	27.4/12.4

*Lock device is ordered separately.