



Upright Freezer Racks

For 96-Well Microtube Boxes



- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.
- ◆ Security lock devices are available for each rack.



Part Number	Dimension L x H x W	Capacity (Box)	Lock Device Option*	Shipping Weight(lb/Kg)
New UFTB-33L	11 ¹ / ₁₆ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 281 x 169 x 140mm	3 x 3 = 9	DP-1107	8.9/4.0
New UFTB-34L	11 ¹ / ₁₆ x 8 ¹³ / ₁₆ x 5 ¹ / ₂ inch 281 x 224 x 140mm	3 x 4 = 12	DP-1109	11.5/5.2
New UFTB-35L	11 ¹ / ₁₆ x 11 x 5 ¹ / ₂ inch 281 x 280 x 140mm	3 x 5 = 15	DP-1111	14.1/6.4
New UFTB-36L	11 ¹ / ₁₆ x 13 ³ / ₁₆ x 5 ¹ / ₂ inch 281 x 335 x 140mm	3 x 6 = 18	DP-1113	16.7/7.6
New UFTB-43L	14 ¹¹ / ₁₆ x 6 ⁷ / ₈ x 5 ¹ / ₂ inch 373 x 169 x 140mm	4 x 3 = 12	DP-1407	11.7/5.3
New UFTB-44L	14 ¹¹ / ₁₆ x 8 ¹³ / ₁₆ x 5 ¹ / ₂ inch 373 x 224 x 140mm	4 x 4 = 16	DP-1409	15.1/6.9
New UFTB-45L	14 ¹¹ / ₁₆ x 11 x 5 ¹ / ₂ inch 373 x 280 x 140mm	4 x 5 = 20	DP-1411	18.5/8.4
New UFTB-46L	14 ¹¹ / ₁₆ x 13 ³ / ₁₆ x 5 ¹ / ₂ inch 373 x 335 x 140mm	4 x 6 = 24	DP-1413	21.9/10.0
New UFTB-53L	18 ³ / ₈ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 466 x 169 x 140mm	5 x 3 = 15	DP-1807	14.4/6.6
UFTB-54L	18 ³ / ₈ x 8 ¹³ / ₁₆ x 5 ¹ / ₂ inch 466 x 224 x 140mm	5 x 4 = 20	DP-1809	18.7/8.5
UFTB-55L	18 ³ / ₈ x 11 x 5 ¹ / ₂ inch 466 x 280 x 140mm	5 x 5 = 25	DP-1811	22.9/10.4
New UFTB-56L	18 ³ / ₈ x 13 ³ / ₁₆ x 5 ¹ / ₂ inch 466 x 335 x 140mm	5 x 6 = 30	DP-1813	17.2/7.8
UFTB-64L	22 x 8 ¹³ / ₁₆ x 5 ¹ / ₂ inch 559 x 224 x 140mm	6 x 4 = 24	DP-2209	22.3/10.1
UFTB-65L	22 x 11 x 5 ¹ / ₂ inch 559 x 280 x 140mm	6 x 5 = 30	DP-2211	27.4/12.4
New UFTB-66L	22 x 13 ³ / ₁₆ x 5 ¹ / ₂ inch 559 x 335 x 140mm	6 x 6 = 36	DP-2213	32.4/14.7
New UFTB-73L	25 ¹¹ / ₁₆ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 652 x 169 x 140mm	7 x 3 = 21	DP-2607	20.0/9.1
New UFTB-74L	25 ¹¹ / ₁₆ x 8 ¹³ / ₁₆ x 5 ¹ / ₂ inch 652 x 224 x 140mm	7 x 4 = 28	DP-2609	25.3/11.5
New UFTB-75L	25 ¹¹ / ₁₆ x 11 x 5 ¹ / ₂ inch 652 x 280 x 140mm	7 x 5 = 35	DP-2611	31.8/14.5
New UFTB-76L	25 ¹¹ / ₁₆ x 13 ³ / ₁₆ x 5 ¹ / ₂ inch 652 x 335 x 140mm	7 x 6 = 42	DP-2613	37.7/17.1

*Lock device is ordered separately.