



Material : Stainless Steel

# Upright Freezer Racks

## For 96-Well & 384-Well Microtiter Plates



Racks are for 96 well and 384 well Microtiter Plates. High quality Upright Freezer Racks, including handles, are constructed of corrosion resistant stainless steel and offer convenient storage of samples in an organized format. Security lock devices are also available for each rack. Using the internal dimensions obtained for your freezer shelf space (see page 1, fig 1) determine the height (H), length (L) and width (W) of each rack you select for your upright freezer. Racks come with locking rods to secure microtiter plates in place.











Rack w/o Locking Rod



Rack w/ Locking Rod



Rack w/DP Lock Device Option

Part Number	Dimension L x H x W	Capacity (Plate)	Lock Device Option*
UFMP-409L	14 <sup>7</sup> / <sub>10</sub> x 9 <sup>3</sup> / <sub>16</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 374 x 233 x 140 mm	56 Plates w/lids 64 Plates w/o lids	 DP-1409
UFMP-411L	14 <sup>7</sup> / <sub>10</sub> x 11 x 5 <sup>1</sup> / <sub>2</sub> inch 374 x 280 x 140 mm	64 Plates w/lids 80 Plates w/o lids	 DP-1411
UFMP-509L	18 <sup>3</sup> / <sub>8</sub> x 9 <sup>3</sup> / <sub>16</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 466 x 233 x 140 mm	70 Plates w/lids 80 Plates w/o lids	 DP-1809
UFMP-511L	18 <sup>3</sup> / <sub>8</sub> x 11 x 5 <sup>1</sup> / <sub>2</sub> inch 466 x 280 x 140 mm	80 Plates w/lids 100 Plates w/o lids	 DP-1911
UFMP-609L	22 x 9 <sup>3</sup> / <sub>16</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 559 x 233 x 140 mm	84 Plates w/lids 96 Plates w/o lids	 DP-2209
UFMP-611L	22 x 11 x 5 <sup>1</sup> / <sub>2</sub> inch 559 x 280 x 140 mm	96 Plates w/lids 120 Plates w/o lids	 DP-2211
UFMP-709L	25 <sup>11</sup> / <sub>16</sub> x 9 <sup>3</sup> / <sub>16</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 652 x 233 x 140 mm	98 Plates w/lids 112 Plates w/o lids	 DP-2609
UFMP-711L	25 <sup>11</sup> / <sub>16</sub> x 11 x 5 <sup>1</sup> / <sub>2</sub> inch 652 x 280 x 140 mm	112 Plates w/lids 140 Plates w/o lids	 DP-2611

\*Lock device is ordered separately.