


# Upright Freezer Rack

**New**



## For Micronic

Loborack-96 with Standard Cover (H=1<sup>1</sup>/<sub>4</sub>inch / 31.7mm)

Loborack-96 with High Cover (H=1<sup>1</sup>/<sub>2</sub>inch / 38.7mm)

## For Fluidx

Low Profile 96-Format Rack for all 0.5ml and 0.75ml Tube with Standard Profile Lid (H=1<sup>5</sup>/<sub>16</sub>inch / 32.9mm)


96-Format Rack for Tube with Protective Jacket with Low Cover (H=1<sup>1</sup>/<sub>2</sub>inch / 38mm)

## For Abgene

0.5ml 2D Barcoded Lock Storage Tube Rack (H=1<sup>1</sup>/<sub>2</sub>inch / 38.5mm)

0.6ml 2D Barcoded Storage Tube Rack (H=1<sup>1</sup>/<sub>2</sub>inch / 38.5mm)

- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.



Part Number	Dimension L x H x W	Capacity	Shipping Weight(μ/κg)
UFML-43L	14 <sup>3</sup> / <sub>16</sub> x 4 <sup>7</sup> / <sub>8</sub> x 5 <sup>7</sup> / <sub>16</sub> inch 361 x 124 x 138 mm	4 x 3 = 12	9.7/4.4
UFML-44L	14 <sup>7</sup> / <sub>16</sub> x 6 <sup>7</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>16</sub> inch 361 x 164 x 138 mm	4 x 4 = 16	12.4/5.7
UFML-45L	14 <sup>7</sup> / <sub>16</sub> x 8 <sup>1</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>16</sub> inch 361 x 205 x 138 mm	4 x 5 = 20	15.2/6.9
UFML-46L	14 <sup>3</sup> / <sub>16</sub> x 9 <sup>7</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>16</sub> inch 361 x 245 x 138 mm	4 x 6 = 24	18.0/8.2
UFML-53L	17 <sup>1</sup> / <sub>4</sub> x 4 <sup>7</sup> / <sub>8</sub> x 5 <sup>7</sup> / <sub>16</sub> inch 451 x 124 x 138 mm	5 x 3 = 15	12.0/5.4
UFML-54L	17 <sup>1</sup> / <sub>4</sub> x 6 <sup>7</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>16</sub> inch 451 x 164 x 138 mm	5 x 4 = 20	15.4/7.0
UFML-55L	17 <sup>1</sup> / <sub>4</sub> x 8 <sup>1</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>16</sub> inch 451 x 205 x 138 mm	5 x 5 = 25	18.9/8.6
UFML-56L	17 <sup>1</sup> / <sub>4</sub> x 9 <sup>1</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>16</sub> inch 451 x 245 x 138 mm	5 x 6 = 30	22.4/10.2
UFML-63L	21 <sup>1</sup> / <sub>4</sub> x 4 <sup>7</sup> / <sub>8</sub> x 5 <sup>7</sup> / <sub>16</sub> inch 540 x 124 x 138 mm	6 x 3 = 18	14.3/6.5
UFML-64L	21 <sup>1</sup> / <sub>4</sub> x 6 <sup>7</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>16</sub> inch 540 x 164 x 138 mm	6 x 4 = 24	18.4/8.4
UFML-65L	21 <sup>1</sup> / <sub>4</sub> x 8 <sup>1</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>16</sub> inch 540 x 205 x 138 mm	6 x 5 = 30	22.5/10.2
UFML-66L	21 <sup>1</sup> / <sub>4</sub> x 9 <sup>7</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>16</sub> inch 540 x 245 x 138 mm	6 x 6 = 36	26.7/12.1