



MW200 Series

MW Series Utility Carts — 12.01

- Consist of Super Erecta wire shelves with plastic split sleeves, two one-piece handles of matching finish, and designated casters with donut bumpers.
- Two- and three-tier models available.
- 375 lbs. (170kg) weight capacity per cart. MW carts use light-duty resilient casters.
- 18" (457mm) wide carts have 4" (102mm) casters. Cart is 38" (965mm) high.
- 21" and 24" (530 and 610mm) wide carts have 5" (127mm) casters. Cart is 39" (990mm) high.

Utility Carts with Stainless Steel Solid Shelves

Width/Length (in.) (mm)		Approx. Pkd. Wt. (lbs.) (kg)		Cat. No. 2-Shelf	List Price Each	Approx. Pkd. Wt. (lbs.) (kg)		Cat. No. 3-Shelf
18x24	457x610	40	18.1	LAB301WM	714.00	48	21.8	LAB302WM
18x30	457x760	45	20.4	LAB401WM	807.50	53	24.0	LAB402WM
18x36	457x914	48	21.8	LAB501WM	807.50	59	26.8	LAB502WM
21x36	530x914	54	24.5	LAB601WM	873.00	74	33.6	LAB602WM
24x36	610x914	60	27.2	LAB801WM	923.00	66	29.9	LAB802WM



MW400 Series

Utility Carts with 1 Solid, 2 Wire Shelves

Width/Length (in.) (mm)		Approx. Pkd. Wt. (lbs.) (kg)		Shelf Material	Handles	Cat. No.
18x24	457x610	43	19.5	Stainless Steel, Chrome	Chrome	LAB104WM
18x30	457x760	47	21.3	Stainless Steel, Chrome	Chrome	LAB204WM
18x36	457x914	52	23.6	Stainless Steel, Chrome	Chrome	LAB304WM
21x36	530x914	59	26.8	Stainless Steel, Chrome	Chrome	LAB404WM
24x36	610x914	65	29.5	Stainless Steel, Chrome	Chrome	LAB604WM



MW600 Series

Utility Carts with 2 Wire Shelves

Width/Length (in.) (mm)		Approx. Pkd. Wt. (lbs.) (kg)		Shelf Material	Handles	Cat. No.
18x24	457x610	34	15.4	Chrome	Chrome	LAB106WM
18x24	457x610	34	15.4	Stainless Steel	Stainless Steel	LAB206WM
18x30	457x760	37	16.8	Chrome	Chrome	LAB306WM
18x30	457x760	37	16.8	Stainless Steel	Stainless Steel	LAB406WM
18x36	457x914	40	18.1	Chrome	Chrome	LAB506WM
18x36	457x914	39	17.7	Stainless Steel	Stainless Steel	LAB606WM
21x36	530x914	44	20.0	Chrome	Chrome	LAB706WM
21x36	530x914	44	20.0	Stainless Steel	Stainless Steel	LAB806WM
24x36	610x914	47	21.3	Chrome	Chrome	LAB116WM
24x36	610x914	46	20.9	Stainless Steel	Stainless Steel	LAB216WM



MW700 Series

Utility Carts with 3 Wire Shelves

Width/Length (in.) (mm)		Approx. Pkd. Wt. (lbs.) (kg)		Shelf Material	Handles	Cat. No.
18x24	457x610	39	17.6	Chrome	Chrome	LAB107WM
18x24	457x610	38	17.2	Stainless Steel	Stainless Steel	LAB207WM
18x30	457x760	43	19.5	Chrome	Chrome	LAB307WM
18x30	457x760	42	19.1	Stainless Steel	Stainless Steel	LAB407WM
18x36	457x914	47	21.3	Chrome	Chrome	LAB507WM
18x36	457x914	45	20.4	Stainless Steel	Stainless Steel	LAB607WM
21x36	530x914	50	22.7	Chrome	Chrome	LAB707WM
21x36	530x914	50	22.7	Stainless Steel	Stainless Steel	LAB807WM
24x36	610x914	60	27.7	Chrome	Chrome	LAB117WM
24x36	610x914	57	25.9	Stainless Steel	Stainless Steel	LAB217WM

Super Adjustable Super Erecta® Stem Caster Carts — 11.01a

- Four- and five-tier models available.
- Consists of Super Adjustable Super Erecta wire shelves.
- 5" (127mm) resilient rubber or polyurethane casters with donut bumpers.
- Super Adjustable Super Erecta shelves allow for quick and easy adjustability.
- Open-wire shelf design minimizes dust and increases air circulation and visibility.
- Strong and versatile. (See Metro Fact below)
- For additional storage levels, Super Adjustable shelves may be added as desired. (See page 37.)



Super Adjustable Super Erecta® Wire Stem Caster Cart



Lift the release at each corner to adjust Super Adjustable shelves at 1" (25mm) increments in seconds.

Metro Fact:

Carts with polyurethane casters are designed to hold up to 900 lbs. (408kg). Carts with rubber casters can hold up to 600 lbs. (272kg) of evenly distributed weight. Heavier weight loads should be stored as low as possible on cart for safe maneuverability.

Use a solid shelf on the lowest levels to protect contents on the bottom of the cart from debris and dripping. (See page 67).

4-Tier Models

With Chrome Shelves and Chrome-Plated Posts

Overall Height 67⁷/₈" (1724mm)

Shelf Width/Length (in.)	Shelf Width/Length (mm)	Approx. Pkd. Wt. (lbs.) (kg)	Catalog Number with Casters	
			Two Swivel Two Brake Resilient Rubber Casters	Two Swivel Two Brake Polyurethane Casters
18x36	457x914	60 27	LABCB633A	LABCE633A
18x48	457x1219	72 32	LABCB653A	LABCE653A
18x60	457x1524	88 40	LABCB663A	LABCE663A
21x36	530x914	69 31	LABCB643A	LABCE634A
21x48	530x1219	81 36	LABCB654A	LABCE654A
21x60	530x1524	98 44	LABCB664A	LABCE664A
24x36	610x914	76 34	LABCB635A	LABCE635A
24x48	610x1219	92 41	LABCB655A	LABCE655A
24x60	610x1524	112 51	LABCB665A	LABCE665A

Note: Models include 63UP posts.

5-Tier Models

With Chrome Shelves and Chrome-Plated Posts

Overall Height 67⁷/₈" (1724mm)

Shelf Width/Length (in.)	Shelf Width/Length (mm)	Approx. Pkd. Wt. (lbs.) (kg)	Catalog Number with Casters	
			Two Swivel Two Brake Resilient Rubber Casters	Two Swivel Two Brake Polyurethane Casters
18x36	457x914	60 27	LABCB633A5	LABCE633A5
18x48	457x1219	72 32	LABCB653A5	LABCE653A5
18x60	457x1524	88 40	LABCB663A5	LABCE663A5
21x36	530x914	69 31	LABCB634A5	LABCE634A5
21x48	530x1219	81 36	LABCB654A5	LABCE654A5
21x60	530x1524	98 44	LABCB664A5	LABCE664A5
24x36	610x914	76 34	LABCB635A5	LABCE635A5
24x48	610x1219	92 41	LABCB655A5	LABCE655A5
24x60	610x1524	112 51	LABCB665A5	LABCE665A5

Note: Models include 63UP posts.