



XT Series

**Extended Time Cryogenic
Small Freezers**

TW 3XTL

LN₂ EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	3XTL
LN ₂ Capacity (L)	3
Static Hold Time (Days)	27
Evap. Rate (Liters/Day)	0.11
DIMENSIONS	
Necktube (in)	2.00
Overall Dia. (in)	7.6
Height (in)	17.2
WEIGHT	
Empty (Lbs)	7.2
Full Weight (Lbs)	12.5
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	N/A
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	N/A



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

TW 8XTL

LN₂ EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid



SPECIFICATIONS	8XTL
LN ₂ Capacity (L)	8
Static Hold Time (Days)	80
Evap. Rate (Liters/Day)	0.1
DIMENSIONS	
Necktube (in)	2.00
Overall Dia. (in)	15.6
Height (in)	19
WEIGHT	
Empty (Lbs)	19.6
Full Weight (Lbs)	33.8
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	N/A
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	N/A



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

TW 10XT

LN₂ EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	10XT
LN ₂ Capacity (L)	10
Static Hold Time (Days)	100
Evap. Rate (Liters/Day)	0.1
DIMENSIONS	
Necktube (in)	2.00
Overall Dia. (in)	11.4
Height (in)	23.8
WEIGHT	
Empty (Lbs)	16.5
Full Weight (Lbs)	34.3
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	540
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	150



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

TW 20XT

LN₂ EXTENDED TIME SMALL FREEZER



For long-term storage at cryogenic temperatures.



SPECIFICATIONS	20XT
LN ₂ Capacity (L)	20.7
Static Hold Time (Days)	230
Evap. Rate (Liters/Day)	0.09
DIMENSIONS	
Necktube (in)	2.00
Overall Dia. (in)	15.6
Height (in)	25.8
WEIGHT	
Empty (Lbs)	26
Full Weight (Lbs)	62.8
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	540
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	150

FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid



800.886.2796 | Fax 317. 297.7988
 sales@ICBiomedical.com
 www.ICBiomedical.com

TW 34XT

LN₂ EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	34XT
LN ₂ Capacity (L)	34
Static Hold Time (Days)	340
Evap. Rate (Liters/Day)	0.1
DIMENSIONS	
Necktube (in)	2.00
Overall Dia. (in)	18.8
Height (in)	26.3
WEIGHT	
Empty (Lbs)	34.75
Full Weight (Lbs)	95.3
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	540
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	150



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com