

TW XT21-AI

AI SMALL FREEZER




STRONGER. STURDIER. AND BUILT TO PERFORM.

Our Extended Time animal insemination small freezer is designed for long-term storage at cryogenic temperatures under the most challenging conditions.

Innovative features that work harder for you:

- Rugged Construction – Ribbed high-strength aluminum body
- Durable Design – Improved magneformed neck tube for demanding field use
- Convenient Storage – Numbered canister index ring and tapered internal spider
- Superior Vacuum Performance – Added super insulation for maximum holding times
- Reliable – Designed and built to stringent specifications, backed by a five-year vacuum warranty



Enhance your XT21-AI with optional accessories:

- Liquid Nitrogen low-level alarm
- 5" canisters available
- Convenient roller base

*See back page
for product
specifications*

TW XT21-AI
AI SMALL FREEZER



MODELS	AI21
Static Holding Time Days	221
Working Time Days	138
Evap. Rate (Liters/Day)	0.95
Liquid Nitrogen Capacity (liters)	21
Weight Empty (lbs./kg.)	26 (11.8)
Weight Full (lbs./kg.)	63.4 (28.8)
Neck Diameter(in./mm.)	2.18 (55.4)
Overall Height (in./mm.)	25.8 (655)
Overall Diameter (in./mm.)	15.6 (396)
Number of Canisters	6
Canister Dimensions (in./mm.)	1.65x11 (41.9x279)
Number of 1.2 ml & 2.0 ml vials - (6/cane)	210
Number of 1/2 cc straws - (10/cane)	720
Number of 1/2 cc straws - Bulk (1 level)	1122 (1.65 x 5 {41.9 x 127} Canister)
OPTIONAL ACCESSORIES	
Roller Base	R018-8C00
Low-Level Alarm	R021-8C15



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
 Made In USA