







## AGCN Advanced Series Analytical

With an outstanding price-to-performance ratio the AGCN is ideal for basic analytical weighing applications. The balance combines a high quality electromagnetic load cell with a complete die-cast metal housing, to deliver unmatched performance and reliability in all settings.

### Functions

					
PARTS COUNTING	TOTALIZING	PERCENT WEIGHING	AUTO CALIBRATION	RS232	DENSITY MEASURING

### Features

- Robust die-cast metal housing
- Automatic Internal Calibration
- Easy access weighing chamber
- Backlit LCD display
- Density Measurement Feature
- Spill proof sealed front panel with dual tare keys
- RS-232 output



Spacious and easily accessible weighing chamber with three sliding doors.



MODEL	AGCN100	AGCN200
CATALOG No.		
Capacity	100g	200g
Readability (d)	0.0001g	
Linearity	+/-0.0002g	
Stabilization	~ 3sec	
Accuracy Class	I	
Pan Size (mm)	Ø 3.5" (90 mm)	
Scale Dimensions	8.5" x 13.5" x 13.6" (216 x 345 x 347 mm)	
RS232 Port	Bidirectional	