


On TARGET 2024 Promotions

Vol. 1-2024
special offers valid
through Oct. 31, 2024



isoPURE™ Mini

Automated Purification System

- Compatible with all popular magnetic bead kits for DNA, RNA purification
- Process up to 16 samples
- Compact, 8x10 in. footprint
- Programmable by PC or Android app




Capacity:	16 samples
Plate Positions:	1 Plate Position
Lysis Temp. Range:	Ambient +5° to 120°C
Elution Temp. Range:	Ambient +5° to 120°C
Mixing Speed:	Variable (10 mixing levels)

AP1016 IsoPure™ Mini, Automated Purification System **\$10,447**



Special Offer \$337 value



Purchase AP1016 and receive a free starter kit, including:

- One case AP1016-MC Tip combs
- One case AP1096-DWP Deep well plates

Also Available

isoPURE™ 96

Automated Purification System

\$47,250


Item:	AP1096
Capacity:	96 samples
Plate Positions:	8 Position carousel
Lysis Temp. Range:	Ambient +5° to 120°C
Elution Temp. Range:	Ambient +5° to 120°C
Mixing Speed:	Variable (10 mixing levels)



Smart-Q™ 100/800

Quantitative Fluorometer *For quantification of DNA, RNA & Protein*

- Sensitive quantitation down to 0.1pg/µl ds DNA
- Large color touch screen for programming and operation
- Two models for single or 8-sample processing
- Built-in reagent calculators for sample preparation
- Compact and portable (11 x 6 inch footprint)



Display:	7.0-inch, color touch screen
Dynamic Range:	5 orders of magnitude
Light Source:	Blue and Red LED
Consumable type:	0.5mL assay tubes

Choose 1 FREE Quantification Kit with purchase of F1001 OR F1008

F1000-HS-100	High Sensitivity dsDNA Quant. Kit
F1000-BR-100	Broad Range dsDNA Quant. Kit

Value: \$131

Use promo code: **OT1.24**

F1001	Smart-Q™ 100 Fluorometer for single samples	\$3,900
F1008	Smart-Q™ 800 Fluorometer for 8-samples	\$6,650

On TARGET



Accuris Analytical Balances

- Readability of 0.0001g
- Available with Quick-Cal™ internal calibration feature
- Multiple mode options (standard weighing, count mode, & percent)

Analytical Balances			Compact Balances		
W3100-120	120g x 0.0001g	\$2,143	W3300-120	120g x 0.01g	\$327
W3100A-120	120g x 0.0001g w/ QuickCal™	\$2,274	W3300-500	500g x 0.01g	\$381 Free!*
W3100-210	210g x 0.0001g	\$2,488	W3300-1200	1200g x 0.01g	\$424
W3100A-210*	210g x 0.0001g w/ QuickCal™	\$2,637			<i>With purchase of W3100A-210</i>

***FREE!**
Compact Balance

With purchase of W3100A-210

Accuris Compact Balances

- Readability of 0.1g to 0.01g
- Easy in-lab calibration feature
- RS232 port for data output to a computer or printer

Use promo code: **OT1.24**



QüadCount™

Automated Cell Counter

- Reduce counting errors, save time and improve efficiency
- Intuitive software and touch-screen interface
- Utilizes cost-effective, disposable, 4-chamber QuadSlides™
- Quick, normal and precise count modes (up to 3.6µl count volume)

E7500	QuadCount™	\$4,275
E7500-S1	QuadSlide™, 50/pk.	\$172

Sale Price: \$4,075

Use promo code: **OT1.24**



Replaces manual cell counting...ask for a demo today!



NEW PRODUCT

Accuris UV Transilluminator

- Large, 21x25cm viewing surface
- Powerful 14 watt bulbs
- Variable intensity
- 302nm midrange UV wavelength
- Optimal for EtBr, SmartGlow™, and SYBR stains
- Hinged, UV-blocking cover



Intro Price \$1,250

E3200 UV Transilluminator, Variable Intensity **\$1,250**