

streamlined

SAMPLE PREP

SHAKE AND STERILIZE



Buy any Lab Line Autoclave
and receive a **FREE**
Hei-SHAKE Rotamax 120



Save time and ensure accurate results! Sample preparation is key to ensuring reproducible and accurate results in your laboratory. Streamline your sample preparation with the combination of Heidolph Tuttnauer's premium lab line sterilizers and the reliability of the Hei-SHAKE Rotamax 120. **Products backed by a 3 year warranty!**



A Powerful Combination

Rotamax 120 Orbital Shaker

Ensure uniform mixing and consistent reactions with the variably adjustable speed control (20-300 rpm) and the orbit of 20 mm of the Hei-SHAKE Rotamax 120 Orbital Shaker.

Lab Line Autoclaves

Sterilize everything in your lab including liquids with the versatile lineup of Heidolph Tuttnauer Lab Line Autoclaves. Choose from benchtop, vertical, free standing configurations as you prevent contamination in your laboratory!

HEIDOLPH TUTTNAUER LAB LINE AUTOCLAVES



Liquid



Plastic



Glass



Sealed bottle



Tools



Wrapped Tools



Clothes



Bioreactors

Benchtop Lab Line

Model	Chamber Dimensions (in)	Volume (L)	P/N
2840ELP	11 x 18	28	023210671
3850ELP	15 x 19.7	65	023210681
3870ELP	15 x 27.2	85	023210672
5075ELP	19.7 x 29.5	160	023210675



Vertical Lab Line

Model	Chamber Dimensions (in)	Volume (L)	P/N
3950V	15.3 x 19.8	60	023210691
3970V	15.3 x 27.6	85	023210692
2840ELVP-WR	11 x 18	28	023210744
3870ELVP-WR	15 x 27.2	85	023210745
5075ELVCP	19.7 x 29.5	160	023210990



Free-Standing Lab Line

Model	Chamber Dimensions (in)	Volume (L)	P/N
3870EL-PVG-WS	14 x 27.2	85	023210680
5075EL-PVG-WS	19.7 x 29.5	160	023210685
2840ELC-PVG-P	11 x 18	28	023210683
3870ELC-PVG-P	15 x 27.2	85	023210678
5075ELC-PVG-P	19.7 x 29.5	160	023260002



Hei-SHAKE ORBITAL SHAKER

Rotamax Orbital Shaker

Model	Platform Dimensions	Orbital Diameter	P/N
Hei-SHAKE Rotamax 120	220 x 220 mm	20 mm	036130160



Cell Cultivation



Fermentation Process



Washing Steps



Staining



Dissolving Samples



Heidolph North America

heidolph

research made easy

