



BeadBug™ and BeadBlaster™ Homogenizers

- Completely lyse, grind & homogenize a variety of sample types
- Available with **REFRIGERATION** to protect delicate samples
- Process up to 24 samples simultaneously
- Wide variety of grinding media in prefilled 2ml & 5ml tubes

Item No.	Description	Price
D1030*	BeadBug 3 Homogenizer	\$998
D1036*	BeadBug 6 Homogenizer	\$2,736
D2400*	BeadBlaster 24 Homogenizer	\$9,460
D2400-R*	BeadBlaster 24R Refrigerated Homogenizer	\$15,098
D1032-05	Prefilled 2ml tubes, 0.5mm zirconium, 50/pk.	\$184
D1032-15	Prefilled 2ml tubes, 1.5mm zirconium, 50/pk.	\$184
D1032-30	Prefilled 2ml tubes, 3mm zirconium, 50/pk.	\$184
D1033-28	Prefilled 2ml tubes, 2.8mm steel, 50/pk.	\$184
D1034-28	Prefilled 5ml tubes, 2.8mm steel, 50/pk.	\$246

Additional grinding media is available in prefilled tubes and loose bulk. Ask for details.



Special Offer:
FREE Prefilled Tube Kit (50pk.)
*Choose 1 free pack of prefilled tubes w/ purchase



Incu-Shaker™ Series Shaking Incubators

- MAGic™** Clamp accessories for instant-magnetic clamp exchange
- Available with **CO₂** control & cooling
- NEW!** Adjustable clamp for 125, 250 and 500ml flasks

Special Offer:
FREE MAGic Clamp Platform

Receive H1000-MR or H1010-MR w/ purchase of Incu-Shaker Shaking Incubator



H1000-MR-CMB
Adjustable clamp for 125, 250 & 500ml flasks

Incu-Shaker Mini

Speed Range:
Temp Range:
Platform Dims:
Universal Platform:
Flask Capacity:

30 to 300rpm
Amb. +5° to 70°C
9.5 x 11.5 in.
H1000-MR
2 x 1L (H1000-MR-1000)
5 x 500ml (H1000-MR-500)
8 x 250ml (H1000-MR-250)
12 x 125ml (H1000-MR-125)

\$3,635



Item: H1001-M

Incu-Shaker 10L & 10LR

Speed Range:
Temp Range:
(H2012, Cooling Model)
Platform Dims:
Universal Platform:
Flask Capacity:

30 to 300rpm
Amb. +5° to 60°C
Amb. -15° to 60°C
17.8 x 17.8 in.
H1010-MR
5 x 2L (H1000-MR-2000)
9 x 1L (H1000-MR-1000)
16 x 500ml (H1000-MR-500)
25 x 250ml (H1000-MR-250)
36 x 125ml (H1000-MR-125)

**\$7,570/
\$10,701**



Item: H2010 & H2012 (Cooling)

Incu-Shaker Mini CO₂

Speed Range:
Temp Range:
CO₂ Range:
Platform Dims:
Universal Platform:
Flask Capacity:

30 to 300rpm
Amb. +5° to 60°C
0 to 20%
9.5 x 11.5 in.
H1000-MR
4 x 1L (H1000-MR-1000)
5 x 500ml (H1000-MR-500)
8 x 250ml (H1000-MR-250)
12 x 125ml (H1000-MR-125)

\$11,120



Item: H3501

Choose one
FREE

With purchase of item below*

(Free with items: BSH6000, H3715, B5200, E7500, H5100-HCT, BT4001/11)



Value: \$462

Includes COMBI-rotor for tubes & strips



Value: \$335

Speed Range: 200-3200 rpm



Value: \$441


Includes 4 tube holders

Your choice of:

MyFuge™ 12
Item: C1012

Benchmixer™ V2
Item: BV1003

RotoMini™
Item: R2020



ISO BLOCK™ Two Independent Block Chambers


- Two isolated chambers with independent temperature control
- Digital control up to 105°C
- Hinged lid improves uniformity
- Quick-Flip universal block accepts 24 x 1.5ml, 48 x 0.2ml, 14 x 0.5ml

Additional blocks available for all common tubes & vials

BSH6000	isoBlock™ Digital Dry Bath, two independent chambers	\$945
BSWCMB	Quick-Flip™ Block, 24 x 1.5ml, 48 x 0.2ml, 14 x 0.5ml	\$203
BSW02	Block, 48 x 0.2ml tubes or 6 PCR strips	\$152

HS-10 Plus Multi-functional Hotplate Stirrer

- Multi-functional design, simultaneously heat & stir up to 5 vessels
- Large, 10 x 10 inch white, ceramic top plate
- Digital speed and temperature control
- Safety LED indicates hot surface



H3715 Benchmark HS-10 Plus Hotplate Stirrer **\$1,049**



SUREAIR™ PCR WORKSTATION WITH UV AND HEPA FILTRATION

- HEPA filter and UV light reduce cross-contamination during PCR preparation
- Compact, 55 x 40 cm footprint
- Cable port allows easy use of centrifuges, vortexers and other instruments
- Includes pipette holder and shelf for storage of pipettes and tips
- Provides clean air (to ISO 5 standards)

B5200	SureAir™ PCR Workstation	\$2,325
Accessories:		
B5200-FIL	SureAir™ Replacement Filter	\$314
B1450-IC	UV-C Indicator Cards, 25 ct.	\$111
B1450-ISS	UV-C Indicator Stickers, 40 ct	\$99

QuadCount™ Automated Cell Counter for mammalian cells

- Fast, accurate and reproducible counting & live/dead determinations
- Intuitive software with auto-focus and touch-screen interface
- Cost-effective QuadSlides™
- Disposable slides with four individual chambers



E7500 QuadCount Cell Counter **\$4,361**
E7500-S1 QuadSlide, 4 chamber slides **\$175**

MULTI-THERM™ Touch COOL - HEAT - SHAKE

- Large color touch screen
- Multiple block options for plates & tubes up to 50ml
- Broad temperature & speed range
- Advanced programming features

Speed range: 200-3,000 rpm (block dependent)
Shaker Orbit: 3mm
Temp. range: Ambient -20° to 105°C



Heated lid prevents condensation & evaporation

ORBi-SHAKER™ CO₂ Remote Controlled, All Environment Shakers

- Designed to fit & maximize capacity in most CO₂ incubators & to withstand temperatures 0° to 65°C
- Horizontally circular 19mm orbit for aeration & mixing
- Adjust settings via the **included** remote control module
- Instantly exchange flask clamps with, **MAGic Clamp™** accessories

Speed range: 30-300 rpm
Max capacity: 2 x 2L or 5 x 1L flasks (BT4001)
5 x 2L or 9 x 1L flasks (BT4011)



BT4001 Orbi-Shaker CO₂, Shaker with remote control, includes 13x12" non-slip platform **\$3,879**
BT4011 Orbi-Shaker CO₂, XL Shaker with remote control, includes 16x16.5" non-slip platform **\$5,379**

H5100-HCT	MultiTherm Touch, Heating/Cooling Vortexer (blocks sold separately)	\$4,231	Available blocks* H5100-05 24 x 0.5ml	\$378
H5100-HL	Optional heated lid	\$811	H5100-MP Block, microplate	\$378
			H5100-384 Block, 384 well plate	\$648

*Additional blocks available, see website for details